dedicio®

# COVID-19 Ag plus Test

#243104D-20

|           |                                |    |
|-----------|--------------------------------|----|
| <b>de</b> | Gebrauchsanweisung .....       | 2  |
| <b>en</b> | Instructions for use.....      | 7  |
| <b>fr</b> | Instructions d'utilisation.... | 12 |
| <b>es</b> | Instrucciones de uso .....     | 17 |
| <b>it</b> | Istruzioni per l'uso .....     | 22 |
| <b>pl</b> | Sposób użycia .....            | 27 |
| <b>pt</b> | Instruções de Utilização ....  | 32 |
| <b>cs</b> | Návod k použití .....          | 37 |
| <b>fi</b> | Käyttöohje.....                | 42 |
| <b>sv</b> | Användarinstruktioner .....    | 47 |
|           | Symbols .....                  | 52 |

**dedicio+**  
dedicated to life

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**15. Referenzen**

1. Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, Nat Rev Microbiol 2019; 17:181-192.
2. Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, TrendsMicrobiol 2016;24:490-502.
3. Weiss SR, Leibowitz JL, Coronavirus pathogenesis, Adv Virus Res 2011; 81:85-164.

Rev. 3, 2021-01-18 OM

### 1. Intended Use

The dedicio® COVID-19 Ag plus Test is a lateral flow chromatographic immunoassay for the qualitative detection of SARS-CoV-2 viral nucleoprotein antigens in human nasal, nasopharyngeal or oropharyngeal specimens (see section 12 'Limitations'). This test is intended for use as an aid in the diagnosis of infections with SARS-CoV-2. Note that the concentration of viral nucleoprotein antigens may vary in the course of the disease and might fall below the detection limit of the test. Possible infectiousness of test subjects cannot be ruled out based on negative test results. The test procedure is not automated and requires no special training or qualification. The dedicio® COVID-19 Ag plus Test is designed for professional use only.

### 2. Introduction and Clinical Significance

COVID-19 (Corona Virus Disease) is the infectious disease caused by the recently discovered coronavirus SARS-CoV-2. The most common symptoms of COVID-19 are fever, dry cough, fatigue, sputum production, shortness of breath, sore throat and headache. Some patients may have myalgia, chills, nausea, nasal congestion and diarrhoea. These symptoms begin gradually and are mild in most of the cases. Some people become infected but do not develop any symptoms and do not feel unwell. Most people (about 80%) recover from the disease without special treatment. Approximately one in six people who get infected with COVID-19 becomes seriously ill and develops difficulty breathing. Elderly people, and those with pre-existing conditions, such as high blood pressure, heart problems or diabetes, are more likely to develop serious illness. So far, about 2% of infected people have died.

COVID-19 is transmitted via respiratory droplets that are exhaled by infected people via coughing, sneezing or talking. These droplets can be inhaled or ingested directly by other people or can contaminate surfaces, which can then be infectious for several days. Most estimates of the incubation period for COVID-19 range from 1 to 14 days, during which people might already be infectious without showing disease symptoms.

### 3. Test Principle

The dedicio® COVID-19 Ag plus Test is a lateral flow chromatographic immunoassay for the qualitative detection of SARS-CoV-2 viral nucleoprotein antigens in human nasal, nasopharyngeal or oropharyngeal specimens.

Anti-SARS-CoV-2 antibodies are immobilised in the test line region (T) of the membrane. A specimen is added to an extraction tube containing buffer in order to release SARS-CoV-2 antigens. During the test, extracted antigens bind to anti-SARS-CoV-2 antibodies conjugated to coloured particles and precoated onto the sample pad of the test cassette. The mixture then migrates along the membrane chromatographically by capillary action and interacts with the reagents on the membrane. The complexes are then captured by anti-SARS-CoV-2 antibodies in the test line region (T). Excess coloured particles are captured in the control line region (C). The presence of a coloured line in the test line region (T) indicates a positive result. The absence of a coloured line in the test line region (T) indicates a negative result.

The formation of a coloured line in the control line region (C) serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

### 4. Reagents and Materials Supplied

- 20 dedicio® COVID-19 Ag plus test cassettes
- Additional material provided according to 93/42/EEC:
  - Due to possible supply shortages of COVID-19 related accessory medical products, the swab manufacturer might change. Therefore, the supplied swabs are from one of the manufacturers listed below.

#### a) 20 sterile swabs, CE0197

Jiangsu Changfeng Medical Industry Co., Ltd  
Touqiao Town, Guangling District, Yangzhou,  
Jiangsu 225109 China (authorised EU  
representative Llins Service & Consulting GmbH,  
Obere Seegasse 34/2, 69124 Heidelberg,  
Germany)

#### b) 20 sterile swabs, CE0197

Jiangsu Rongye Technology Co., LTD, Touqiao  
Town, Yangzhou City, Jiangsu Province, China  
(authorised EU representative Riomavix S.L., Calle  
de Almansa 55, 1D, Madrid 28039 Spain)

#### c) 20 sterile swabs, CE0197

CITOTEST LABWARE MANUFACTURING CO., LTD  
No.48, Xinxu Road, Haimen, Jiangsu province  
(authorised EU representative WellKang Ltd,  
Enterprise Hub, NW Business Complex,  
1 Beraghmore Rd., Derry, BT48 8SE, Northern  
Ireland)

- 20 extraction tubes incl. dropper caps
- 20 buffer ampoules for single use (300 µL each)\*
- 1 reagent holder
- 1 package insert

\*Buffer containing the following preservative: ProClin™ 300: <0.03%.

The detergents contained in the buffer lyse and neutralise the virus.

No hazard labelling is required according to Regulation (EC) N° 1272/2008 CLP. Concentrations are below the exemption threshold.

### 5. Additional Materials Required

- Timer

### 6. Storage & Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze the test kit. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

## 7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not use the buffer if it is discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear appropriate personal protective equipment, such as face mask, isolation gown, gloves and eye protection during specimen collection, preparation and test procedure.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- Further specimen processing and patient management should follow local COVID-19 guidelines and regulations.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Temperature can adversely affect test results.
- Used testing materials should be disposed of according to local regulations.

## 8. Specimen Collection and Preparation

### Nasal specimen:

- It is important to obtain as much secretion as possible. Insert the swab into the nostril.
- Gently push the swab until resistance is met at the level of the turbinate (up to 2.5 cm into the nostril).
- Rotate the swab 5 times against nasal wall to ensure that both mucus and cells are collected.
- Slowly withdraw the swab while continuing to rotate it.
- Repeat this process for the other nostril using the same swab in order to ensure that sufficient specimen volume is collected from both nasal cavities.

### Nasopharyngeal specimen:

- Insert the swab into the nostril, parallel to the palate (not upwards) until resistance is encountered or the distance is equivalent to that from the ear to the nostril of the patient, indicating contact with the nasopharynx.

- Gently rub and roll the swab. Leave swab in place for several seconds to absorb secretions.
- Slowly remove the swab while rotating it. Specimens can be collected from both nostrils using the same swab, but it is not necessary to collect specimens from both sides if the tip is saturated with fluid from the first collection.

### Oropharyngeal specimen:

- Gently insert a sterile swab into the pharynx and collect secretions by brushing the swab several times against the reddened posterior pharyngeal wall and both tonsillar pillars. Avoid touching the tongue, teeth and gums.


### Note:

- Use only synthetic fibre swabs with plastic shafts. Do not use calcium alginate swabs or swabs with wooden shafts, as they may contain substances that inactivate some viruses and inhibit further testing.
- Swab specimens should be tested immediately after collection. Use freshly collected specimens for best test performance. If the test cannot be carried out immediately after specimen collection, the swab can be stored for up to a maximum of one hour in a clean, dry and sealed tube.
- Do not use specimens that are obviously contaminated with blood, as it may interfere with the flow of specimens and lead to inaccurate test results.

## 9. Test Procedure

Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.

1. Place a clean extraction tube labeled with the patient or control identification into the designated area of the reagent holder.
2. Detach one buffer ampoule.
3. Open the ampoule by twisting the tip off.
4. Hold the buffer ampoule vertically over the tube to ensure that the entire buffer solution flows into the bottom part. Squeeze the buffer ampoule and add the entire buffer solution to the extraction tube without touching the edge of the tube.
5. Insert the swab with the collected specimen into the tube. Swirl the swab and squeeze it at least 5 times by compressing the wall of the extraction tube against the swab to extract the antigens contained in the swab.
6. Remove the swab, pressing it firmly against the wall of the tube to release as much liquid as possible. Dispose of the swab in accordance with guidelines for the handling of infectious agents.


**Note:** If the dropper cap is attached, the extracted specimen can be stored in the extraction tube for up to a maximum of 1 hour until the removal of the test cassette from its foil pouch.

7. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.

8. Place the test cassette on a clean and level surface.

9. Attach a dropper cap to the extraction tube, invert the tube and transfer 3 drops (approximately 60 µL) of the extracted solution to the specimen well (S) of the test cassette.

15 min


10. Start the timer.

11. Wait for the coloured line(s) to appear. Read the test result after 15 minutes. Do not interpret the result after more than 20 minutes.

### 10. Result Interpretation

#### Positive:

Two coloured lines appear in the result area. One line appears in the control line region (C) and the other line appears in the test line region (T).

**Note:** The colour intensity in the test line region (T) may vary depending on the concentration of SARS-CoV-2 viral nucleoprotein antigens in the specimen. Any shade of colour in the test line region (T) should be considered a positive result. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.

#### Negative:

Only one coloured line appears in the control line region (C). No coloured line appears in the test line region (T).

#### Invalid

The control line (C) fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

### 11. Quality Control

An internal procedural control is included in the test cassette:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

*Good laboratory practice (GLP)* recommends the use of external control materials to ensure proper test kit performance.

### 12. Limitations

- The dedicio® COVID-19 Ag plus Test is for professional *in-vitro* diagnostic use only. It should be used for the qualitative detection of SARS-CoV-2 viral nucleoprotein antigens in human nasal, nasopharyngeal or oropharyngeal specimens only. Neither the quantitative value nor the rate of increase/decrease in the concentration of SARS-CoV-2 viral nucleoprotein antigens can be determined with this qualitative test.
- The dedicio® COVID-19 Ag plus Test only detects the presence of SARS-CoV-2 viral nucleoprotein antigens in specimens and should not be used as the sole criterion for a diagnosis of COVID-19.
- Both viable and non-viable SARS-CoV-2 viruses can be detected using the dedicio® COVID-19 Ag plus Test.
- The sections 'Specimen Collection and Preparation' as well as 'Test Procedure' must be followed closely while testing. Failure to follow them may lead to inaccurate test results because the antigen concentration in the swab is highly dependent on the correct procedure.
- As with all diagnostic tests, all results should be interpreted in conjunction with other clinical information available to the physician.
- In the course of SARS-CoV-2 infection, the concentration of viral nucleoprotein antigens may fall below the detection limit of the test.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of a SARS-CoV-2 infection and should be confirmed via molecular assay.
- Positive and negative predictive values are highly dependent on prevalence. The local prevalence should be taken into consideration when interpreting diagnostic test results.
- Positive results do not preclude co-infections with other pathogens (e.g. influenza virus A/B).
- The test does not differentiate between SARS-CoV and SARS-CoV-2.

### 13. Expected Values

SARS-CoV-2 viral particles are normally present in the respiratory tracts of COVID-19 patients. A positive test result can indicate an acute infection. Virus concentrations in nasal, nasopharyngeal or oropharyngeal swab specimens may vary in the course of the disease and might fall below the detection limit of rapid tests, even though patients are still showing symptoms. Conversely, the virus might continue to be detectable over long periods of time even in convalescent patients. Possible infectiousness of test subjects cannot be ruled out based on negative test results.

**14. Performance Characteristics**

**Clinical performance**

**Diagnostic sensitivity and specificity**

The dedicio® COVID-19 Ag plus Test was evaluated with clinical nasal swab specimens whose status was confirmed using RT-PCR.

The results are presented in the following tables.

| dedicio® COVID-19 Ag plus Test |          | RT-PCR   |          |       |
|--------------------------------|----------|----------|----------|-------|
|                                |          | Positive | Negative | Total |
|                                | Positive | 167      | 0        | 167   |
|                                | Negative | 5        | 431      | 436   |
|                                | Total    | 172      | 431      | 603   |

Diagnostic sensitivity: 97.09% (93.38% - 98.75%)\*

Diagnostic specificity: >99.99% (99.12% - 100%)\*

Overall agreement: 99.17% (98.07% - 99.65%)\*

\*95% confidence interval

**Detection limit**

The 95% detection limit of the dedicio® COVID-19 Ag plus Test is 75.5 TCID<sub>50</sub>/mL and was determined via dilution of inactivated SARS-CoV-2 virus.

**Measuring range**

No prozone effect was observed when testing inactivated SARS-CoV-2 viruses with concentrations of up to 1.51 x 10<sup>6</sup> TCID<sub>50</sub>/mL.

**Interfering substances**

The following substances, normally present in respiratory specimens or artificially introduced into the respiratory tract, were evaluated at the concentrations listed below. All substances were mixed with a 1/20000 SARS-CoV-2 virus dilution (1.51 x 10<sup>6</sup> TCID<sub>50</sub>/mL) and showed no interference with the dedicio® COVID-19 Ag plus Test.

| Substance                                       | Concentration |
|---|---------------|
| Zanamivir                                       | 5 mg/mL       |
| Oseltamivir                                     | 10 mg/mL      |
| Artemether/Lumefantrine                         | 50 µM         |
| Doxycycline hyclate                             | 70 µM         |
| Quinine   | 150 µM        |
| Lamivudine                                      | 1 mg/mL       |
| Ribavirin                                       | 1 mg/mL       |
| Daclatasvir                                     | 1 mg/mL       |
| Mucin from bovine submaxillary glands, type I-S | 100 µg/mL     |
| EDTA-anticoagulated human blood                 | 5% (v/v)      |
| Biotin  | 100 µg/mL     |
| Neosynephrin® (phenylephrine)                   | 10% (v/v)     |
| Afrin® Nasal Spray (oxymetazoline)              | 10% (v/v)     |
| Saline nasal spray                              | 10% (v/v)     |
| homeopathic Zicam® Allergy Relief Nasal Gel     | 5% (v/v)      |
| Sodium cromoglycate                             | 20 mg/mL      |
| Olopatadine hydrochloride                       | 10 mg/mL      |
| Paracetamol                                     | 199 µM        |
| Acetylsalicylic acid                            | 3.62 mM       |
| Ibuprofen                                       | 2.425 mM      |
| Mupirocin                                       | 10 mg/mL      |

| Substance     | Concentration |
|---------------|---------------|
| Tobramycin    | 5 µg/mL       |
| Erythromycin  | 81.6 µM       |
| Ciprofloxacin | 30.2 µM       |

**Cross-reactivity**

Specimens spiked with the following pathogens were tested using the dedicio® COVID-19 Ag plus Test:

| Pathogen  | Concentration                            |
|---|--|
| Influenza A H1N1 virus                                | 6x10 <sup>5</sup> TCID <sub>50</sub> /mL |
| Influenza A H3N2 virus                                | 9x10 <sup>5</sup> TCID <sub>50</sub> /mL |
| Influenza A H5N1 virus                                | 8x10 <sup>4</sup> TCID <sub>50</sub> /mL |
| Influenza A H7N9 virus                                | 9x10 <sup>4</sup> TCID <sub>50</sub> /mL |
| Influenza B   | 3x10 <sup>5</sup> TCID <sub>50</sub> /mL |
| Adenovirus type 1                                     | 3x10 <sup>5</sup> TCID <sub>50</sub> /mL |
| Adenovirus type 2                                     | 2x10 <sup>5</sup> TCID <sub>50</sub> /mL |
| Adenovirus type 3                                     | 1x10 <sup>6</sup> TCID <sub>50</sub> /mL |
| Adenovirus type 5                                     | 4x10 <sup>5</sup> TCID <sub>50</sub> /mL |
| Adenovirus type 7                                     | 1x10 <sup>6</sup> TCID <sub>50</sub> /mL |
| Adenovirus type 55                                    | 4x10 <sup>5</sup> TCID <sub>50</sub> /mL |
| Respiratory syncytial virus type A/B                  | 3x10 <sup>5</sup> TCID <sub>50</sub> /mL |
| Human coronavirus 229E                                | 3x10 <sup>5</sup> TCID <sub>50</sub> /mL |
| Human coronavirus OC43                                | 1x10 <sup>5</sup> TCID <sub>50</sub> /mL |
| Human coronavirus NL63                                | 1x10 <sup>4</sup> TCID <sub>50</sub> /mL |
| Human coronavirus HKU1                                | 1x10 <sup>3</sup> ng/mL                  |
| MERS-Coronavirus Florida/USA-2/Saudi Arabia.2014      | 4x10 <sup>4</sup> TCID <sub>50</sub> /mL |
| Parainfluenza virus type 1/2/3/4                      | 1x10 <sup>5</sup> TCID <sub>50</sub> /mL |
| Rhinovirus type A16                                   | 1x10 <sup>5</sup> TCID <sub>50</sub> /mL |
| Legionella pneumophila Bloomington-2                  | 1x10 <sup>5</sup> cells/mL               |
| Legionella pneumophila 82A3105                        | 1x10 <sup>5</sup> cells/mL               |
| Mycobacterium tuberculosis K                          | 1x10 <sup>5</sup> cells/mL               |
| Mycobacterium tuberculosis Erdman                     | 1x10 <sup>5</sup> cells/mL               |
| Mycobacterium tuberculosis HN878                      | 1x10 <sup>5</sup> cells/mL               |
| Mycobacterium tuberculosis CDC1551                    | 1x10 <sup>5</sup> cells/mL               |
| Mycobacterium tuberculosis H37Rv                      | 1x10 <sup>5</sup> cells/mL               |
| Streptococcus pneumoniae 475298 [Maryland (D1) 6B-17] | 1x10 <sup>5</sup> cells/mL               |
| Streptococcus pneumoniae 178 [Poland 23F-16]          | 1x10 <sup>5</sup> cells/mL               |
| Streptococcus pneumoniae 262 [CIP 104340]             | 1x10 <sup>5</sup> cells/mL               |
| Streptococcus pneumoniae Slovakia 14-10 [29055]       | 1x10 <sup>5</sup> cells/mL               |
| Streptococcus pyogenes T1                             | 1x10 <sup>5</sup> cells/mL               |
| Mycoplasma pneumoniae mutant 22                       | 1x10 <sup>5</sup> cells/mL               |
| Mycoplasma pneumoniae FH strain of Eaton Agent        | 1x10 <sup>5</sup> cells/mL               |
| Mycoplasma pneumoniae M129-B7                         | 1x10 <sup>5</sup> cells/mL               |

No cross-reactivity with the specimens was observed when tested using the dedicio® COVID-19 Ag plus Test. SARS-CoV positive specimens (Urbani strain, 1 x 10<sup>6</sup> PFU/mL) show cross-reactivity with the dedicio® COVID-19 Ag plus Test.

## Precision

### Repeatability and reproducibility

Precision was established by testing 20 replicates of negative, low positive, medium positive and high positive controls. Repeatability was established within the reproducibility study. Testing was performed by 5 operators using 3 independent dedicio® COVID-19 Ag plus test lots at 2 different sites on 3 separate days.

The dedicio® COVID-19 Ag plus Test demonstrated acceptable repeatability and reproducibility. The negative and positive values were correctly identified >99% of the time.

## 15. References

1. Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, Nat Rev Microbiol 2019; 17:181-192.
2. Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, TrendsMicrobiol 2016;24:490-502.
3. Weiss SR, Leibowitz JL, Coronavirus pathogenesis, Adv Virus Res 2011; 81:85-164.

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#### 15. Bibliographie

1. Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, Nat Rev Microbiol 2019; 17:181-192.
2. Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, Trends Microbiol 2016;24:490-502.
3. Weiss SR, Leibowitz JL, Coronavirus pathogenesis, Adv Virus Res 2011; 81:85-164.

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a) 20 hisopos estériles, CE0197

Jiangsu Changfeng Medical Industry Co., Ltd  
Touqiao Town, Guangling District, Yangzhou,  
Jiangsu 225109 China (representante EU  
autorizado Llians Service & Consulting GmbH,  
Obere Seegasse 34/2, 69124 Heidelberg,  
Germany)

b) 20 hisopos estériles, CE0197

Jiangsu Rongye Technology Co., LTD, Touqiao  
Town, Yangzhou City, Jiangsu Province, China  
(representante EU autorizado Rimovavix S.L., Calle  
de Almansa 55, 1D, Madrid 28039 Spain)

c) 20 hisopos estériles, CE0197

CITOTEST LABWARE MANUFACTURING CO., LTD  
No.48, Xinxiu Road, Haimen, Jiangsu province  
(representante EU autorizado WellKang Ltd,  
Enterprise Hub, NW Business Complex,  
1 Beraghmore Rd., Derry, BT48 8SE, Northern  
Ireland)

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a) 20 tamponi sterili, CE0197

Jiangsu Changfeng Medical Industry Co., Ltd  
Touqiao Town, Guangling District, Yangzhou,  
Jiangsu 225109 China (rappresentante autorizzato  
EU Llins Service & Consulting GmbH, Obere  
Seegasse 34/2, 69124 Heidelberg, Germany)

b) 20 tamponi sterili, CE0197

Jiangsu Rongye Technology Co., LTD, Touqiao  
Town, Yangzhou City, Jiangsu Province, China  
(rappresentante autorizzato EU Riomavix S.L.,  
Calle de Almansa 55, 1D, Madrid 28039 Spain)

c) 20 tamponi sterili, CE0197

CITOTEST LABWARE MANUFACTURING CO., LTD  
No.48, Xinxu Road, Haimen, Jiangsu province  
(rappresentante autorizzato EU WellKang Ltd,  
Enterprise Hub, NW Business Complex,  
1 Beraghmore Rd., Derry, BT48 8SE, Northern  
Ireland)

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nal von minden GmbH