
HemoCue Hb 201+

Operating Manual

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GB

Table of Contents

| | |
|--|----|
| HemoCue Hb 201 ⁺ Operating Manual | 5 |
| Components | 6 |
| Start up | 8 |
| Set-up..... | 10 |
| Audio signal, time and date..... | 10 |
| Measuring..... | 12 |
| Capillary blood..... | 12 |
| Venous or arterial blood | 20 |
| Set-up..... | 22 |
| QC-test..... | 22 |
| Memory function – Scroll | 24 |
| Delete results | 26 |
| PC/Printer function..... | 28 |
| Maintenance | 30 |
| Daily maintenance | 30 |
| Optical parts | 32 |
| Trouble Shooting Guide | 34 |
| Specifications..... | 46 |

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GB The HemoCue Hb 201+ Analyzer provides determination of haemoglobin quickly, easily and with lab quality results. Capillary, venous or arterial whole blood may be used. Being light and portable, as well as having a mains power option, the HemoCue Hb 201 system is equally suited for both static and mobile use. The analyzer stores test results, date and time for up to 600 measurements. This manual provides the basic instructions for routine use as well as technical specifications. Further information may be obtained from your HemoCue distributor or directly from HemoCue AB.

 All system components are designed and manufactured to provide maximum safety. Any other use of the system may impair the safety.

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
Components
Bestandteile
Composants
Onderdelen

GB

1. The HemoCue Hb 201+ Analyzer.*
2. Mains adapter.**
3. 4 type AA batteries, 1,5 V.***
4. A vial of HemoCue Hb 201 Microcuvettes.***
Individually packaged HemoCue Hb 201 Microcuvettes.***
5. HemoCue Hb 201+ Operating Manual, HemoCue Hb 201+ Quick Reference Guide and Instruction CD.
6. HemoCue Cleaner.

The HemoCue Hb 201+ and its mains adapter are delivered in a carton. Open the carton on a stable surface and lift out the instrument and accessories. On the bottom of the analyzer there is a lid covering the battery compartment. Press the flap to remove the lid. Place 4 type AA batteries in the battery compartment, observing the indication of polarity in the battery holders. Replace the lid. Consult local environmental authorities for proper disposal of batteries.

* Do not open the cover of the analyzer. Note: The warranty is voided if the analyzer has been opened.

**△ Only use adapters listed under specifications.

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
GB

Make sure the analyzer is placed horizontally and on a stable surface.

- 1a. If mains power is available, connect the adapter to the socket on the back.
- 1b. If no mains power is available, insert the 4 type AA batteries, 1.5 V, into the battery compartment.
- 1c. If a battery symbol appears on the display the batteries are running low. The analyzer will continue to give accurate results but the batteries should be replaced as soon as possible.
2. The analyzer can be connected to a PC (for further information, see separate manual for PC connection) or directly to a printer (see *Set-up Printer function*).
3. Pull the cuvette holder out to its loading position.
4. Press and hold the left button until the display is activated (all symbols appear on the display).
5. The display shows the version number of the programme, after which it will show "Σ", "Hb" and "🔊" (audio signal). During this time the analyzer will automatically verify the performance of the optronic unit.
6. After 10 seconds the display will show three flashing dashes and the HemoCue symbol. This indicates that the HemoCue Hb 201+ is ready for use.

To turn off the analyzer press and hold the left button until the display reads OFF and becomes blank. Disconnect the equipment from the power source by removing the AC adapter from the outlet.

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
⑥

GB

The time and date function within the set-up menu must be activated before data can be transferred to a PC.

1. Press both buttons at the same time.
2. The display now shows a flashing QC-symbol.
- 3ab. Use the right button to scroll until the display shows an audio symbol in the upper right corner. The signal can be activated or deactivated by pressing the left button.
4. Continue to scroll using the right button until the display shows characters for time, date and year. The hour figure will be flashing.
5. Use the right button to change the hours, minutes, day, month and year. To change a flashing figure hold the left button down for quick advancement or advance slowly by quickly pressing and releasing the button.
6. When all the settings are completed, hold the right button down for approximately 5 seconds. The analyzer will then automatically return to the measuring position.

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
Meting Capillair bloed

GB

Note: Always handle blood specimens with care, as they might be infectious. Always wear protective gloves when handling blood specimens.

1. After start-up, the cuvette holder should be in its loading position. The display will show three flashing dashes and the HemoCue symbol.
2. Make sure the patient's hand is warm and relaxed. Use only the middle or ring finger for sampling. Avoid fingers with rings on.
3. Clean with disinfectant and allow to dry.
4. Using your thumb, lightly press the finger from the top of the knuckle towards the tip. This stimulates the blood flow towards the sampling point.

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

GB

5. For best blood flow and least pain, sample at the side of the fingertip, not in the centre.
6. Whilst lightly pressing towards the fingertip, prick the finger using a lancet.
7. Wipe away the first 2 or 3 drops of blood.
8. Re-apply light pressure towards the fingertip until another drop of blood appears

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

GB

9. When the blood drop is large enough, fill the microcuvette in one continuous process. Do NOT refill! Note: The microcuvette should be filled within 3 minutes after the microcuvette has been taken out of its package.
10. Wipe off excess blood on the outside of the microcuvette tip. Make sure that no blood is drawn out of the microcuvette during this procedure.
11. Look for air bubbles in the filled microcuvette. If present, take a new sample. Small bubbles around the edge can be ignored.

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
GB

- Place the microcuvette into the cuvette holder and start the measurement as soon as possible but no later than 10 minutes after filling the microcuvette by gently pushing the cuvette holder to its measuring position.
- During the measurement "∞" will be shown on the display.
- The measuring time is 15-60 seconds for Hb values below 20 g/dL. The result will remain on the display as long as the cuvette holder is in the measuring position. When operating on battery power the analyzer will automatically turn off after approximately 5 minutes.
- Although the reagents are present in the microcuvette in extremely low quantities, consult local environmental authorities for proper disposal. Always handle blood specimens with care, as they might be infectious.

All results are stored in the analyzer and can be viewed in the display using the scroll function (see section *Set-up Memory function - Scroll*).

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Measuring Venous or arterial blood

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


4

GB

Note: Always handle blood specimens with care, as they might be infectious. Always wear protective gloves when handling blood specimens.

1. Mix the blood well before performing the measurement. Note, allow the blood to reach room temperature before mixing.
- 2ab. Place a drop of blood onto a hydrophobic surface using a transfer device.
3. Fill the microcuvette in one continuous process. Do NOT refill! Wipe off excess blood on the outside of the microcuvette tip. Make sure that no blood is drawn out of the microcuvette during this procedure. Note: The microcuvette should be filled within 3 minutes after the microcuvette has been taken out of its package.
4. Look for air bubbles in the filled microcuvette. If present, take a new sample. Small bubbles around the edge can be ignored.

Perform the analysis as per steps 12–15 in section *Measuring Capillary blood*.

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1. Press both buttons at the same time.
2. The display shows a flashing QC-symbol.
3. Select QC-test by pressing the left button.
4. The analyzer automatically returns to the measuring position and the QC-symbol appears in the display.
5. Mix the control solution well before performing the measurement. Note, allow the control solution to reach room temperature before mixing. Fill the microcuvette with the control solution recommended by HemoCue and carry out the measurement in accordance with steps 12–15 in section *Measuring Capillary blood*.
6. Once the measurement has been made the analyzer automatically returns to the measuring position and the QC-symbol disappears from the display. Patient sampling can begin.

If more QC-tests are to be performed, proceed in the same way following steps 1–6. To deactivate the QC-test, follow steps 1–2 and scroll using the right button until another set-up activity is shown on the display. Hold the right button down for approximately 5 seconds. The analyzer automatically returns to its measuring position and the QC-symbol is no longer visible.

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
Set-up Geheugenfunctie – Scrollen

GB

The analyzer automatically stores up to 600 results. When the memory is full the analyzer will automatically overwrite the oldest result. To view results stored in the memory use the following procedure.

- 1ab. The cuvette holder may be in the loading position or in the measuring position. If in the measuring position the most recent measurement is displayed.
2. Press the left or right button to scroll backwards or forwards between the results. The stored values are now shown in the display. When the button is held down the user can see the order of the results (the most recent measurement is always number 1) and when the button is released the result is shown.
- 3ab. If no buttons are pressed whilst in the loading position the analyzer automatically resets after 5 seconds and new measurements can be made. If in the measuring position pull out the cuvette holder and wait for the HemoCue symbol after which a new measurement can be made.

Note: The analyzer will not respond to other commands while a measurement is being made. If the error code E35 is shown while scrolling through the results, the measurement has not been stored in the memory, see the *Troubleshooting Guide*.

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
Set-up Resultaten wissen

GB

1. Press both buttons at the same time.
2. The display now shows a flashing QC-symbol.
3. Scroll using the right button until the display shows a flashing waste bin in the left corner.
- 4ab. To delete the most recent result press and hold the left button. A warning triangle together with a "1" appears in the display for approximately 3 seconds. The most recent result has now been deleted and the display shows a "0" to confirm that this has been done. The most recent result can only be deleted immediately after the measurement has been made.
- 5ab. To delete all stored results hold both buttons down at the same time for 5 seconds. The display shows the waste bin, the warning triangle and the total number of results stored. Continue to hold both buttons down until the number changes to "000" and the warning triangle disappears. All stored results have now been deleted.

Note: Deleted results cannot be recovered. If results are to be saved they must be printed out or downloaded to a PC before starting the deleting process.

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


Set-up Pc/printerfunctie

GB

The analyzer can be on or off when connected to the printer. Only the most recent result can be printed out directly after the measurement.

1. Connect the cable* to the analyzer and printer*.
2. Perform the analysis by following the steps 12–15 in section *Measuring Capillary blood*.
3. When the result is shown on the display, the printer will automatically print the result, date and time.

* Not included.

Note: Only use ASCII printers recommended by HemoCue. In order to connect the analyzer to a PC, a special PC program is required.

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Maintenance Maintenance quotidienne

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GB

The cuvette holder should be cleaned after each day of use.

1. Pull the cuvette holder out to the loading position.
2. While pressing the catch, carefully rotate the cuvette holder sideways in open position as far as possible to the left.
3. Remove the cuvette holder from the analyzer, it will come off the stainless steel pin it rotates on.
4. Clean the cuvette holder with alcohol (20-70 %) or mild detergent

Wait 15 minutes before replacing the cuvette holder and using the analyzer. Make sure the cuvette holder is dry before inserting.

The optical parts should be cleaned when directed to do so in the Troubleshooting Guide, see section *Maintenance Optical parts*.

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

Maintenance Optical parts

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GB

Dirty optical parts may cause an error code. Follow step 1-3 under Daily maintenance and then clean optical parts as follows. Note! Make sure that the HemoCue® Cleaner reaches both upper and lower cover glasses, see pictures 5a and 5c.

5. With the cuvette holder removed from the analyzer push the Cleaner into the opening of the optic unit, as far in as possible. Move from side to side 5-10 times, see picture 5b and thereafter push in and pull out the Cleaner 5-10 times, cleaning the cover glasses, placed to the left, see picture 5c. If the Cleaner is stained, repeat with a new Cleaner.
6. Wait 15 minutes before replacing the cuvette holder and using the analyzer. Make sure the cuvette holder is dry before inserting.

The cover may be cleaned with alcohol (20-70 %) or mild detergent.

Note! As an alternative to the HemoCue® Cleaner, a lint free cotton swab, non-pretreated, moistened with alcohol (20-70 % without additive) or water may also be used. If a cotton swab is used make sure it is not too wet and not too dry. Use a dry swab to wipe away excess liquid in the optic house after cleaning with a moistened swab. To avoid scratches on the cover glasses, only the cotton part of the swab should come in contact with the cover glasses.

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GB Trouble Shooting Guide

If you are unable to resolve the problem by following this Trouble Shooting Guide, please contact your local HemoCue distributor or HemoCue AB. Note: Do not open the cover of the analyzer. The warranty is voided if the analyzer has been opened.

| Symptom | Explanation | Action |
|-----------------------------------|--|---|
| The analyzer shows an error code. | May be a temporary fault. | Turn off the analyzer and turn it on again after 30 seconds. Take a new microcuvette and repeat the measurement. If the problem continues, see specific error code below. |
| E00 | No stable endpoint is found within the time range. 1. The cuvette is faulty. 2. The circuit board is out of order. | 1a. Check the expiry date for the microcuvettes. 1b. Take a new microcuvette and repeat the measurement. 2. The analyzer needs service. Contact your distributor. |
| E01–E05 | 1. Fault or dirt in the optronic or electronic units. | 1a. Turn off the analyzer and clean the optronic unit using the HemoCue Cleaner. 1b. The analyzer needs service. Contact your distributor. |
| E06 | 1. Unstable blank value. The analyzer might be cold. | 1. Turn off the analyzer and allow it to reach room temperature. If the problem continues, the analyzer needs service. Contact your distributor. |
| E07 | 1. The battery power is too low. | 1a. The batteries need to be replaced. Turn off the analyzer and replace the batteries, 4 type AA, 1.5 V. 1b. Use the mains power supply. |
| E08 | The absorbance is too high. 1. Light blocking item in the cuvette holder. | 1a. Check that the analyzer and microcuvettes are being used according to the HemoCue Hb 201+ operating manual and instructions for use. 1b. The analyzer needs service. Contact your distributor. |
| E09–E30 | 1. Fault or dirt in the optronic or electronic units. | 1a. Turn off the analyzer and clean the optronic unit using the HemoCue Cleaner. 1b. The analyzer needs service. Contact your distributor. |
| E35 | 1. The current measurement has been incorrectly stored and will remain as error code E35 until the entire memory has been deleted. | 1a. No action is normally required. 1b. If storage problems continue, the analyzer needs service. Contact your distributor. |

| Symptom | Explanation | Action |
|---|--|--|
| HHH | 1. Measured value exceeds 256 g/L (25.6 g/dL, 15.9 mmol/L). | |
| No characters on the display. | <ol style="list-style-type: none"> 1. The analyzer is not receiving power. 2. If on battery power, the batteries need to be replaced. 3. The display is out of order. | <ol style="list-style-type: none"> 1a. Check that the mains adapter is connected to the mains power supply. 1b. Check that the mains adapter is securely connected to the analyzer. 1c. Check that the cable is not damaged. 2. Turn off the analyzer and replace the batteries, 4 type AA, 1.5 V. 3. The analyzer needs service. Contact your distributor. |
| The display gives erroneous characters. | <ol style="list-style-type: none"> 1. The display is out of order. 2. The microprocessor is out of order. | <ol style="list-style-type: none"> 1. The analyzer needs service. Contact your distributor. 2. The analyzer needs service. Contact your distributor. |
| The display shows "▣". | <ol style="list-style-type: none"> 1. The batteries need to be replaced. 2. If on mains power, the mains adapter or the circuit board is out of order. | <ol style="list-style-type: none"> 1. Turn off the analyzer and replace the batteries, 4 type AA, 1.5 V. 2a. Check that the mains adapter is properly connected and working. 2b. The analyzer needs service. Contact your distributor. |
| The display does not switch from "Hb" and "Hb" to three flashing dashes and "🔴"(ready for measuring). | <ol style="list-style-type: none"> 1. The magnet in the cuvette holder may be missing. 2. The magnetic sensor is out of order. | <ol style="list-style-type: none"> 1. The analyzer needs service. Contact your distributor. 2. The analyzer needs service. Contact your distributor. |

| Symptom | Explanation | Action |
|---|---|---|
| <p>Measurements on control materials are out of range – either too HIGH or too LOW.</p> | <ol style="list-style-type: none"> 1. The microcuvettes are beyond their expiry date, damaged or have been improperly stored. 2. The optical eye of the microcuvette is contaminated. 3. The control has not been mixed properly and/or is not at room temperature. 4. Air bubbles in the microcuvette. 5. The optronic unit is dirty. 6. The control is not suitable for use with the HemoCue Hb 201+ system. 7. The calibration of the analyzer has been changed. 8. The controls are beyond their expiry dates or have been improperly stored. | <ol style="list-style-type: none"> 1. Check the expiry date and the storage conditions of the microcuvettes. 2. Remeasure the sample with a new microcuvette. 3. Make sure that the control is mixed properly and at room temperature. 4. Check the microcuvette for air bubbles. Remeasure the sample with a new microcuvette. 5. Clean the optronic unit using the HemoCue Cleaner. 6. Only use controls intended for the HemoCue Hb 201+ system recommended by HemoCue. 7. The analyzer needs service. Contact your distributor. 8. Check the expiry date and the storage conditions of the control. Take a new microcuvette and repeat the measurement. |
| <p>Measurements on patient samples are higher or lower than anticipated.</p> | <ol style="list-style-type: none"> 1. The microcuvettes are beyond their expiry date, damaged or have been improperly stored. 2. The optical eye of the microcuvette is contaminated. 3. Air bubbles in the microcuvette. 4. The optronic unit is dirty. 5. The calibration of the analyzer has been changed. | <ol style="list-style-type: none"> 1. Check the expiry date and the storage conditions of the microcuvettes. 2. Remeasure the sample with a new microcuvette. 3. Check the microcuvette for air bubbles. Remeasure the sample with a new microcuvette. 4. Clean the optronic unit using the HemoCue Cleaner. 5. The analyzer needs service. Contact your distributor. |

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

Specifications

General

The HemoCue Hb 201⁺ is a system used for the determination of the total amount of haemoglobin in whole blood. The system consists of a specially designed analyzer with specially designed microcuvettes containing dried reagents. The microcuvette serves as pipette, reaction vessel and as a measuring microcuvette. No dilution is required. The haemoglobin measurement takes place in the analyzer, which follows the progress of the reaction until the end point has been reached. The system is factory calibrated against the haemiglobincyanide (HiCN) method, the international reference method for the determination of the haemoglobin concentration in blood, and needs no further calibration.

Intended purpose

Quantitative determination of haemoglobin in capillary, venous or arterial blood. The HemoCue Hb 201⁺ system is used for quantitative determination of haemoglobin in whole blood using a specially designed analyzer, the HemoCue Hb 201⁺ Analyzer, and specially designed microcuvettes, the HemoCue Hb 201 Microcuvettes. HemoCue Hb 201 Microcuvettes are for *In Vitro* Diagnostic use only. The HemoCue Hb 201⁺ Analyzer is only to be used with HemoCue Hb 201 Microcuvettes.

IVD Medical Device Directive

The HemoCue Hb 201⁺ complies with the IVD Medical Device Directive 98/79/EC and carries the CE mark.

Theory

Sodium deoxycholate haemolyses the erythrocytes and haemoglobin is released. Sodium nitrite converts haemoglobin to methaemoglobin which, together with sodium azide, gives azidemethaemoglobin. The absorbance is measured at two wavelengths (570 nm and 880 nm) in order to compensate for turbidity in the sample.

Storage and handling

Operating temperature for the system is 15-30 °C (59-86 °F). The system should not be operated at high (i.e. >90 % non-condensing) humidity. Use the HemoCue Hb 201 Microcuvettes prior to their expiry date. The expiry date is printed on each package.

Storage for microcuvettes kept in a vial

The microcuvettes are to be stored at room temperature (15–30 °C, 59–86 °F) and in a dry place. The expiry date is printed on each vial. Once the seal is broken the microcuvettes are stable for three months. Always keep the vial properly closed.

Storage for individually packaged cuvettes

The microcuvettes are to be stored at room temperature (15–30 °C, 59–86 °F) and in a dry place. The expiry date is printed on each package.

HemoCue Hb 201⁺ Analyzer

The analyzer can be stored and transported at 0–50 °C (32–122 °F). Operating temperature is 15–30 °C (59–86 °F). Allow the analyzer to reach ambient temperature before use. The analyzer should not be stored at high (i.e. > 90 % non-condensing) humidity.

Specimen collection and preparation

Capillary, venous or arterial blood may be used. Appropriate anticoagulants (e.g. EDTA or heparin) may be used, preferably in solid form to avoid dilutional effects. Mix all specimen tubes thoroughly on a mechanical mixer for at least 2 minutes or invert the tube 8–10 times by hand. Hemoglobin remains unchanged for days, provided that the blood does not become infected. If the specimen has been stored in a refrigerator, it will be viscous and the blood should be allowed to warm up to room temperature before mixing².

In case of peripheral circulatory failure, hemoglobin measurements from capillary sampling may be misleading. Confirmation of an unacceptable or unexpected laboratory result is considered as a common, well established routine in the laboratory. Extra cellular fluid present in capillary samples due to individual physiology or the affects of unacceptable capillary sampling technique can cause false hemoglobin results.

One of several pre analytical factors associated with technique revolves around the size and action mode of the lancet and how it is used. This and other factors associated with capillary sampling technique can have an impact on the puncture, blood flow and ultimately on the results obtained. Additional factors unrelated to sampling technique, such as an individual's capillary circulation status, must also be recognized as potential causes of false results. Therefore when testing for hemoglobin using capillary blood as the sample, confirmation of an unacceptable or unexpected result is needed to rule out potential pre analytical factors as the cause. This is particularly helpful when a cut off value is used, i.e. when screening for nutritional status based on hemoglobin. Since the individual performing the test may not always be aware of all the potential pre analytical factors or even those that are directly or indirectly associated with sampling techniques, a second finger stick should be performed when wanting to rule out pre analytical factors as a cause of an unacceptable or unexpected result. Performance of the second finger stick by a different individual may aid in lessening the affect of the factors that could have been technique related.

Always handle blood specimens with care, as they might be infectious. Consult local environmental authorities for proper disposal. Always wear protective gloves when handling blood specimens. The microcuvette is for single use only.

Materials required

- HemoCue Hb 201⁺ Analyzer
- HemoCue Hb 201 Microcuvettes
- Lancet (for capillary samples)
- Pipette or other transfer device (for venous, arterial or control material samples)
- Lint-free wipe (non-fraying)
- Hydrophobic surface (for venous, arterial or control material samples)

Quality Control

The HemoCue Hb 201⁺ Analyzer has an internal quality control, the "self test". Every time the analyzer is turned on, it will automatically verify the performance of the optronic unit of the analyzer. This test is performed every second hour if the analyzer remains switched on.

Follow local guidelines regarding quality control procedures. If a quality control test is required by local or other regulations and therefore should be performed, only use controls recommended by HemoCue, see relevant package insert for more information.

Expected values²

Adult Males 130–170 g/L (13.0–17.0 g/dL, 8.1–10.5 mmol/L)

Adult Females 120–150 g/L (12.0–15.0 g/dL, 7.4–9.3 mmol/L)

Infants, after neonatal period 110–140 g/L (11.0–14.0 g/dL, 6.8–8.7 mmol/L)

Children, two years to teenage: gradual increase to adult normals.

Due to a wide range of conditions (dietary, geographical, etc) which affect normal values, it is recommended that each laboratory establish its own normal range.

Measuring range

0–256 g/L (0–25.6 g/dL, 0–15.9 mmol/L).

Results above 256 g/L (25.6 g/dL, 15.9 mmol/L) will be displayed as HHH.

Limitations

HemoCue Hb 201 Microcuvettes are for *In Vitro* Diagnostic use only. The HemoCue Hb 201+ Analyzer is only to be used together with HemoCue Hb 201 Microcuvettes. For further limitations of the procedure, see the HemoCue Hb 201 Microcuvettes instructions for use.

Technical Specifications

Dimensions: 85×160×43 mm

Weight: 350 g (batteries included)

4 batteries type AA, 1.5 V

Mains adapter: CE marked

Only use adapters as recommended in the table, last in Specifications after section *Manufacturer*.

Pollution degree: 2

Overvoltage category: II

The instrument is tested according to IEC 61010-1, Second edition: 2001 and EN 61010-1: 2001, IEC/EN 61010-2-101: 2002, IEC 61326-1: 2005 and EN 61326-1: 2006, IEC 61326-2-6, First edition: 2005 and EN 61326-2-6: 2006, IEC 60601-1-2, Third edition: 2007 and EN 60601-1-2: 2007 and complies with the IVD Medical Device Directive 98/79/EC. The instrument is made for continuous mode.

Essential Performance

The essential performance is quantitative determination of hemoglobin in capillary, venous or arterial whole blood.

Warranty

The analyzer carries a 24-month warranty from the day of receipt.

After the warranty period service/repair is carried out at fixed prices.

Any other use of the system than recommended by the manufacturer will void the warranty.

Service and Disposal

The analyzer should be cleaned as recommended under Maintenance prior to service or disposal. Consult local environmental authorities for proper disposal.

Spare parts – Accessories

The following spare parts and accessories are available:

Mains adapter

Battery Lid

PC programme

Cables for PC and printer connection

HemoCue Cleaners

Hard Carrying Case

HemoCue Lancets

HemoCue Hb 201+ Control Cuvette

Recommended separation distance between Portable and mobile RF communications equipment and HemoCue Hb 201⁺ Analyzer

The HemoCue systems are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of HemoCue systems can help prevent electromagnetic interference by maintaining a minimum distance between portable and RF communications equipment (transmitters) and HemoCue systems as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter (W) | Separation distance according to frequency of transmitter (m) | | |
|---|---|--------------------------------------|---------------------------------------|
| | 150 kHz to 80 MHz $d=1.2\sqrt{P}$ | 80 MHz to 800 MHz $d=1.2\sqrt{P}$ | 800 MHz to 2.5 GHz $d=2.3\sqrt{P}$ |
| 0.01 | 0.12 | 0.12 | 0.23 |
| 0.1 | 0.38 | 0.38 | 0.73 |
| 1 | 1.2 | 1.2 | 2.3 |
| 10 | 3.8 | 3.8 | 7.3 |
| 100 | 12 | 12 | 23 |

For transmitters rated at maximum output power not listed above, the recommended separation distances (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where (P) is the maximum output power rating of the transmitter in watts (W) according to transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – Electromagnetic immunity

The HemoCue systems are intended for use in the electromagnetic environment specified below. The customer or user of the HemoCue systems should assure that it is used in such an environment.

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment – guidance |
|--|--|---|---|
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±6 kV contact ±8 kV air | ±6 kV contact ±8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Electrical fast transient/ burst IEC 61000-4-4 | ±2 kV for power supply lines ±1 kV for input/output lines | 2 kV for power supply lines ±1 kV for input/output lines | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | ±1 kV differential mode | ±1 kV differential mode | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11 | <5 %U (>95 % dip in U) for a 0.5 cycle 40 %U (60 % dip in U) for 5 cycles 70 %U (30 % dip in U) for 25 cycles <5 %U (>95 % dip in U) for 5 seconds For explanation of U see NOTE 1 | <5 %U (>95 % dip in U) for a 0.5 cycle 40 %U (60 % dip in U) for 5 cycles 70 %U (30 % dip in U) for 25 cycles <5 %U (>95 % dip in U) for 5 seconds | Mains power quality should be that of a typical commercial or hospital environment. If the user of the HemoCue systems requires continued operation during power mains interruptions, it is recommended that the HemoCue systems be powered from an uninterruptible power supply or a battery. |

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment – guidance |
|--|--|--|--|
| <p data-bbox="302 256 430 303">Conducted RF IEC 61000-4-6</p> <p data-bbox="302 338 430 385">Radiated RF IEC 61000-4-3</p> | <p data-bbox="540 256 706 303">3 Vrms 150 kHz to 80 MHz</p> <p data-bbox="540 338 706 385">3V/m 80 MHz to 2.5 GHz</p> <p data-bbox="540 421 755 438">See NOTE 2 and NOTE 3</p> | <p data-bbox="777 256 846 274">3 Vrms</p> <p data-bbox="777 338 836 356">3 V/m</p> | <p data-bbox="950 256 1338 415">Portable and mobile RF communications equipment should be used no closer to any part of the HemoCue systems, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter</p> <p data-bbox="950 421 1230 438">Recommended separation distance</p> <p data-bbox="950 473 1024 491">$d=1.2\sqrt{P}$</p> <p data-bbox="950 509 1190 526">$d=1.2\sqrt{P}$ 80 MHz to 800 MHz</p> <p data-bbox="950 544 1190 562">$d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p data-bbox="950 597 1338 697">Where (P) is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and (d) is the recommended separation distance in meters (m).</p> <p data-bbox="950 732 1338 832">Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b).</p> <p data-bbox="950 867 1338 909">Interference may occur in the vicinity of equipment marked with the following symbol.</p> <div data-bbox="1117 938 1187 1000" style="text-align: center;">  </div> |

NOTE 1 U is the a.c. mains voltage prior to application of the test level.

NOTE 2 At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 3 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the HemoCue systems are used exceeds the applicable RF compliance level above, the HemoCue systems should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the systems.
- b) Over the frequency range 150 KHz to 80 Mhz, field strength should be less than 3 V/m.

Technical specifications (EMC-RF)**Use only cables with the following specification:**

USB shielded maximum 2 m

Serial shielded maximum 1.5 m

| Guidance and manufacturer's declaration – electromagnetic emissions | | |
|--|------------|---|
| The HemoCue systems are intended for use in the electromagnetic environment specified below. The customer or the user of the HemoCue systems should assure that it is used in such an environment. | | |
| Emission test | Compliance | Electromagnetic environment – guidance |
| RF emissions | Group 1 | The HemoCue systems uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment |
| RF emissions | Class B | The HemoCue systems are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Voltage fluctuations/flicker emissions IEC 61000-3-3 | Complies | |

NOTE 1 It is the manufacturer's responsibility to provide equipment electromagnetic compatibility information to the customer or user.

NOTE 2 It is the user's responsibility to ensure that a compatible electromagnetic environment for the equipment can be maintained in order that the device will perform as intended.

Symbols Used

Caution

CE mark

Class II equipment

Only valid within the European Community. Indicates separate collection for waste of electrical and electronic equipment.

Temperature limitation

Efficiency Level

Serial port

DC inlet

To maintain safety use only adapter marked HCA01

Consult instructions for use

Relative humidity, non-condensing

References

1. Makarem, A. In clinical Chemistry: Principles and Techniques, 2nd ed., Henry, R.J., Cannon, D.C., and Winkelman, J.W., Harper and Row, Hagerstown, M.D., 1974, pp. 1125–1147.
2. Dacie and Lewis, Practical Haematology, Ninth Edition, 2001.
3. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard NCCLS Document H15-A.
4. Evaluation of precision performance of clinical chemistry devices; approved Guideline NCCLS Document EP5-A.

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Literatur

1. Makarem, A. In *Clinical Chemistry: Principles and Techniques*, 2nd ed., Henry, R.J., Cannon, D.C., and Winkelman, J.W., Harper and Row, Hagerstown, M.D., 1974, pp. 1125–1147.
2. Dacie and Lewis, *Practical Haematology*, Ninth Edition, 2001.
3. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard NCCLS Document H15-A.
4. Evaluation of precision performance of clinical chemistry devices; approved Guideline NCCLS Document EP5-A.

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HemoCue AB, Box 1204, SE-262 23 Ängelholm, Sweden. Phone +46 77 570 02 10. Fax +46 77 570 02 12.