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
Content

Page

1. Important information	12
1.1 Warranty terms	12
1.2 Stretchers identifikation	12
1.3 User groups	13
1.4 Declaration of conformity	13
2. Information on the safety and protection of the product	13
3. Purpose and use of the scoop stretcher	13
3.1 Adjustment of the scoop stretcher's length	14
3.2 Opening of the scoop stretcher	15
3.3 Patient pick-up and closing of the scoop stretcher	15
4. Cleaning and Care	17
5. Maintenance	18
6. Storage	18
7. Service life	19
8. Optional accessories	20
9. Identification and Warranty Certificate	21

1. Important information

Before using the scoop stretcher, please read the operating manual and safety instructions carefully. Keep the operating manual in a safe place to avoid operating errors.

Should you require further information or encounter difficulties in its use or maintenance, please contact your local dealership or ultraMEDIC GmbH directly. Always specify the serial number to identify the scoop stretcher.

1.1 Warranty terms

The manufacturer grants a warranty period of 24 months on this product. The warranty only encompasses the costs of the spare parts.

The following points are not covered by the warranty: a) possibly omitted shipping or labour costs incurred on the purchaser's side; b) defects in the product due to normal wear and tear; c) damages caused by the user; d) defects or damages caused by non-compliance with the instructions in this document; e) any modifications or repairs to the product without prior approval by ultraMEDIC GmbH.

1.2 Stretchers identification

The stretcher is identified by an autoadhesive badge fixed on the main structure (Figure 1). The Identification badge must never be removed. In case of damaging please require a duplicate. Any form of warranty will decay if the Identification badge is removed from the stretcher.


Figure1 Identification badge

1.3 User groups

This product may only be used by qualified, expert personnel. A prerequisite is theoretical and practical experience in handling the scoop stretcher.

1.4. Declaration of conformity

The scoop stretchers conform to the DIN EN 1865-1 standard and to the Medical Device Regulation (EU) 2017/745. A declaration of conformity for these products can be acquired from the local dealership or the manufacturer.

2. Information on the safety and protection of the product

1. Before and after every use, check the components and patient restraining belts for damage and completeness.
2. Ensure that the patient restraining belts are correctly installed.
3. Before initial use, conduct a test run with an empty scoop stretcher to become aware of the handling of the product.
4. In the event of damage or suspected damage to the straps or the scoop stretcher, remove the product from operation immediately until any occurring irregularities have been remedied.
5. All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State where the user and/or the patient is established.

3. Purpose and use of the scoop stretcher

The ultraSCOOP STRETCHER (X-RAY) is a scoop stretcher for motionless repositioning, immobilisation and transport of injured persons who may have suffered a spinal fracture.

3.1 Adjustment of the scoop stretchers length


Figure 2

The stretcher's length can be adjusted to 4 different lengths depending on the patient's size by using the profile extension of the foot area (Figure 3a,b). The device can also be used to fold the scoop stretcher for storage.

To adjust the length of the stretcher, follow these steps:

1. Place the closed scoop stretcher in a horizontal position.
2. Open the right and left blocking device (Figure 4a,b,c,d).
3. Pull out the two foot profiles and adjust to the desired position.
4. After setting the desired length, close the previously opened blocking devices again. Check whether the blockings of the two profile pull-outs are audibly engaged in one of the marked positions after pulling them out. The stretcher's length can be adjusted to 4 different lengths. (Figure 3a und 3b using the example of SAN-0125)
5. Check that the foot profile is properly engaged. The foot profile must not be able to be pulled out.


Figure 3a


Figure 3b


Figure 4a
SAN-0125
locked


Figure 4b
SAN-0125
unlocked


Figure 4c
SAN-0288
locked


Figure 4d
SAN-0288
unlocked

3.2. Opening of the scoop stretcher

To place the stretcher under the patient without moving him, the stretcher must be opened. To do this, please open the locks on the head and/or foot area. (Figure 2). The locks allow operators to open the scoop stretcher into the desired position (in 2 parts or "A" formed).


Figure 5 Press to unlock
(by using the SAN-0125 as an example)

Disassembling the scoop stretcher into 2 parts:

1. Place the closed scoop stretcher in a horizontal position.
2. Hold the right and left sections at the same time and press the central lock (Figure 5). Remove both sections from each other until they are fully open.
3. Repeat point 2 to open the foot side (Figure 6a).

To open the scoop stretcher in "A" position, please unlock only the lock of the foot side (Figure 6b).


Figure 6a
Scoop stretcher open on both sides


Figure 6b
Scoop stretcher opened in "A" position

Safety belts

The equipment includes - for a safe and correct resting position of the patient - 3 quick-release safety belts with a storage bag.

3.3 Patient pick-up and closing of the scoop stretcher

1. To position the ultraSCOOP STRETCHER (X-RAY) under the patient without moving him, unlock the locks on the foot and head side (point 3.2).

There are various application techniques for the scoop stretcher. You can place the halves under the patient from the sides or only open the scoop stretcher at the bottom end, which then forms an "A". Carefully position the top end under the patient's head and close the bottom end (with a so-called "scooping" motion).

The type of application must be individually determined, depending on the nature of the injury.

Make sure that while the patient is being placed on the stretcher, their body parts are not pinched and their clothes or hair do not get caught in the scoop stretcher.

2. Adjust the two sections to the length and position next to the patient to be transported. If the patient's legs protrude beyond the edge of the scoop stretcher, check to see whether the patient's legs require additional immobilisation. If the injured person is shorter than 1.20 m, we recommend if by using of ultraSCOOP STRETCHER X-RAY additional using the straps SAN-0281 or SAN-0286.

3. Slide the two sections of the stretcher under the patient.

4. Bring the two sections together until the corresponding blocking devices are audibly locked (lock head and foot side).

5. Ensure that both snap locks on head and foot side of the ultraSCOOP STRETCHER (X-RAY) are locked and only then secure the patient with the straps.

6. Secure all 3 patient safety straps. The latches must always audibly lock into place. For safety reasons, subject the latches to a test by pulling at the appropriate point. Ensure that the injured party's breathing is not encumbered by the traction of the straps.

5. Check all buckles and straps before lifting and transporting the patient.

4. Cleaning and Care

Never use products with aggressive agents (bleach, ammonia, etc.) for cleaning, avoid the use of scratching materials to remove stains such as brushes, steel wool and knife blades.

For **disinfection**, you can use commercially available alcohol-based agents. Other solutions or cleaning chemicals may damage the material.

Safety belts: The belts that have gotten wet should only dry naturally. Avoid direct heat radiation from heat sources. These factors could damage the material.

Metallic parts: Simply clean with a sponge drenched with water and normal neuter soap, or, if necessary, with proper disinfectant. Rinse and dry carefully. After that apply a thin protective wax coating.

Plastic parts: Clean with a sponge drenched with water diluted with disinfectant and then rinse and dry completely.

Lubrication: Lubrication can be simply executed by applying in the indicated points of the figure beneath one or two drops of medium viscosity oil. (Figure 7).


Figure 7
Lubrication

5.

In order to ensure the safety of the user and the injured party during a rescue mission, we recommend that you have the ultraSCOOP STRETCHER (X-RAY) and straps visually inspected by an expert on an **annual** basis.

For German users, we recommend consideration to the following regulations and guidelines of MPBetreibV.

For users outside of Germany, a regular check according to country-specific requirements is required.

The following parts must be checked and the following defects must be excluded:

Base plate	deformation, tears, breaks, loosened screws, broken parts
Safety belts	Cutting and abrasion points, intact seams, tissue-damaging stains
Metal parts	Rust, oxidation, deformation
Plastic parts	Deformation, tears, breaks
other damages	

In the event of damage or suspected damages, remove the product from operation immediately until any occurring irregularities have been remedied.

6. Storage

Store the ultraSCOOP STRETCHER in a dry and dark place. Avoid direct sunlight and direct heat radiation from heat sources. These factors can damage the material.

Protect the blocking devices and moving parts from dust and dirt deposits to maintain their maximum effectiveness during use.

The scoop can be folded to reduce the required space and making storing easier. (Figure 8abc).

Figure 8a

Place the closed scoop stretcher in a horizontal position.


Figure 8b

Open the right and left locking devices and unthread the two feet profiles at the same time until the folding joints are completely got out

Figure 8c

Taking hold of the frame at feet extremity turning all this part until complete superposition with the bust zone profile

7. Service life

The service life of the ultraSCOOP STRETCHER (X-RAY) depends on the individual conditions of use. Depending on how frequently and intensely the scoop stretcher is used, wear and tear on the ultraSCOOP STRETCHER (X-RAY) is normal. In this regard, correct use is assumed. The material is subject to an ageing process even if handled carefully. The product should be discarded depending on its condition.

Service life of the patient safety belts is set at max. **10 years**.

For the service life of the optional accessories, please refer to the relevant operating instructions.

8. Optional accessories

SAN-0284	ultraHEAD-LOCK I	Head restraining set
SAN-0285	ultraHEAD-LOCK II	Head restraining set
SAN-0286	ultraSTRAPS Xtreme	Safety strap system with clip locks
SAN-0281	ultraSTRAPS	Safety strap system with Velcro fastener
SAN-0280-T	ultraSPINE BAG	Storage bag
SAN-0084-S	ultraBELT	Safety strap set
SAN-0260-1	X-BOON ONE	Vacuum mattress
SAN-0260-2	X-BOON TWO	Vacuum mattress
SAN-0278	ultraVACPUMP I	Vacuum pump

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