**Abbott**

Panbio™

# COVID-19 IgG/IgM RAPID TEST DEVICE

(FINGERSTICK WHOLE BLOOD/VENOUS WHOLE  
BLOOD/SERUM/PLASMA)

---

REF **ICO-T402/ICO-T40203**

IVD

**EN: Instructions for Use**

Page 1

*In vitro* diagnostic rapid test for the qualitative detection of IgG and IgM antibodies (Ab) to SARS-CoV-2.

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## **INTENDED USE**

The Panbio™ COVID-19 IgG/IgM Rapid Test Device (Fingerstick Whole Blood/Venous Whole Blood/Serum/Plasma) is an *in vitro* diagnostic rapid test for the qualitative detection of IgG and IgM antibodies to SARS-CoV-2 in human serum, plasma, venous and fingerstick whole blood. The Panbio™ COVID-19 IgG/IgM Rapid Test Device (Fingerstick Whole Blood/Venous Whole Blood/Serum/Plasma) is for professional use only and is intended to be used as an aid in the diagnosis of SARS-CoV-2 infection. The product may be used in any laboratory and non-laboratory environment that meets the requirements specified in the Instructions for Use and local regulation.

The test provides preliminary test results. Negative results will not preclude SARS-CoV-2 infection and they cannot be used as the sole basis for treatment or other management decision.

The test is not intended to be used as a donor screening test for SARS-CoV-2.

## **SUMMARY**

The novel coronaviruses belong to the  $\beta$  genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

Both molecular and serological tests are needed to definitively confirm a virus carrier. Molecular tests (NAT assays for SARS-CoV-2) played an important role in identifying SARS-CoV-2 infection. The antibodies detected by serology test indicate that a person has an immune response to SARS-CoV-2, whether symptoms developed from infection or the infection was asymptomatic. Antibody tests are important to confirm that an infection has occurred.

It has been shown that the severity of COVID-19 symptoms is correlated to IgM and IgG antibody levels. The increased IgG response is an attribute of the immune response to high viral loads. IgM levels are lower in more severely affected patients at a later stage of the disease whereas IgG levels are higher in the same group. Due to patient variability in the immune response to COVID-19, antibody testing may be more beneficial at >14 days after a positive diagnosis or suspected exposure.

The Panbio™ COVID-19 IgG/IgM Rapid Test Device (Fingerstick Whole Blood/Venous Whole Blood/Serum/Plasma) is suitable for use in the detection of a wide range of clinical presentations when seroconversion is reached. Results must be combined with clinical observations, patient history, and epidemiological information. Antibodies may not be detected in the first few days of infection. If acute infection is suspected, direct testing for SARS-CoV-2 is necessary.

## **PRINCIPLE**

The Panbio™ COVID-19 IgG/IgM Rapid Test Device (Fingerstick Whole Blood/Venous Whole Blood/Serum/Plasma) is an immunochromatographic test for the qualitative detection of IgG and IgM antibodies to SARS-CoV-2. The test device may be used for near-patient testing.

The test contains a membrane strip and a plastic housing. The test device shows the letters C, G and M at the right side of the reading window and the letter S above the Specimen well of device where C stands for Control line, G for IgG test line, M for IgM test line and S for Specimen well, respectively. To use the test, the whole blood/serum/plasma is applied into the Specimen well (S) first and then two drops of buffer are applied. The mixture of specimen and buffer migrates along the membrane strip to the reading window. On the nitrocellulose membrane within the reading window anti-human IgG and anti-human IgM antibodies are precoated at the G area and M area separately and a goat anti-rabbit antibody is precoated at the C area. If the specimen is SARS-CoV-2 IgG antibodies positive, the G line will become visible. If the specimen is SARS-CoV-2 IgM antibodies positive, the M line will become visible. If the specimen is SARS-CoV-2 antibodies negative, only the C line will become visible. The Control line (C) must always be visible if the test has been performed correctly.

The visible Control line indicates that the result is valid. If the Control line does not

appear, the test result is invalid. When the Control line, IgG test line and/or IgM test line are visible this indicates a positive result. When only the Control line is visible this indicates a negative result.

## MATERIALS

### Materials Provided

COMPONENTS	ICO-T402	ICO-T40203
1. Test devices individually foil pouched with a desiccant	x25	x25
2. Buffer (3mL/vial)	x1	x1
3. Specimen dropper (for fingerstick whole blood only)	x25	x25
4. Lancet (for fingerstick whole blood only)	-	x25
5. Plaster	-	x25
6. Alcohol pad	-	x26
7. Instructions for use	x1	x1
8. Quick reference guide (for fingerstick whole blood)	-	x1

### Materials Required but Not Provided

- Lancet (for fingerstick whole blood only) [Applies to ICO-T402 only]
- Specimen collection equipment and containers
- Micropipette
- Cotton wool or gauze pad (for fingerstick whole blood only)
- Centrifuge
- Timer
- Biohazard waste containers for sharps and non-sharps

## PRECAUTION/KIT STORAGE AND STABILITY

- The test kit should be stored at 2-30°C (storage in refrigerator is permitted). Higher temperatures up to 55°C for a limited period of time (i.e. up to 48 accumulated hours) have no impact on test performance. Do not store the kit in the freezer.

- The test device and buffer are stable until the expiration date printed on the outer package. Do not use it beyond the expiration date.
- The test device must remain in the sealed pouch until use.
- Do not use it if the pouch is damaged or the seal is broken.
- The test device is recommended to be used at room temperature (15-30°C).
- Perform the test as soon as possible after removing the test device from its foil pouch (within one hour).
- After opening the pouch, the test is sensitive to relative humidity above 70%.
- Do not re-use the test device.
- EDTA, sodium citrate and heparin have been validated and can be used as anticoagulants. Other anticoagulants have not been validated. Do not use specimen treated with other anticoagulants.
- Do not mix buffers from different lots of test kits.
- Do not touch the Specimen well or reading window of the test device directly with the finger as this can cause incorrect results.

## **WARNINGS**

- For *in vitro* diagnostic professional use only.
- Do not eat, drink or smoke in the area where handling specimens or test kits.
- Read the instructions carefully before performing the test.
- Wear PPE such as gloves, laboratory coats, and eye protection when handling specimens. Follow proper infection control guidelines for handling all blood specimens and related items.
- Avoid splashing or aerosol formation and clean up spills thoroughly using an appropriate disinfectant, such as alcohol, isopropanol, or sodium hypochlorite.
- Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials, as if they were infectious waste, in a biohazard container.
- The buffer contains 0.09% sodium azide as a preservative which may be toxic if ingested. When eye contamination happens, rinse with running water for at least 15 minutes.

## SPECIMEN COLLECTION, STORAGE AND PREPARATION

### FINGERSTICK WHOLE BLOOD

- Clean entire fingertip (3rd or 4th finger from non-dominant hand) with an alcohol pad or an alternative cleaning material. Allow fingertip to dry.
- Puncture the side of the finger with a lancet.
- Wipe off the first drop of blood with a dry cloth or gauze.
- Squeeze the middle part of the specimen dropper and bring the open tip to the blood. Slowly release the pressure to draw the blood into the specimen dropper. Avoid air bubbles. Stop filling when the blood flows up to the marked black line (20  $\mu$ L).
- After collecting the finger blood specimen, apply plaster or place a gauze pad or cotton ball on the finger until the bleeding stops.
- Whole blood collected by fingerstick should be tested immediately but no longer than 2 minutes after collection.

### VENIPUNCTURE WHOLE BLOOD

- Using venipuncture, collect whole blood into an EDTA, sodium citrate or heparin collection tube and mix well.
- If the blood specimen is not immediately tested, it should be refrigerated at 2-8°C.
- If stored at 2-8°C, the blood specimen should be tested within 2 days.
- Do not freeze whole blood specimens.

### SERUM OR PLASMA

- Serum: Collect the whole blood into the collection tube (NOT containing anticoagulants) by venipuncture, leave to settle for 30 minutes to allow for blood coagulation and then centrifuge the blood to get a serum specimen.
- Plasma: Collect the whole blood into an EDTA, sodium citrate or heparin collection tube by venipuncture and then centrifuge the blood to get a plasma specimen.
- The centrifuge setting 1,000-1,500 x g for approximately 5 minutes is required and refrigeration is not required.
- If serum or plasma specimens are not tested immediately, they should be refrigerated at 2-8°C for up to 3 days. If testing is delayed more than 3 days, the specimen should be frozen at -20°C $\pm$ 10°C. Three (3) freeze/thaw cycles show no impact on test result.

## TEST PROCEDURE

**Allow the test device, buffer and specimen to reach room temperature (15-30°C) prior to testing.**

1. Remove the test device from the foil pouch and use it as soon as possible (within one hour).
2. Place the test device on a clean and flat surface. Write specimen ID on the test device.
3. Add specimen and buffer to the Specimen well (S) of the test device.

### For Fingerstick Whole blood

- Fill the specimen dropper and transfer approximately 20  $\mu\text{L}$  of fingerstick whole blood specimen to the Specimen well (S) by touching the tip of the Specimen Dropper to the bottom of the Specimen well (S).
- Next, hold the buffer bottle vertically and add 2 drops of buffer (approximately 60  $\mu\text{L}$ ) to the Specimen well (S). Start the timer and wait for 10 minutes.


### For Venipuncture Whole blood

- A micropipette may be used to apply 20  $\mu\text{L}$  of whole blood to the Specimen well (S).
- Next, hold the buffer bottle vertically and add 2 drops of buffer (approximately 60  $\mu\text{L}$ ) to the Specimen well (S). Start the timer and wait for 10 minutes.

### For Serum or Plasma

- A micropipette may be used to apply 10  $\mu\text{L}$  of serum or plasma to the Specimen well (S).
- Next, hold the buffer bottle vertically and add 2 drops of buffer (approximately 60  $\mu\text{L}$ ) to the Specimen well (S). Start the timer and wait for 10 minutes.

Please refer to the illustration below.

4. As the test begins to work, a red color moves across the reading window of the test device.
5. Wait for the red line/lines to appear. Read test results at 10 minutes. Do not read results after 20 minutes.

## INTERPRETATION OF RESULTS


A red line will appear at the C area of the reading window to show that the test is working properly. This line is the Control line.

A red line that appears at the G area of the reading window is the IgG test line.

A red line that appears at the M area of the reading window is the IgM test line.

\*NOTE: The intensity of the color in the IgG and IgM test lines may vary depending on the concentration of COVID-19 antibodies present in the specimen. Therefore, any visible line at the G area or the M area should be considered positive.

Please refer to the illustration below.

<p>(+)</p> 	<p><b>IgG POSITIVE:</b> * The presence of both Control line and IgG test line within the reading window indicates an IgG positive result.</p>
<p>(+)</p> 	<p><b>IgM POSITIVE:</b> * The presence of both Control line and IgM test line within the reading window indicates an IgM positive result.</p>
<p>(+)</p> 	<p><b>IgG and IgM POSITIVE:</b> * The presence of Control line and both IgM and IgG test lines within the reading window indicates a positive result for both, IgG and IgM.</p>
	<p><b>NEGATIVE:</b> The presence of only the Control line and no IgM and IgG test line within the reading window indicates a negative result.</p>
	<p><b>INVALID:</b> If the Control line is not visible in the reading window the test result is invalid. If this occurs, it is recommended to read the IFU again before re-testing the specimen with a new test device.</p>

## QUALITY CONTROL

A Control line is visible within the reading window after the test is performed. The Control line is used in the test as a procedural control. A visible Control line confirms that the lateral flow of the test is successful but is not the confirmation that the specimen and buffer have been applied properly. Quality control specimens are not supplied in this kit; however, it is recommended that quality control specimens can be tested as a good laboratory practice.

## LIMITATIONS

1. The Panbio™ COVID-19 IgG/IgM Rapid Test Device (Fingerstick Whole Blood/Venous Whole Blood/Serum/Plasma) is for *in vitro* diagnostic use only. This test is used for the detection of antibodies to SARS-CoV-2 in human serum, plasma and fingerstick and venipuncture whole blood. Other body fluids or diluted specimens may not give accurate results and should not be used.
2. For venipuncture whole blood and plasma, EDTA, sodium citrate or heparin may be used as the anticoagulants. Other anticoagulants have not been validated and may give incorrect results.
3. As with all diagnostic tests, all results must be interpreted together with other clinical information available to the medical doctor.
4. The test detects SARS-CoV-2 antibodies that are formed as part of the body's immune response to the infection caused by the virus rather than detecting the virus itself.
5. A negative result does not eliminate the possibility of a SARS-CoV-2 infection, particularly in those who have been in contact with the virus less than 14 days before testing. Follow-up testing with a molecular diagnostic should be considered to rule out infection in these individuals.
6. The test will show negative results under the following conditions: The titer of the novel coronavirus antibodies in the sample is lower than the minimum detection limit of the test, or the novel coronavirus antibody has not appeared at the time of sample collection.
7. There is no positive correlation between the intensity of a red line at G and/or M area and the titer of antibody in the specimen.
8. Positive results may be due to past or present infection with non-SARS-CoV-2 coronavirus strains, such as coronavirus HKU1, NL63, OC43, or 229E.
9. Reading test results earlier than 10 minutes or later than 20 minutes may give incorrect results.

## PERFORMANCE CHARACTERISTICS

### 1. Diagnostic Sensitivity and Specificity

#### **Plasma**

Clinical performance of the Panbio™ COVID-19 IgG/IgM Rapid Test Device (Fingerstick Whole Blood/Venous Whole Blood/Serum/Plasma) was determined by testing 89 positive frozen EDTA and heparin plasma specimens (confirmed positive with a leading commercial PCR test) and 250 negative frozen EDTA plasma specimens collected prior to September 2019.

Overall study results were as follows:

Sensitivity of 97.8% (87/89) with 95% confidence interval of [92.1%-99.7%],

Specificity of 92.8% (232/250) with 95% confidence interval of [88.9%-95.7%],

Percent agreement of 94.1% (319/339) with 95% confidence interval of [91.0%-96.4%].

#### **Panbio™ COVID-19 IgG/IgM Rapid Test Device Performance - Overall**

		<b>PCR Positive</b>	<b>Clinical Negative</b>	<b>Total</b>
Panbio™ COVID-19 IgG/IgM Rapid Test Device Result	IgG or IgM Positive	87	18	105
	IgG and IgM Negative	2	232	234
	Total	89	250	339

The results for IgG demonstrate a specificity of 240/250 and a sensitivity of 86/89.

The results for IgM demonstrate a specificity of 240/250 and sensitivity of 46/89.

**Panbio™ COVID-19 IgG/IgM Rapid Test Device Performance – Elapsed Time of Sample Collection Post-Positive PCR Result**

		Elapsed Time of Sample Collection Post Positive PCR Result		
		5 Days	10 Days	Unknown
Panbio™ COVID-19 IgG/IgM Rapid Test Device Result	IgG or IgM Positive	47	36	4
	IgG and IgM Negative	1	1	0
	Total	48	37	4
	Sensitivity [95% CI]	97.9% [88.9%-99.9%]	97.3% [85.8%-99.9%]	100.0% [39.8%-100.0%]

***Fingerstick/Venous Whole Blood***

The clinical performance of the Panbio™ COVID-19 IgG/IgM Rapid Test Device (Fingerstick Whole Blood/Venous Whole Blood/Serum/Plasma) was evaluated with 26 fingerstick whole blood and 25 venous blood specimens that were confirmed to be SARS-CoV-2 positive with leading commercial PCR methods and the Abbott ARCHITECT SARS-CoV IgG test. In addition, 24 fingerstick whole blood specimens and 24 venous whole blood specimens were collected from presumed COVID-19 negative individuals as confirmed with the Abbott ARCHITECT SARS-CoV-2 IgG test and/or PCR and subjected to testing on the device. Fingerstick specimens were tested by non-laboratory health care personnel and venous whole blood specimens were tested by laboratory trained operators.

Overall study results are shown in the table below.

**Panbio™ COVID-19 IgG/IgM Rapid Test Device Performance – Fingerstick and Venous Whole Blood**

	Panbio™ COVID-19 IgG/IgM Rapid Test Device			
	Fingerstick Whole Blood		Venous Whole Blood	
Positive	25/26		24/25	
Presumed negative	24/24		23/24	
	Sensitivity	Specificity	Sensitivity	Specificity
	96.2% [80.4%-99.9%]	100.0% [85.8%-100.0%]	96.0% [79.7%-99.9%]	95.8% [78.9%-99.9%]

Matrix equivalency was also demonstrated in separate studies using serum, EDTA plasma, EDTA venous whole blood and fingerstick whole blood specimen from 35 healthy donors and heparin/ sodium citrate plasma and heparin/ sodium citrate venous whole blood from 10 healthy donors. Negative specimens were tested and in addition, spiked positive specimens were prepared by dilution of anti-SARS-CoV-2 IgG and IgM strong positive control specimens and tested on the Panbio™ COVID-19 IgG/IgM Rapid Test Device. The study demonstrated consistent results across all specimen types.

## 2. Precision

### ***Intra-Assay (same lot)***

Within-run precision has been determined by testing five specimens in 3 replicates: a negative serum, a negative plasma, a negative venipuncture whole blood, a negative fingerstick whole blood, and a SARS-CoV-2 low positive plasma. All above values were correctly identified >99% of the time.

### ***Inter-Assay***

Between-run precision has been determined by testing the same five specimens in 3 replicates: a negative serum, a negative plasma, a negative venipuncture whole blood, a negative fingerstick whole blood, and a SARS-CoV-2 low positive plasma. Three different

lots of the Panbio™ COVID-19 IgG/IgM Rapid Test Device (Fingerstick Whole Blood/Venous Whole Blood/Serum/Plasma) have been tested using above specimens. The specimens were correctly identified >99% of the time.

### **3. Interfering studies**

To ensure that other medical conditions (potentially interfering substances) did not affect the performance of the Panbio™ COVID-19 IgG/IgM Rapid Test Device (Fingerstick Whole Blood/Venous Whole Blood/Serum/Plasma), samples of SARS-CoV-2 negative blood were tested from people who had other conditions. These included (brackets show number of correct results/number of samples): HAMA (3/3); Antinuclear antibody (3/3); Rheumatoid factor (3/3); Influenza A and B IgG (3/3).

The Panbio™ COVID-19 IgG/IgM Rapid Test Device (Fingerstick Whole Blood/Venous Whole Blood/Serum/Plasma) was also tested with 6 inhibition interference substances. These included (tested concentration shown in brackets): Ascorbic acid (0.2mg/mL); Hemoglobin (10mg/mL); Bilirubin (0.6mg/mL); Oxalic acid (0.6mg/mL); Human Serum Albumin (20mg/mL); Triglyceride (50mg/mL). These substances were spiked with SARS-CoV-2 negative and SARS-CoV-2 low positive plasma and the test results indicated these interference substances gave no impact on the performance of the Panbio™ COVID-19 IgG/IgM Rapid Test Device (Fingerstick Whole Blood/Venous Whole Blood/Serum/Plasma).

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© 2020 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.
**Abbott Rapid Diagnostics Jena GmbH**