REF 02FK10CE, 02FK16CE, 02FK17CE

Bioline™
HCV

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About the test

[Introduction] Hepatitis C virus (HCV) is recognized as a major agent of chronic hepatitis, transfusion acquired non-A, non-B hepatitis and liver disease throughout the world. HCV is an enveloped positive-sense, single-stranded RNA virus. Testing for HCV infection begins serology testing with either a rapid or a laboratory-conducted assay for HCV antibody in blood. A reactive result indicates presumptive HCV infection. When to confirm the current HCV infection it is recommended nucleic acid testing (NAT) for the detection of HCV RNA be performed following HCV antibody reactive test result. Abbott Diagnostics Korea has constructed HCV genes for the expression of recombinant antigens in bacterium systems such as *E. coli* focused on structural and non-structural immunogenic regions of the HCV-encoded polyprotein. The major immunoreactive antigens of these proteins have been reported as core, NS3, NS4 and NS5 regions of HCV genome, which are known to be highly immunodominant. For detection of HCV infection, these recombinant proteins were used as capture materials in this immunochromatographic test. Improving on the first generation HCV serologic antibody test using a single recombinant antigen, recombinant proteins with multiple antigens have been used to minimize non-specific cross-reactivity and to increase the sensitivity in this assay.

[Test principle] The Bioline™ HCV test contains a nitrocellulose membrane strip, which is pre-coated with recombinant HCV capture antigen (core, NS3, NS4 and NS5) at the test line region (T). The protein A-colloid gold

conjugate and the specimen moves along the membrane chromatographically to the test region. There the antigen-antibody protein A gold particle complex forms into a visible line with high degree of sensitivity and specificity. This test device has letter “T” and “C” representing “Test Line” and “Control Line” on the surface of the case. Both the test line and control line in result window are not visible before applying the specimen. The control line is a procedural control. The control line should always appear if the test procedure is performed properly and the reagents in the control line are working.

[Intended use] The Bioline™ HCV test is an *in vitro* immunochromatographic, rapid assay designed for the qualitative detection of antibodies specific to HCV, in human serum, plasma, venous whole blood or finger prick whole blood. The Bioline™ HCV test kit is intended only for professional use as the initial test, as an aid to diagnosis. Reactive specimens should be followed up with nucleic acid testing (NAT) for the detection of HCV RNA together with a HCV immunoblot testing to identify current HCV infection. This test may not be suitable for diagnosis of early infection or blood donation screening.

Materials provided and active ingredients of main components

1. The Bioline™ HCV test kit contains the following items to perform the assay:

Catalog No.	Contents
02FK10CE	1. 30 Test devices with desiccant in individual foil pouches 2. Assay diluent (1 x 5 ml/vial) 3. 1 Instructions for use
02FK16CE	1. 25 Test devices with desiccant in individual foil pouches 2. Assay diluent (1 x 5 ml/vial) 3. 25 Capillary pipettes (10 µl), 25 Sterile lancets, 25 Alcohol swabs 4. 1 Instructions for use
02FK17CE	1. 25 Test devices with desiccant in individual foil pouches 2. Assay diluent (1 x 5 ml/vial) 3. 25 Capillary pipettes (10 µl), 25 Safety lancets, 25 Alcohol swabs 4. 1 Instructions for use

2. Active ingredients of main components
 - 1 Test device includes:
 - Gold conjugates: Protein A – gold colloid ($1.0 \pm 0.2 \mu\text{g}$)
 - Test line: Recombinant HCV antigen (core, NS3, NS4, NS5) ($1.5 \pm 0.3 \mu\text{g}$)
 - Control line: Goat anti-human Immunoglobulin ($2.0 \pm 0.4 \mu\text{g}$)

- Assay diluent includes: 50 mM Tris-HCl Buffer (5 ml), Sodium azide (0.02 %)

Materials required but not provided

- Micropipette, Protective gloves, Timer, Biohazard container

Kit storage and stability

1. The test kit should be stored at a temperature between 1 °C and 30 °C. Do not freeze the kit or its components.
2. Assay diluent may be opened and resealed for each assay. Cap should be firmly sealed between each use. Assay diluent is stable until expiration date if kept at 1 - 30 °C.
3. The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from foil pouch.
4. Do not use the test kit beyond its expiration date. The shelf life of the kit is as indicated on the outer package.
5. Do not use the test kit if the pouch is damaged or the seal is broken.

Warnings

1. The test devices are for *in vitro* diagnostic use only. Do not reuse the test device.
2. The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be proficient.
3. Do not use the pipette by mouth, smoke, drink, eat, apply cosmetics, or handle contact lenses in areas where specimens or kit components are being handled.
4. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
5. Clean up spills thoroughly using an appropriate disinfectant.
6. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials in a biohazard container as if they were infectious waste.
7. Do not mix or interchange different specimens.
8. Do not eat the desiccant from the foil pouch.
9. Avoid splashing or aerosol formation of specimen and assay diluent.
10. Do not mix or interchange components among different lots or those for other products.
11. Do not drink assay diluent.
12. Care should be taken to avoid contamination of the bottle nozzle when dropping assay diluent into the specimen well.
13. The assay diluent contains a proprietary anti-microbial

agent, sodium azide, which presents no hazard to the user if normal laboratory safety precautions are followed. If contact with assay diluent to the eyes and/or skin, wash affected area with soap and water immediately. If irritation or signs of toxicity occur, seek medical attention.

14. The assay diluent contains sodium azide, which may react with lead or copper plumbing to form highly explosive metal azide compounds. When disposing of these reagents through plumbing fixtures, flush with a large volume of water to prevent azide build-up in drains.
15. Even though the color of assay diluent is changed to faint yellow, it doesn't affect the performance and stability of the kit.

Specimen collection and storage

1. Whole blood

[Collection by venipuncture]

- Using venipuncture, draw the whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate).
- If the blood specimen is not immediately tested, it must be refrigerated at 2 - 8 °C.
- If stored at 2 - 8 °C, the blood specimen must be tested within 3 days of refrigeration.
- Do not use a blood specimen stored for more than 3 days; it can cause a nonspecific reaction.
- Bring blood specimens to room temperature (15 - 30 °C) prior to use.

[Collection using a lancet]

- Clean the area to be lanced with an alcohol swab.
- Squeeze the fingertip then prick the lateral side of the finger with a lancet provided. Wipe away the first blood drop. Then, safely dispose of the lancet immediately after.
- Immerse the open end of a new capillary pipette (10 µl) in the next blood drop and release the pressure to draw blood into the capillary pipette up to the fill line.

2. Plasma or serum

- [Plasma] Using venipuncture, draw the whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate) and then centrifuge blood to obtain the plasma specimen.
- [Serum] Using venipuncture, draw the whole blood into the collection tube (NOT containing anticoagulants) then leave for 30 minutes to allow blood coagulation to occur. Centrifuge the tube to generate a serum specimen.
- If plasma or serum specimens are not tested immediately, they must be refrigerated at 2 - 8 °C. For storage period longer than 2 weeks, freezing (below -20 °C) is required. Bring plasma or serum specimens to room temperature (15 - 30 °C) prior to use.
- Plasma or serum specimens containing a precipitate may yield inconsistent test results. Such specimens must be clarified by centrifugation prior to assaying.

3. Precautions

- Repeated frozen-thawed cycle for specimen should be avoided.

- Anticoagulants including heparin, EDTA and sodium citrate do not affect the test result. Use of other anticoagulants has not been validated. Their use may affect the test result.
- Use new pipette tips for each specimen in order to avoid cross-contamination of other specimens which could cause erroneous results.

Test procedure

1. Bring all kit components and specimens to reach a temperature between 15 °C and 30 °C prior to testing.
2. Remove the test device from foil pouch and place it on a flat, dry surface. Label the test device with a patient identifier.
3. **[Using a micropipette]**
Dispense 10 µl of plasma, serum or whole blood specimen into the specimen well marked "S".
Or,
[Using a capillary pipette]
Dispense 10 µl of drawn whole blood specimen into the specimen well marked "S"
4. Dispense 4 drops of assay diluent into the specimen well "S".
 **Caution:** If you do not hold the bottle vertically, it can lead to inaccurate results. Exactly, 4 drops should be added. Adding more than 4 drops may result in reddish color background or an invalid result.
5. As the test begins to work, you will see purple color move across the result window in the center of the test device.

- Interpret test results 5 - 20 minutes after adding assay diluent. Do not read after 20 minutes.

 **Caution:** If the test result is not legible after 5 minutes due to high background color, read again later but within 20 minutes of adding the diluent. Reading outside of this time frame (before 5 min or after 20 min) may result in false results.

Test interpretation

- A colored control line will appear at “C” in the result window to show that the test is working properly.
 - The “T” section of the result window indicates the test result.
 - Non-reactive result:** The presence of only the control line (C) within the result window indicates a non-reactive result.
 - Reactive result:** The presence of the test line (T) and the control line (C) within the result window, regardless of which line appears first, indicates a reactive result.
-  **Caution:** The presence of any test line, no matter how faint, the result is considered reactive.
- Invalid result:** If the control line (C) is not visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly or the test may have deteriorated. It is recommended that the specimen be retested using a new test device.

Test limitations

- A non-reactive result does not preclude the possibility of infection with HCV. Other clinically available tests

are required if questionable results are obtained. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

- Due to the inherent design of qualitative IVD tests, a faint or absent test line (false non-reactive) may occur in specimens containing high antibody densities; the prozone effect. In order to obtain a definitive result, all clinical and laboratory findings should be evaluated.

Internal quality control

The Bioline™ HCV test device has 2 pre-coated lines on the surface of the test: “T” (test line) and “C” (control line). Neither the test line nor the control line is visible in the result window before applying a specimen. The control line is used for procedural control and shows only that the diluent has been applied successfully and that the active ingredients of the main components on the strip are functional, but it is not an assurance that the specimen has been properly applied; it is not a reactive specimen control.

Performance characteristics

Bioline™ HCV test kit has been evaluated in 3 different sites as below. The results of individual laboratories may vary from these data because the results can be unique to the population it serves depending upon geographical, patient, dietary, environmental and other factors.

1. Study 1 : European performance evaluation according to the common technical specification (2009/886/EC)
- 1) Diagnostic sensitivity for anti-HCV detection
 - Anti-HCV positive specimens
410 anti-HCV positive specimens were tested by Sanquin in the Netherlands and the German Red Cross. The diagnostic sensitivity for anti-HCV antibody detection, calculated on 410 positive specimens, was 99.3 %.
Table 1. Result obtained with the Bioline™ HCV on all Anti-HCV positive specimens

	Bioline™ HCV	
	Reactive	Non-reactive
Anti-HCV positive serum/plasma (n=213)	211	2
Anti-HCV positive whole blood/plasma couples (n=100)	100	0
Anti-HCV with known genotype (n=97)	96	1
Total (n=410)	407	3
Sensitivity (95 % CI)	99.3 % (97.9 - 99.8 %)	

- Anti-HCV positive serum/plasma
213 anti-HCV positive specimens confirmed with CHIRON HCV RIBA™ 3.0 SIA or INNO-LIA HCV Score without genotype information were tested at Sanquin. In two specimens the Bioline™ HCV obtained a non-reactive result.
- Anti-HCV positive whole blood/plasma couples (Paired specimens)

100 paired whole blood and plasma specimens anti-HCV positive with Abbott Architect anti-HCV, taken from HCV infected patients were tested at the German Red Cross. The Bioline™ HCV test was reactive on all 100 paired specimens.

- Anti-HCV with known genotype
94 specimens from patients with a known HCV genotype 1 to 6, based on the VERSANT HCV Genotype 2.0 Assay (LiPA) were tested by Sanquin, and 3 specimens of genotype 5 were tested by the German Red Cross.

Table 2. Anti-HCV positive specimens with known genotype

HCV genotype		# Specimens	Bioline™ HCV	
			Reactive	Non-reactive
Genotype 1 (n=23)	1	1	1	0
	1a	10	10	0
	1b	12	11	1*
Genotype 2 (n=22)	2a/2c	13	13	0
	2b	9	9	0
Genotype 3 (n=25)	3	1	1	0
	3a	22	22	0
	3b	1	1	0
	3k	1	1	0
Genotype 4 (including non 4a) (n=20)	4a/4c/4d	4	4	0
	4c/4d	14	14	0
	4h	2	2	0
Genotype 5 (n=5)	5a	2	2	0
	5	3	3	0
Genotype 6 (n=2)	6	1	1	0
	6a	1	1	0

*One genotype 1b specimen was non-reactive on the Bioline™ HCV.

- Sensitivity on seroconversion panels
20 commercially available seroconversion panels (SeraCare Life Sciences and Zeptoemix) were tested with Bioline™ HCV by Sanquin.
In total 54 specimens were tested reactive with Bioline™ HCV and 47 with the competitor tests.

Panel ID	Bioline™ HCV	Ortho HCV 3.0 or Ortho Enhanced SAvE Anti- HCV 3.0 (*)
PHV904	3/7	3/7
PHV905	3/9	4/9
PHV911	3/5	3/5
PHV913	0/4	0/4
PHV914	5/9	3/9
PHV915	2/4	2/4
PHV918	1/8	2/8
PHV919 (*)	7/7	3/7
PHV920	6/10	6/10
6213	2/12	2/12
6214	6/13	4/13
6224	2/6	0/6
6226	4/12	3/12
6228	3/12	3/12
6229	3/8	3/8
9044	2/6	2/6
9045	2/8	1/8
9047	4/10	4/10
9054	1/10	1/10 (*)
9058	2/5	1/5 (*)
Total reactive bleeds (*)	54	47

(³): The Bioline™ HCV test obtained faint reactive lines on all panel members of panel PHV919. Since the first bleeds are non-reactive with all the competitor tests, these initial faint reactions with Bioline™ HCV were considered as nonspecific and therefore this panel was not taken into account for the evaluation.

- Bioline™ HCV was evaluated using additional 12 commercial HCV seroconversion panels (SeraCare Life Sciences, ZeptoMetrix and Biomex GmbH) and showed comparable performance to 5 CE-marked current Anti-HCV screening tests.

No.	Panel	Bioline™ HCV	Bioline™ HCV	Test 1	Test 2	Test 3	Test 4	Test 5
	Cat. no. #	No. reactive	No. indeterminate	No. reactive				
1	PHV 906	7	0	7	n.a.	6	7	n.a.
2	PHV 908	10	0	10	n.a.	7	8	8
3	PHV 910	3	0	3	3	3	3	3
4	PHV 912	0	1	n.a.	1	1	1	1
5	PHV 916	4	0	3	4	0	2	3
6	PHV 921	10	0	9	10	n.a.	7	9
7	6211*	3	0	3	3	2	2	3
8	6212	8	0	8	8	6	7	8
9	6222	1	0	1	n.a.	n.a.	1	1
10	6227	2	0	2	n.a.	2	2	2
11	9055	0	0	1	1	1	0	0
12	RP006**	4	0	n.a.	n.a.	n.a.	6	n.a.

n.a.: information not available

*ZeptoMetrix 6211: only panel members 35 to 40 tested.

**Biomex RP006: panel members 8 and 9 not available.

2) Diagnostic specificity

In total 1,000 plasma specimens and 500 whole blood specimens from blood donors negative for anti-HCV with ABBOTT PRISM were tested. The specimens originated from 2 collection sites in Germany, Frankfurt and Ulm. The diagnostic specificity calculated on 1,500 negative specimens was 100 %.

Table 3. Diagnostic specificity on anti-HCV negative blood donors

	Bioline™ HCV	
	Reactive	Non-reactive
Blood donors plasma specimens	0	1,000
Blood donors whole blood specimens	0	500
Specificity (95 % CI)	100% (99.7 - 100 %)	

2. Study 2 : WHO laboratory evaluation

• WHO HCV specimen reference panels

The 483 specimens that clinically derived serum/plasma specimens from Australian, European, African, Latin American and Asian origin was tested on the Bioline™ HCV. There were 163 anti-HCV positive specimens and 320 anti-HCV negative specimens. Each specimen was confirmed by commercial EIA tests, Murex anti-HCV EIA version 4.0 (DiaSorin S.A. Italy) and Monolisa anti-HCV PLUS version 2.0 (Biorad Laboratories).

Reference assay		Bioline™ HCV		Total Results
		Reactive	Non-reactive	
EIA	Positive	163	0	163
	Negative	0	320	320
Sensitivity (95 % CI)		100 % (97.8 - 100 %)		
Specificity (95 % CI)		100 % (98.8 - 100 %)		

• HCV seroconversion panels

In four different seroconversion panels (PHV913, PHV919, PHV920 and PHV922), the Bioline™ HCV detected antibodies to HCV on average 3 specimens later than the reference assay (DiaSorin Anti-HCV Murex EIA, Version 4.0).

• HCV performance panels

The Bioline™ HCV assay detected 15/15 HCV antibody positive specimens in the HCV mix titre performance panel (0810-0175) and 8/10 HCV antibody positive specimens in the HCV low titre performance panel (0810-0192).

3. Study 3

Specimens were collected from 243 persons in Bangladesh with four (4) types of blood (serum, plasma, venous whole blood, and finger prick blood). Bioline™ HCV kit showed 100% (243/243) equivalent result to the 3 types of specimen corresponding serum, plasma, and venous whole blood. In case of finger prick blood, there was only one discordance specimen, it was given 99.6% (242/243) agreement with other specimen types. This specimen was confirmed as weak positive by a commercial anti-HCV assay.

Specimen type		Finger prick blood		
		Positive	Negative	Total
Serum, plasma, venous whole blood	Positive	85	1*	86
	Negative	0	157	157
	Total	85	158	243
Agreement (95 % CI)		99.6 % (97.7 - 99.9 %)		

*The specimen was weak positive of Hepatitis C confirmed by commercial anti-HCV assay.

4. Analytical performance

1) Analytical specificity

- 200/200 hospitalized patients were non-reactive on the Bioline™ HCV.
- 200/200 specimens of pregnant women (incl. 20 specimens of multipara) were non-reactive on the Bioline™ HCV.
- 100/100 specimens containing potentially cross reactive substances were non-reactive on the Bioline™ HCV (20 Anti-HBs positive, 20 Anti-HBc positive, 20 Anti-HIV positive, 10 Anti-HTLV I/II positive, 20 Anti-HEV positive, 10 Rheumatoid factor positive).
- The following potential cross-reacting pathogens had no effect on test results of Bioline™ HCV.

Anti-HBs	CMV	HIV	Syphilis
<i>Borrelia burgdorferi</i>	EBV	HTLV	Toxoplasma
Chlamydia	HBsAg	Influenza	<i>Trypanosoma cruzi</i> I / II

- Bioline™ HCV was tested with potential interfering substances present in the patient's specimen and showed no interference by pregnant women, high cholesterol (up to 327 mg/dL), high bilirubin (up to 2.8 mg/dL), rheumatoid factor (up to 189.4 IU/ml), lipemic, hemolyzed, autoimmune, alcoholic cirrhosis, multiparous pregnancy.
- Pharmaceutical substances
There was no significant interference with the following 25 drugs on Bioline™ HCV. All drugs were tested at concentrations of 250 µg/mL.

Abacavir	Ergocalciferol	Naproxen
Acetaminophen	Folic acid	Nevirapine
Acetylsalicylic acid	Hydrochlorothiazide	Pantoprazole
Amoxicillin	Ibuprofen	Pyrazinamide
Aspirin	Iron Chloride	Rifampicin
Cholecalciferol	Isoniazid	Ritonavir
Cyclobenzaprine	L-ascorbic acid	Salicylic acid
Darunavir	Magnesium sulfate	
Diclofenac	Metformin	

2) Prozone effect

Bioline™ HCV may exhibit prozone effect (false non-reactive result) in specimens which have higher than S/CO of approx. 11.0 in Abbott Architect and S/CO approx. 3.0 in Ortho HCV version 3.0 ELISA anti HCV assays.

3) Specimen matrix

Validation on whole blood was done by testing 500 negative and 100 positive anti-HCV specimens. The performance of Bioline™ HCV on whole blood was comparable to the performance on plasma specimens. (See table 1 and table 3 presented above.) Validation on serum and different plasma specimen types (EDTA / Heparin / Sodium citrate) was performed by German Red Cross. The results obtained on negative and positive specimens are identical in serum, EDTA plasma, heparin plasma and citrate plasma.

Table 4. Serum / plasma (EDTA / Heparin / Sodium citrate) equivalence

Specimen type	No. of Bioline™ HCV reactive/No. of positive specimens	No. of Bioline™ HCV non-reactive/No. of negative specimens
Serum	25/25	25/25
EDTA plasma	25/25	25/25
Heparin plasma	25/25	25/25
Na-Citrate plasma	25/25	25/25

4) Complement factors interference in fresh serum specimens

In total 25 negative specimens, spiked with an anti-HCV positive specimen were tested within 24 hours after collection and retested after being stored at 4 °C for 1, 2, 3 and 4 days. No differences were observed on the results obtained on the fresh specimens and the same specimen stored for 1 to 4 days at 4 °C.

Table 5. Result obtained on spiked fresh specimens

		Bioline™ HCV	
		Reactive	Non-reactive
25 negative specimens	day 0	0	25
25 negative specimens spiked with anti-HCV	day 0	25	0
	day 1	25	0
	day 2	25	0
	day 3	25	0
	day 4	25	0

- 5) Reproducibility of the Bioline™ HCV has been demonstrated using in-house reference panels for within-run, between-run and batch-to-batch studies. All values were identical to reference panel acceptance criteria.

PREPARATION

1

Open the package and look for the following:

1. Test device with desiccant in individual foil pouch
2. Assay diluent
3. Instructions for use

Including only for Catalog No. 02FK16CE and 02FK17CE

1. Capillary pipette (10 µl)
2. Lancet
3. Alcohol swab

2

Carefully read the instructions on how to use the Bioline™ HCV test kit.

3


Look at the expiration date at the back of the foil pouch. If the expiration date has passed, use another kit. To avoid false results, ensure that the assay diluent used is from the same kit as the new test device.

4

Open the foil pouch and look for the following:

1. Test device
2. Desiccant

Then, label the device with the patient identifier.

● 10 µl : 10 µl of serum, plasma or whole blood specimen


□▶●●●● : Assay diluent 4 drops

TEST PROCEDURE

I. Blood (by venipuncture), Plasma or Serum specimen


- 1** Take 10 μ l of serum, plasma or whole blood specimen using a micropipette.

Dispense 10 μ l of serum, plasma or whole blood specimen into the specimen well "S".

- 2** Hold assay diluent bottle vertically and dispense 4 drops of assay diluent into the specimen well "S". Exactly, 4 drops should be added. Do not let bottle tip touch device in order to avoid cross-contamination.

4 x 


- 3** Interpret test results 5 - 20 minutes after adding assay diluent. Reading outside of this time frame (before 5 min or after 20 min) may provide false results.

TEST PROCEDURE


II. Blood specimen (with a lancet)

- 1** Clean the area to be lanced with an alcohol swab.


- 2** Squeeze the fingertip then prick the lateral side of the finger with a lancet provided. Wipe away the first blood drop. Then, safely dispose of the lancet immediately after.


[02FK16CE]


[02FK17CE]


- 3** Immerse the open end of a new capillary pipette (10 μ l) in the next blood drop and release the pressure to draw blood into the capillary pipette up to the fill line.

- 5** Hold assay diluent bottle vertically and dispense 4 drops of assay diluent into the specimen well marked "S". Exactly, 4 drops should be added. Do not let bottle tip touch device in order to avoid cross-contamination.

- 4** Dispense 10 μ l of drawn whole blood specimen in the specimen well marked "S". Lightly touch the capillary pipette to the specimen pad while dispensing.


- 6** Interpret test results 5 - 20 minutes after adding assay diluent. Reading outside of this time frame (before 5 min or after 20 min) may provide false results.

TEST INTERPRETATION

NON-REACTIVE


The presence of only the control line (C) within the result window indicates a non-reactive result.

REACTIVE


The presence of the test line (T) and the control line (C) within the result window, regardless of which line appears first, indicates a reactive result.

⚠ Caution: The presence of any test line, no matter how faint, the result is considered reactive.

INVALID

If the control line (C) is not visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly or the test may have deteriorated. It is recommended that the specimen be retested using a new test device.

Product Disclaimer:

While every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The subject of the diagnosis should consult a doctor for further confirmation of the test result.

Warning:

The manufacturers and distributors of this product shall not be liable for any direct, indirect, or consequential losses, liability, claims, costs or damages arising from or related to an incorrect reactive or non-reactive result using this product.

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1. Arash G., Czeslaw W., Chao Lin, Stephen M. Feinstone, and Charles M. Rice : Expression and Identification of Hepatitis C Virus Polyprotein Cleavage Products. *Journal of Virology*, March. 1993, p.1385-1395
2. Young Gyu Cho, Min Kyung Yi, Kyung Lib jang, Chang Min Kim and Young Chul Sung : Cloning and Overexpression of the Highly Immunogenic Region of HCV Genome from Korean Patients. *Mol. Cells*, Vol. 3, 4-7 - 416
3. S. Osborne, E. Ceconato, S. Griva, F. Garetto, R. Calogero, C. Rosa and F. Bonelli : Expression in *E. coli* and purification of a chimeric p22-NS3 recombinant antigen of Hepatitis C Virus(HCV). *Federation of European Biochemical Societies*, Volume 324, number 3, 253-257
4. A. Yoshikawa, K. Takahashi, S. Kishimoto : Serodiagnosis of hepatitis C virus infection by ELISA for antibodies against the putative core protein (p20C) expressed in *Escherichia coli*. *Journal of Immunological Methods*, 148 (1992) 143-150
5. Mitchell B.L. et al, Impact freeze-thaw cycles and storage time on plasma samples used in mass spectrometry based biomarker discovery projects. *Cancer Informatics*, 2005;1:98-104
6. Cuhadar S. et al, The effect of storage time and freeze-thaw cycles on the stability of serum samples. *Biochemia Medica*, 2013;23(1):70-77

 **Lancet**
Manufacturer
Shandong Lianfa Medical Plastic Products Co., Ltd.
No.1 Shuangshan Sanjian Road 250200
Zhangjiu City, Jinan, Shandong PEOPLE'S
REPUBLIC OF CHINA

 **Lancet**
Manufacturer
HTL-STREFA S.A.
ul. Adamówek 7, 95-035 Ozorków, Poland

 **Alcohol swab**
Manufacturer
FA INC.
442-48, Chengyeon-ro, Yeondong-
myeon, Sejong-si, 339-862 Korea

Authorized Representative
Shanghai International Trading Corp.
GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany

Authorized Representative
MT Promedt Consulting GmbH
Altenhofstrasse 80, 66386 St. Ingbert, Germany
Phone: +49 6894 581020

Authorized Representative
MT Promedt Consulting GmbH
Altenhofstrasse 80, 66386 St. Ingbert, Germany
Phone : +49 6894 581020, Fax : +49 6894 581021

Abbott Diagnostics Korea Inc.
65, Borahagal-ro, Giheung-gu, Yongin-si, Gyeonggi-do,
17099, Republic of Korea
Tel : 82-31-899-2800 | Fax : 82-31-899-2840
abbott.com/poct

Date issued : 2019. 11
02FK10CE/02FK16CE/02FK17CE-06-A0

© 2019 Abbott. All rights reserved.
All trademarks referenced are trademarks of either the Abbott group of
companies or their respective owners.