IVD **CE** German

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[Referenzen]

- Zhu N, Zhang D, Wang W, Li X, Yang B, Wu G, Gao GF, et al. A Novel Coronavirus from Patients with Pneumonia in China, 2019. *N Engl J Med*. 2020 Jan 24.
- Ji W, Wang W, Zhao X, Zai J, Li X. Homologous recombination within the spike glycoprotein of the newly identified coronavirus may boost cross-species transmission from snake to human. *J Med Virol*. 2020 Jan 22.
- Novel Coronavirus primer and Probe sequences, National Institute of Viral Disease, Chinese Center for Disease Control and Prevention, 2020.01.21
- Clinical management of severe acute respiratory infection when novel coronavirus (nCoV) infection is suspected—Interim guidance. 2020.

[Hersteller]

Name des Registrierenden/Herstellers: Guangzhou Pluslife Biotech Co., Ltd.
Adresse: Room 402, 6 Lianhuan Road, Huangpu District, Guangzhou, Guangdong, China
Postleitzahl: 510700
Kontakt: +86 20 84156831

[EU-Vertreter]

Name: Kingsmead Service B.V.
Adresse: Zonnehof 36, 2632 BE, Nootdorp, Netherlands
Tel.: +31 (0) 646571005
E-Mail: office@kingsmead-service.com

Dokument-Nr.: PB/CE-COV2-06 rev.A/0
Version: A/0
Datum: 28. Dez. 2021

IVD CE English

SARS-CoV-2 Nucleic Acid Testing Card FOR PROFESSIONAL USE ONLY

[Product Name]

Product Name: SARS-CoV-2 Nucleic Acid Testing Card

Trade name: Pluslife SARS-CoV-2 Card

[Intended Use]

This assay kit is used for in vitro qualitative detection of the N gene and ORF1ab gene of novel coronavirus SARS-CoV-2 in nasal swab samples from suspected pneumonia cases of SARS-CoV-2 infection and suspected cluster cases and other persons requiring diagnosis or differential diagnosis of SARS-CoV-2 infection.

Test results obtained with this kit are for clinical reference only and shall not be used as the sole standard for clinical diagnosis. A comprehensive analysis of the patient's clinical manifestations and other laboratory tests are recommended.

[Testing Principle]

This assay kit is based on isothermal amplification technology and enzyme digestion probe technology. The SARS-CoV-2 N gene and ORF1ab gene were selected as the amplification target region in designing the specific probe and primers. A large number of copies of the target sequence were generated in the reaction system during the isothermal amplification process. When the probe hybridizes to the complementary sequence, it is cleaved and fluorescence is emitted. The Integrated Nucleic Acid Testing Device detects and analyzes the fluorescence signal automatically, reporting a negative, positive or invalid result. The assay includes internal control for monitoring during sample collection, processing, and amplification, in order to reduce false negative results.

[Composition and Article Number]

Article No. and specifications	RM10102 02-1 1 Test	RM10102 02-2 2 Tests	RM10102 02-5 5 Tests	RM10102 02-10 10 Tests	RM10102 02-20 20 Tests	RM10102 02-50 50 Tests
Component name						
SARS-CoV-2 Reaction Card (piece)	1	2	5	10	20	50
Nucleic Acid Releasing Agent 01 (1 tube)	1	2	5	10	20	50
Disposable Sampling Swab (branch)	1	2	5	10	20	50
Waste Bag (piece)	1	2	5	10	20	50

Note:

1. The above components of different batches of kits must not be used interchangeably.
2. The card holder is in the packaging box of the Integrated Nucleic Acid Testing Device.

[Storage Conditions and Shelf Life]

1. 2–28°C storage, validity period 9 months.
2. The production date and expiration date are shown on the packaging label.

[Applicable Instruments]

Integrated Nucleic Acid Testing Device (PM001)

[Sample Requirements]

Nasal Swab

[Testing Method]

The room temperature should be between 15–28°C. Please read all the instructions carefully before you begin.

Be careful of the operations printed in red font. These operations are easily omitted. Omitting these operations will probably lead to an incorrect/invalid result.

STEP 1: Instrument Preparation

1. Take out the device, power adapter, and card holder from the Integrated Nucleic Acid Testing Device packaging.

2. Put the Integrated Nucleic Acid Testing Device on a flat surface, connect the power supply, press the button on the front of the Instrument to start the warm-up process (power light flashes red). After about 2 minutes, warm-up is completed and the device is in standby (power light is blue).

STEP 2: Sample Collection

1. Wash or disinfect your hands before performing the test procedure.
2. Bring the kit to room temperature (15 to 28°C) before use, then take out the components required for one test.


NOTE: Make sure there is no damage to the SARS-CoV-2 Reaction Card packaging, and no liquid leakage of the nucleic acid releasing agent 01. Do not use if any leakage has occurred.

NOTE: Before sample collection, maintain proper infection control and use recommended personal protective equipment (PPE), including an N95 or higher-level respirator (or face mask if a respirator is not available), eye protection, gloves, and a gown.

3. Have the test subject wipe away any excess mucus, if present.
4. Take out the disposable sampling swab from the notch and hold at the handle end.
5. Tilt the test subject's head back (approximately 70 degrees).
6. Gently insert the entire absorbent tip of the disposable sampling swab into one nostril until you meet resistance (approximately 1.5 to 2 cm).
7. Rotate it 3 times against the inner wall of one nostril, so that the absorbent surface of the swab is moistened.
8. Slowly remove the disposable sampling swab from the first nostril.
9. Repeat steps 6–8 on the other nostril using the same disposable sampling swab.

NOTE:


- 1) Avoid scratching the test subject.
- 2) Samples should be tested immediately after collection.

STEP 3: Sample Processing

1. Open the seal of the nucleic acid releasing agent 01 tube carefully to avoid spilling the liquid. Discard the seal into the waste bag.
2. Insert the used disposable sampling swab into the vial, making sure the absorbent tip is in the liquid. Then twist the disposable sampling swab tip against the bottom and sides of the nucleic acid releasing agent 01 tube 10 times and pinch the disposable sampling swab tip.
3. Discard the disposable sampling swab into the waste bag.
4. Screw the cap.

NOTE: Prevent the nucleic acid releasing agent 01 from contacting eyes or skin. In the event that this happens, wipe off the liquid immediately and rinse with plenty of water.

STEP 4: Sample Testing

Make sure the Integrated Nucleic Acid Testing Device is in standby (power light is blue). 1. Tear open the aluminum foil bag of one SARS-CoV-2 Reaction Card and take it out, place the SARS-CoV-2 Reaction Card on the card holder and unscrew the cap of the sample tube on the SARS-CoV-2 Reaction Card.

NOTE: Proceed to the following steps with the SARS-CoV-2 Reaction Card as soon as possible after the aluminum foil bag has been opened and the cap of the tube has been unscrewed.


2. Open the top cap of the nucleic acid releasing agent 01 tube from STEP 3, use one hand to stabilize the card holder, and use the other hand to pour the nucleic acid releasing agent 01 solution to the SARS-CoV-2 Reaction Card sample tube, between the two liquid injection lines, by squeezing the wall of the nucleic acid releasing agent 01 tube.

NOTE: There are two liquid injection lines marked on the SARS-CoV-2 Reaction Card sample tube. Add nucleic acid releasing agent 01 solution into the SARS-CoV-2 Reaction Card sample tube **until the liquid level is between the two lines**.

3. Place the nucleic acid releasing agent 01 tube in the waste bag for disposal.
4. Tightly screw down the cap of the SARS-CoV-2 Reaction Card sample tube.


5. Firmly press the protruding arc-shaped air bag on the sample tube of the SARS-CoV-2 Reaction Card to deform it and recess it into the tube.

NOTE: Wait for 10 seconds, DO NOT move on to the next step **until the nucleic acid releasing agent 01 fully fills all of the chambers** in the center of the card.


6. Hold the card and shake it up and down 10 times within about 5 seconds. The SARS-CoV-2 Reaction Card is ready to be tested.

NOTE: The SARS-CoV-2 Reaction Card must be tested immediately.

7. Open the top cover of the Integrated Nucleic Acid Testing Device, insert the SARS-CoV-2 Reaction Card into the device according to the direction indicated on the SARS-CoV-2 Reaction Card, pushing it to the fixed position of the bottom card slot, then close the top cover.

8. Press the button to start operation. The light flashes blue during operation.
9. Wait 15-35 minutes.
10. When the operation is completed, the indicator light on the device will display the test result. Record the result immediately. The assay is finished.

11. Open the top cover, **take out the SARS-CoV-2 Reaction Card immediately** and put it into the waste bag, seal the waste bag, and dispose of it following local regulations.

12. If proceeding with another test, return to STEP 2. If not, press the button for at least 3 seconds to turn off the device.

STEP 5: Interpretation of Test Results and Suggestions

The test results are determined as follows:

Phenomenon	Description	Result determination	Suggestions
	Positive indicator light on	The sample was determined to be SARS-CoV-2 positive.	In case of a positive result: a) Perform an RT-PCR confirmation test in the laboratory. b) Remain isolated until the result of confirmatory test. Contact a doctor if symptoms progress or in the event of any unusual symptoms.
	Negative indicator light on	The sample was determined to be SARS-CoV-2 negative.	In case of a negative result: a) If symptoms of COVID-19 appear, do a new test. b) If you have symptoms and are over 60, please contact your doctor to assess the need for further examinations.
	All lights on at the same time	Invalid result. Test should be repeated. Possible reasons: 1. The sample quantity is insufficient. The reaction is inhibited. 2. Operation error. 3. The sample is contaminated.	In case of invalid result: a) No conclusion can be made with this result. b) Perform a new test. c) If the problem persists, please contact the local distributor for assistance.

[Limitations of Detection Methods]

- Test results obtained with this kit are for clinical reference only. Comprehensive analysis and interpretation should be carried out based on the patient's symptoms/signs, medical history, and other laboratory diagnosis results. Test results obtained with this kit should not be used as the sole basis for clinical diagnosis, treatment or management of patients.
- If the amount of virus in the sample is insufficient, false negative results may occur.
- False positive results may occur if cross-contamination happens. Please dispose of the waste correctly following local regulations.
- Mutation of the target sequence during the virus epidemic or sequence changes caused by other reasons may lead to false negative results.

[Product Performance Indicators]

- Positive coincidence rate: the positive coincidence rate is 100% when testing the company's positive reference substances.

- Negative coincidence rate: the negative coincidence rate is 100% when testing the company's negative reference substances.
- Sensitivity (limit of detection): 400 copies/mL.
- Repeatability: the intra-assay repeatability detection rate reaches 100%; the inter-assay repeatability detection rate reaches 100%.
- Specificity: other pathogens that are similar to SARS-CoV-2 or cause similar symptoms have no cross-reactivity, such as coronavirus (HKU1, OC43, NL63, 229E), SARS coronavirus, MERS coronavirus, influenza A H1N1, H5N1, H7N9 and H9N2, influenza B virus, respiratory syncytial virus type A and B, human parainfluenza virus type II, adenovirus type 3 and 7, enterovirus EV71, Mycoplasma pneumoniae, EB virus, human cytomegalovirus and Mycobacterium tuberculosis.
- Exogenous interfering substances: therapeutic drugs for SARS-CoV-2 presented in the sample do not interfere with the test results of the kit, such as benflin, hydroxymezolin, sodium chloride, beclomethasone, dexamethasone, fluronium, triamcinolide, budesonide, mometasone, fluticasone, histamine, interferon, zanamivir, ribavirin, oseltamivir, peramivir, lopinavir, mupirocin, levofloxacin, azithromycin, tobramycin, ritonavir, meropenem, and ceftriaxone.
- Possible endogenous substances in the nasal swab sample, such as whole blood and mucin, do not interfere with the test results of the kit.

[Precautions]

- This kit is for in vitro diagnostic uses only, please read this instruction carefully before use, and use strictly in accordance with the instructions.
- Correct collection of swab samples and accurate operation according to the test method are critical to the accuracy of the test results.
- The test environment temperature should be between 15–28°C. If the kit is stored at a lower temperature, it must be returned to room temperature before opening, in order to avoid moisture condensation.
- When opening the top cap of the SARS-CoV-2 Reaction Card tube, do not touch the reaction tube and the inside of the tube cap.
- If the level of N gene and ORF1ab gene of novel coronavirus SARS-CoV-2 in the sample is below the limit of detection of the product, it may result in a false negative result.
- If the sample is tested more than 30 minutes after collection, it may result in a false negative result. Therefore, samples should be tested immediately after collection.
- If samples are not collected, transported or handled properly, it may lead to false negative results.
- The use of extracts other than those provided with this product can lead to incorrect results.

- A positive result only indicates the presence of N gene and ORF1ab gene of novel coronavirus SARS-CoV-2. Positive results do not exclude the possibility of co-infection with other pathogens.
- Negative results do not exclude SARS-CoV-2 infection and should not be used as the sole basis for decisions about treatment and patient management (including infection control measures).
- In the event of a negative result, the patient's recent exposure, medical history, and clinical signs and symptoms similar to those of COVID-19 should be taken into account.

- Patients whose symptoms persist for more than 7 days should be treated with caution if the test result is negative and, if necessary, confirmed by molecular testing for patient management.

- The length of the anterior nasal cavity may vary for users in different regions, 2 cm is for reference only. It is recommended that the user insert the swab until significant resistance is felt.
- All parts of this kit are for external use and should be kept out of reach of children.
- The test kit shall not be used after the expiry date indicated on the outer packaging.
- Exposure of the contents of the SARS-CoV-2 Reaction Card is likely to cause contamination. Do not damage the SARS-CoV-2 Reaction Card at any time.
- Disposal: all parts used have a potential risk of infection. Please use the provided waste bag for disposal.

- If collection and sample processing are not well controlled, it may result in cross-contamination, and false positives may occur.
- A variety of factors during the storage, transportation and use of reagents may cause changes in performance, such as sample collection, sample processing, and non-standard operation in the testing process. Please strictly follow the instructions. Due to the characteristics of the nasal swab and other sample collection process and the virus infection process itself, there may be false-negative results caused by insufficient sample size. They should be combined with other clinical diagnosis and treatment information to make a comprehensive judgment, and retests performed if necessary.

[Explanation of Symbols]

	CE Mark		Keep dry
	Consult instructions for use		Batch code
	Use by		In vitro diagnostic medical device
	Temperature limits		Date of manufacture
	Manufacturer		Do not use if package is damaged
	Authorized representative in the European Community		Keep away from sunlight
	Do not reuse		Sufficient for <n> tests
	Catalog number		

[References]

- Zhu N, Zhang D, Wang W, Li X, Yang B, Wu G, Gao GF, et al. A Novel Coronavirus from Patients with Pneumonia in China, 2019.N Engl J Med. 2020 Jan 24.
- Ji W, Wang W, Zhao X, Zai J, Li X. Homologous recombination within the spike glycoprotein of the newly identified coronavirus may boost cross-species transmission from snake to human. J Med Virol. 2020 Jan 22.
- Novel Coronavirus primer and Probe sequences, National Institute of Viral Disease, Chinese Center for Disease Control and Prevention, 2020.01.21
- Clinical management of severe acute respiratory infection when novel coronavirus (nCoV) infection is suspected— Interim guidance. 2020.

[Manufacturer]

Registrant/manufacturer name: Guangzhou Pluslife Biotech Co., Ltd.

Address: Room 402, 6 Lianhuanyuan Road, Huangpu District, Guangzhou, Guangdong, China

Zip: 510700

Contact: +86-20-84156831

[EU Representative]

Name: Kingsmead Service B.V.

Address: Zonnehof 36, 2632 BE, Nootdorp, Netherlands

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[Referencias]

1. Zhu N, Zhang D, Wang W, Li X, Yang B, Wu G, Gao GF, et al. A Novel Coronavirus from Patients with Pneumonia in China, 2019.N Engl J Med. 2020 Jan 24.
2. Ji W, Wang W, Zhao X, Zai J, Li X. Homologous recombination within the spike glycoprotein of the newly identified coronavirus may boost cross-species transmission from snake to human.J Med Virol. 2020 Jan 22.
3. Novel Coronavirus primer and Probe sequences, National Institute of Viral Disease, Chinese Center for Disease Control and Prevention, 2020.01.21
4. Clinical management of severe acute respiratory infection when novel coronavirus (nCoV) infection is suspected— Interim guidance.2020.

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[Riferimenti]

- [1] Zhu N, Zhang D, Wang W, Li X, Yang B, Wu G, Gao GF, et al. A Novel Coronavirus from Patients with Pneumonia in China, 2019.N Engl J Med. 2020 Jan 24.
- [2] Ji W, Wang W, Zhao X, Zai J, Li X. Homologous recombination within the spike glycoprotein of the newly identified coronavirus may boost cross-species transmission from snake to human. J Med Virol. 2020 Jan 22.
- [3] Novel Coronavirus primer and Probe sequences, National Institute of Viral Disease, Chinese Center for Disease Control and Prevention, 2020.01.21
- [4] Clinical management of severe acute respiratory infection when novel coronavirus (nCoV) infection is suspected— Interim guidance. 2020.

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[Referenties]

- Zhu N, Zhang D, Wang W, Li X, Yang B, Wu G, Gao GF, et al. A Novel Coronavirus from Patients with Pneumonia in China, 2019. *N Engl J Med*. 2020 Jan 24.
- Ji W, Wang W, Zhao X, Zai J, Li X. Homologous recombination within the spike glycoprotein of the newly identified coronavirus may boost cross-species transmission from snake to human. *J Med Virol*. 2020 Jan 22.
- Novel Coronavirus primer and Probe sequences, National Institute of Viral Disease, Chinese Center for Disease Control and Prevention, 2020.01.21
- Clinical management of severe acute respiratory infection when novel coronavirus (nCoV) infection is suspected— Interim guidance. 2020.

[Fabrikant]

Naam registrant/fabrikant: Guangzhou Pluslife Biotech Co., Ltd.

Adres: Room 402, 6 Lianhuayan Road, Huangpu District, Guangzhou, Guangdong, China
Postcode: 510700
Contact: +86-20-84156831

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[2]. Ji W, Wang W, Zhao X, Zai J, Li X. Homologous recombination within the spike glycoprotein of

the newly identifiedcoronavirus may boost cross-species transmission from snake to human.J Med Virol. 2020 Jan 22.

[3]. *Novel Coronavirus primer and Probe sequences, National Institute of Viral Disease, Chinese Center for Disease Control and Prevention, 2020.01.21*

[4]. *Clinical management of severe acute respiratory infection when novel coronavirus (nCoV) infection is suspected— Interim guidance.2020.*

[Fabricante]

Nome do fabricante registrado: Guangzhou Pluslife Biotech Co., Ltd.

Endereço: Room 402, 6 Lianhuayan Road, Huangpu District, Guangzhou, Guangdong, China

Código: 510700

Contacto: +86-20-84156831

[Representante UE]

Nome: Kingsmead Service B.V.

Endereço: Zonnehof 36, 2632 BE, Nootdorp, Netherland

Tel: +31(0)646571005

Email: office@kingsmead-service.com

Documento No.: PB/CE-COV2-06 rev.A/0

Versão: A/0

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