thermocult digital

D

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thermocult

Operating Instructions

Small-Size Incubator with permanently set incubation temperature

- permanently set incubation temperature of 35°C +/-2°C
- economical space configuration with horizontal or vertical mounting or placement. Suitable for 20 dip-slide cultures, 8 Petri- dishes or 6 blood culture bottles.
- 2 separate control lights for power and heating system
- built in digital-thermometer
- Shelf for tray cultures and staining racks/dishes (optional)
- All-around heating elements for uniform heat distribution
- Clear door-visual control, rapid check
- Radio interference suppression Degree N
- Casing of high quality poly-material
- Common cleaners and disinfectants can be used

Inhaltsverzeichnis	Seite
Control elements	9
Placement of the unit	9
Start up	10
Working temperature	10
Incubation temperature adjustment	10
Hints for the user	10
Cleaning the incubator	10
Troubleshooting	10
Change of the batteries (thermometer)	11
Technical specs.	12
Statement of conformity	13

Contents of the package upon delivery:

1 Instrument
4 foam rubberfeet
1 drilling template
1 batterie-changingtool
1 operating instruction incl.
CE-statement of conformity
1 Warranty-card

Shelf is not included

Control Elements

- 1 Master switch with Service control light (on top: ON)
- 2 Heating System control light
- 3 Digital-Thermometer
- 4 Cover for the temperature adjustment screw
- 5 Shelve (Accessories, please order separately)

1. Placement of the unit

The Small-size incubator thermocult can be used in a horizontal (Fig.1) or a vertical (Fig. 2) Position.

When setting the incubator upright, make sure to attach the foam rubber feet. The incubator **should not sit directly on the pad**

Fig. 1

Fig. 2

2. Start-Up

- a) Check to ensure that the voltage indicated on the rating plate and the Power supply voltage match.
- b) Plug the line plug into an electrical outlet (230V ~)
- c) Turn on the master switch (1): The red Service control light (1) and the yellow heating System control light (2) are on.
- d) The constant chamber temperature is reached after approx. **60 minutes**. The yellow heating System control light (2) will turn on and off at irregular intervals.

3. Working Temperature

The built-in temperature controller maintains a constant Incubation temperature of 35°C with a variation of $\pm 2^\circ\text{C}$ at various spots of the chamber (ambient temperature: 20°C - 25°C).

4. Incubation Temperature Adjustment

An adjustment is required if, after prolonged use of the Instrument, the permissible fluctuation span is no longer maintained. In order to proceed, remove the cover(4). If you turn the screw behind it to the right, you can adjust a higher incubation temperature. You can lower it by turning the screw to the left.

A one-quarter rotation corresponds to approx. 15°C.

5. Hints for the user

To close properly, the door must be pressed shut until impact. The thermocult unit must be installed in a Position where the door is not exposed to direct sun light.

The shelf should be used for all vessels or plates **so that they do not sit directly on the floor**.

It is recommended to switch off the device at longer disuse.

6. Cleaning the Incubator

Both the inner and the outer casing as well as the door of the incubator consist of high quality poly-materials. The surfaces will tolerate all commonly used cleaners and disinfectants.

Solvents (e.g. methanol, acetone, and the like) should not be used to clean the unit.

7. Trouble-shooting

Have the unit returned to the manufacturer by your dealer. If the Problem falls within the warranty, make sure to include the warranty certificate and add a short description of the fault.

Batteriewechsel / Change of Batteries

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ATTENTION: Pull out the mains plug first!

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Push change tool in the left recess beside the

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Thermometer lever up with the tool and take out

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
Reinsert flush the thermometer

Technical Specs.

Type	thermocult
Temperature constancy	$\pm 0.3^{\circ}\text{C}$
Incubation temperature	$35^{\circ}\text{C} \pm 2^{\circ}\text{C}$
Rated voltage	230V ~
Current and frequency	A.C. ~50/60 Hz
Rated current	0.14 A at 230 V
Rated capacity	30W
Class of protection	II (protective insulation) Not "ex" proof
Radio interference suppressed	degree N (VDE0875)
Technical safety compliance	
Capacity	20 dip-slide cultures, or 8 Petri dishes or 6 blood culture bottles
Dimensions:	
Casing	30 x 17 x 15cm (l x d x h)
Work chamber	21 x 14.5 x 11cm (l x d x h)
Weight	approx. 1.3kg
Accessories (optional)	shelves
Batteries:	2 x LR44/AG13 Button Cells 1,5 V

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