Spengler

AUTOTENSIO®

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Safety Notice	44
Unit Illustration	48
Important Testing Guidelines	51
Quick Start	52
Unit Operation	54
Battery Installation	54
System Settings	55
Applying The Wrist Monitor	56
Testing	58
Power Off	62
Memory Check	63
Memory Deletion	64
Low Battery Indicator	64
Troubleshooting	65
Blood Pressure Information	66
Blood Pressure Q&A	70
Maintenance	72
Specifications	76
Warranty	80

Safety Notice

Thank you for purchasing the Autotensio® Blood Pressure Monitor.

The unit has been constructed using reliable circuitry and durable materials. Used properly, this unit will provide years of satisfactory use.

This device is intended for non-invasive measuring of an adult individuals' systolic, diastolic blood pressure and heart rate using the oscillometric method. The device is not intended for use on infants and children. The device is designed for home or clinical use. All values can be read out in one LCD DISPLAY. Measurement position is on adult wrist only.

Please read this manual thoroughly before using the unit. Please retain this manual for future reference. For specific information about your blood pressure, please CONSULT YOUR DOCTOR.

To avoid risk and damage follow all warning precautions. Operate unit only as intended. Read all instructions prior to use.

WARNING SIGNS AND SYMBOLS USED	
	Caution
	Mandatory
	Prohibited
	Type BF Equipment
	Consult Instructions For Use
	Serial Number
	Discard the used product to the recycling collection point according to local regulations
	The product conforms to the requirements of the EC medical devices directive 93/42/EEC
	Manufacturer

Safety Notice

Caution

Individuals with serious circulation problems may experience discomfort. Consult your physician prior to use.

Contact your physician if test results regularly indicate abnormal readings. Do not attempt to self-treat these symptoms without consulting your physician first.

Make sure the electronic blood pressure doesn't stop the blood circulation.

Product is designed for its intended use only.
Do not misuse in any way.

Product is not intended for infants or individuals who cannot express their intentions.

Do not disassemble or attempt to repair.

Do not use cell phones and other devices, which generate strong electrical or electromagnetic fields, near the device, as they may cause incorrect readings and interference or become interference source to the device.

Caution

This device must not be used in presence of electrosurgery.

This device must not be used during patient transport.

Do not apply cuff on a wound, risk of more injuries.

Battery Precautions

Do not mix new and old batteries simultaneously.

Replace batteries when Low Battery Indicator “” appears on screen.


Be sure battery polarity is correct.

Do not mix battery types. Long-life alkaline batteries are recommended.


Remove batteries from device when not in operation for more than 3 months.

Dispose batteries properly; observe local laws and regulations.


Monitor Unit

Display


Contents

1. Monitor Unit

2. Owner's Manual

3. Plastic Storage Case

1. Avoid eating, exercising, and bathing for 30 minutes prior to testing.
2. Sit in a calm environment for at least 5 minutes prior to testing.
3. Do not stand while testing. Sit in a relaxed position while keeping your wrist level with your heart.
4. Avoid speaking or moving body parts while testing.
5. While testing, avoid strong electromagnetic interference such as microwave ovens and cell phones.
6. Wait 3 minutes or longer before re-testing.
7. Try to measure your blood pressure at the same time each day for consistency.
8. Test comparisons should only be made when monitor is used on the same wrist, in the same position, and at the same time of day.
9. This blood pressure monitor is not recommended for people with severe arrhythmia.

Quick Start

1. Install batteries. (See Figure A)


Figure A

2. Remove clothing from the wrist area. (See Figure B)
3. Rest for several minutes prior to testing. Wrap cuff around left wrist. (See Figure C)


Figure B


Figure C

4. Sit in a comfortable position and place wrist level with heart.
(See Figure D)
5. Press “

Figure D



Figure E

Battery Installation

Slide battery cover off as indicated by arrow.

Install 2 new AAA alkaline batteries according to polarity.

Close battery cover.

Note:

- 1) Replace batteries when Low Battery Indicator "☒" appears on screen.
- 2) Batteries should be removed from device when not in operation for an extended period of time.

1. Time /Date Setting

With power off, Press “

mode. Set the month first by adjusting the “M” button. Press “

Continue setting the day, hour, and minute in the same fashion. Every time the “A digital clock display showing "0:00 AM" in a small font at the bottom. Above the display is a large, empty rectangular frame with a thin black border, likely representing the screen of a device.

2. Save Settings


While in any setting mode, press “

Note: Unit will automatically save all information and shut off if left it for 3 minutes.

Applying The Wrist Monitor

Do not apply over clothing. If wearing a long sleeved shirt, be sure to roll sleeve back to forearm.

Apply monitor to wrist as illustrated. Tighten cuff firmly as not to wiggle.

Do not stand while testing. Sit in a comfortable position and place wrist level with heart.


Testing

1. Power On

Press button “” to turn the unit on.


The LCD screen will appear for one second as unit performs a quick diagnosis. A voice tone will indicate when unit is ready for testing.

Note: Unit will not function if residual air from previous testing is present in cuff. The LCD will flash “” until pressure is stabilized.

2. Pressurization


Initial pressure is first pumped to 190 mmHg. If the current user's systolic blood pressure is over 190mmHg, the unit will automatically re-inflate to the proper shelf.

Note: Pressurization will gradually subside and ultimately stop when cuff is not properly applied to the wrist.
If this occurs, pressing the “” button to turn the unit off.

3. Testing


After cuff inflation, air will slowly subside as indicated by the corresponding cuff pressure value. A flashing “” will appear simultaneously on screen signaling heart beat detection.

Note: Remain relaxed during testing. Avoid speaking or moving body parts.

4. Result Display

The screen will display measurements for systolic and diastolic blood pressure with voice broadcast. A indicator representing the current measurement will appear next to the corresponding WHO Classification.

Note: Refer to Page 27-28 for detail WHO Blood Pressure Classification Information.

Power Off

The “” button can be pressed to turn off the unit in any mode. The unit can turn off the power itself about 3 minutes no operation in any mode.


Safety Precaution: If pressure in cuff becomes too extreme while testing, press the “” button to turn power off.

The cuff pressure will rapidly dissipate once the unit is off.

Memory Check


With power off, you may check past test results by using the "M" buttons.

Upon activating test results, you can press the "M" buttons to scroll through all test results stored in memory. The LCD will display the last measurement memory with voice broadcast as NO:1 reading.

Memory Deletion

Memory may be deleted while in Memory Check mode. Press and hold the “” button for approximately 3 seconds to delete all memory records. With voice broadcast “Memory Clear” and then transfer into testing mode. Then press the “” button to turn the unit off.

Note: Memory cannot be recovered once it has been deleted.

Low Battery Indicator

The unit will broadcast “Low Battery” when battery life is depleting and unable to inflate cuff for testing. The “” appears simultaneously for approximately 5 seconds prior to shutting down. Replace batteries at this time. No memory loss will occur throughout this process.

Troubleshooting


Problem	Possible Cause	Solution
Blood pressure results are not within typical range	Cuff is too tight or not properly positioned on the wrist	Firmly reposition cuff on wrist making sure no wiggle is present. (See Page 16)
	Inaccurate test results due to body movement or monitor movement	Sit in a relaxed position placing wrist level with heart. Avoid speaking or moving body parts while testing. (See Page 9)
“Err” displayed	Cuff fails to inflate properly	Make sure hose is properly fastened to cuff and monitor unit.
	Improper operation	Read user manual carefully and re-test properly.
	Pressurization is over 300mmHg	Read user manual carefully and re-test properly.

Blood Pressure

Blood pressure is the force of blood pushing against the walls of arteries. It is typically measured in millimeters of mercury (mmHg.) Systolic blood pressure is the maximum force exerted against blood vessel walls each time the heart beats. Diastolic blood pressure is the force exerted on blood vessels when the heart is resting between beats.

An individual's blood pressure frequently changes throughout the course of a day. Excitement and tension can cause blood pressure to rise, while drinking alcohol and bathing can lower blood pressure. Certain hormones like adrenaline (which your body releases under stress) can cause blood vessels to constrict, leading to a rise in blood pressure.


If these measuring numbers become too high, it means the heart is working harder than it should.


WHO Blood Pressure Classification Indicator

The Autotensio® is equipped with a classification indicator based on established guidelines from the World Health Organization.

The chart below (color coded on monitor unit) indicates test results.


Mild Hypertension

■ : Blood Pressure Classification Indicator

Health Reminder

Hypertension is a dangerous disease that can affect the quality of life. It can lead to a lot of problems including heart failure, kidney failure, and cerebral hemorrhaging. By maintaining a healthy lifestyle and visiting your physician on a regular basis, hypertension and relative diseases are much easier to control when diagnosed in the early stages.

Note: Do not be alarmed if an abnormal reading occurs.

A better indication of an individual's blood pressure occurs after 2-3 readings are taken at the same time each day over an extended period of time. Consult your physician if test results remain abnormal.

Q: What is the difference between measuring blood pressure at home or at a professional healthcare clinic?

A: Blood pressure readings taken at home are now seen to give a more accurate account as they better reflect your daily life. Readings can be elevated when taken in a clinical or medical environment. This is known as White Coat Hypertension and may be caused by feeling anxious or nervous.

Note: Abnormal test results may be caused by:

1. Improper cuff placement
Make sure cuff is snug-not too tight or too loose.
2. Improper body position
Make sure to keep your body in an upright position.
3. Feeling anxious or nervous
Take 2-3 deep breaths, wait a few minutes and resume testing.

Q: What causes different readings?

A: Blood pressure varies throughout the course of a day. Many factors including diet, stress, cuff placement, etc. may affect an individual's blood pressure.


Q: Should I apply the cuff to the left or right wrist? What is the difference?

A: Either wrist can be used when testing, however, when comparing results, the same wrist should be used.
Testing on your left wrist may provide more accurate results as it is located closer to your heart.


Q: What is the best time of day for testing?

A: Morning time or any time you feel relaxed and stress free.


1. Avoid dropping, slamming, or throwing the unit.


2. Avoid extreme temperatures. Do not use outdoors.


- When cleaning the unit, use a soft fabric and lightly wipe with mild detergent. Use a damp cloth to remove dirt and excess detergent.


- Cuff Cleaning: Do not soak cuff in water! Apply a small amount of rubbing alcohol to a soft cloth to clean cuff's surface. Use a damp cloth (water-based) to wipe clean. Allow cuff to dry naturally at room temperature.


5. Do not use petrol, thinners or similar solvents.

6. Remove batteries when not in operation for an extended period of time.

7. Do not disassemble product.

8. It is recommended the performance should be checked every 2 years.
9. Expected service life: Approximately three years at 10 tests per day.

Specifications

Product Description	Wrist-type Fully Automatic Blood Pressure Monitor	
Model	Autotensio®	
Display	LCD Digital Display - Size:45mmx30mm	
Measurement Method	Oscillometric Method	
Measurement Range	Pressure	0mmHg ~ 300mmHg
	Pulse	30 to 180 Beats/Minute
Measurement Accuracy	Pressure	±3mmHg
	Pulse	±5%
Pressurization	Automatic Pressurization	
Memory	120 Memories with Date and Time	
Function	WHO Classification Indicator	
	Low Battery Detection	
	Automatic Power-Off	

Power Source	2 Alkaline Batteries Size AAA	
Battery Life	Approximately 2 months at 3 tests per day	
Unit Weight	Approx. 228g (8 oz.) (Excluding Battery)	
Unit Dimensions	Approx. 79 x 66 x 28mm	
Cuff Circumference	Approx.75 (W)×300(L) mm Fits wrist circumference 13.5-21.5 cm(5.3"-8.5")	
Operating Environment	Temperature	10°C ~ 40°C (50°F~104°F)
	Humidity	15% ~ 93% RH
	Pressure	700hPa~1060hPa
Storage Environment	Temperature	-25°C~70°C (-13°F~158°F)
	Humidity	≤93 % HR
Classification	Internal Powered Equipment, Type BF	

Specifications are subject to change without notice.

This Blood Pressure Monitor complies with the European regulations and bears the CE mark "CE 0459". This blood pressure monitor also complies with mainly following standards (included but not limited):

Safety standard:

EN 60601-1 Medical electrical equipment part 1: General requirements for safety

EMC standard:

EN 60601-1-2 Medical electrical equipment part 1-2: General requirements for safety- Collateral standard: Electromagnetic compatibility- Requirements and tests

Performance standards:

EN 1060-1 Non-invasive sphygmomanometers - General requirements

EN 1060-3 Non-invasive sphygmomanometers - Supplementary requirements for electromechanical blood pressure measuring systems.

EN 1060-4 Non-invasive sphygmomanometers - Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers.

Correct Disposal of This Product
(Waste Electrical & Electronic Equipment)

This marking shown on the product indicates that it should not be disposed with other household waste at the end of its life. To prevent potential harm to the environment or to human health, please separate this product from other types of wastes and recycle it responsibly. When disposing this type of product, contact the retailer where product was purchased or contact your local government office for details regarding how this item can be disposed in an environmentally safe recycling center. Business users should contact their supplier and check the terms and conditions of the purchasing agreement. This product should not be mixed with other commercial wastes for disposal. This product is free of hazardous materials.

Warranty

The Autotensio® Blood Pressure Monitor is guaranteed for 2-years from the date of purchase. If the Autotensio® does not function properly due to defective components or poor workmanship, we will repair or replace it freely. The warranty does not cover damages to your Autotensio® due to improper handling.

Please contact local retailer for details.

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SPENGLER SAS