EN

DE

FR

IT

ES

bloom lab

instructions for use

download the app to get started


www.bloomdiagnostics.com/app

What's in the box

REF BL-REF-2

1 instructions for use

1 Bloom Lab



1 USB-C cable

1 power supply

Quick start guide

1

2

[EN] See operating steps in the Bloom Test quick start guide

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

Table of contents

1. General information	3	8. Store-use operating steps	8
1.1. Intended use	3	9. Cleaning the Bloom Lab	9
1.2. Who should use the Bloom System?	3	10. Maintenance	9
1.3. Standards and legal requirements	3	10.1. Updating a Bloom Lab	9
2. Safety regulations	3	10.2. Setting up an administrator for a Bloom Lab	9
2.1. Disclaimer	3	10.3. Changing the administrator of a Bloom Lab	9
2.2. General notes on safety	4	11. Storage and transport	9
2.3. How to handle the Bloom Lab correctly	4	12. Disposal	9
2.4. Handling of bodily fluids	4	13. LED indications and error handling	9
2.5. Use of second-hand Bloom products	4	13.1. Proper function	10
2.6. Use of non-Bloom parts with the Bloom System	4	13.2. Error case	11
2.7. Damaged products or packaging	5	14. Troubleshooting	11
3. Terms and definitions	5	15. Technical details	13
4. Product information	5	15.1. Technical specification	13
4.1. Overview of the Bloom System	5	15.2. Operating principle	15
4.2. Bloom Lab	5	15.3. Performance	15
4.3. Bloom Test	5	15.4. Trademarks	15
4.4. Bloom App	6	16. Warranty	15
5. Bloom Lab components	6	17. Symbol information	15
6. Using the Bloom System for the first time	6	18. Contact information	15
7. General operating steps	6		

1. General information

1.1. Intended use

The Bloom Lab is intended for optically measuring Bloom Test strips and sharing the results with other Bloom devices, such as the Bloom App, for display and further analysis.

1.2. Who should use the Bloom System?

The Bloom System is intended to be used by those interested in preventative health tracking, self-optimisation, and tracking of medical conditions.

The Bloom System should only be used and operated by persons aged 18 years and over. No prior training is necessary.

1.3. Standards and legal requirements

The Bloom Lab is classified as an IVD according to the EC Council Directive 98/79/EC (IVDD), and conforms to the relevant provisions of the IVDD in accordance with Annex III. It further complies with the European directives 2011/65/EU (RoHS), 2012/19/EU (WEEE) and 2014/53/EU (RED). The Bloom Lab has been developed in accordance with the IEC 61010-1:2010+AMD1:2016, IEC 61010-2-101:2018, EN 61326-1:2013, EN 61326-2-6:2013, EN 60601-1-2:2014, EN 60601-1-11:2015 and EN 13532:2002.

2. Safety regulations

2.1. Disclaimer

The results obtained from using the Bloom System should not be regarded as a medical assessment.

The Bloom System provides the user with health guidance. It is not a diagnostic device. Despite all efforts taken to avoid such situations, false results may nevertheless be received. Therefore, the user should not make any decision of medical relevance based on the results received from the Bloom System without first consulting their medical practitioner.

When the Bloom System is used for tracking of existing medical conditions, the user should only adapt their treatment of that condition based on Bloom Test results if they have been trained to do so or in conjunction with a medical practitioner. Any treatment prescribed by a medical practitioner might require additional independent testing by a medical practitioner and/or laboratories.

2.2. General notes on safety

Do not take this test if you are under 18 years of age, if you are, or may be, pregnant, if you are intoxicated, or if you find yourself in an acute situation of emergency or distress. In case of doubt, always consult with your medical practitioner.

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Do not use the device after rapid temperature or humidity change in the environment as this could lead to an incorrect interpretation of your results. Please leave the device without power after rapid external temperature/humidity change.

The Bloom Lab is limited to indoor use.

If the Bloom Lab is used in a manner not specified by the manufacturer the degree of protection provided by the Bloom Lab may be impaired.

2.3. How to handle the Bloom Lab correctly

In order to avoid hazardous situations, the Bloom Lab should be handled with care at all times, this includes:

- The Bloom Lab should not be:
 - exposed to direct sunlight, physical instability, humidity, moisture, rapid temperature changes, strong electric/magnetic fields, strong power surges, dirt or dust;
 - exposed to high impact, as a result of falling or otherwise;
 - compressed or have weight placed upon it from any direction;
 - handled with force or pressure;
 - exposed to vibrations, movement or extreme changes in air pressure while in usage;
 - disassembled, except by trained and certified Bloom Diagnostics' technical support staff;
 - modified.
- The Bloom Lab should be used in the temperature range between +15°C to +40°C;
- The Bloom System should not be used with any hardware or software that is not manufactured or developed by Bloom Diagnostics for the explicit purpose of being used with the Bloom System. This includes using generic power cables or uncertified or generic strips;
- Bluetooth® connection or power supply should not be disconnected or lost while measurement, pairing, or usage of the Bloom Lab is ongoing.

2.4. Handling of bodily fluids

All bodily fluids should be treated as potentially infectious. If the Bloom Lab comes into contact with bodily fluids either from you, or from other people using the Bloom Lab, including family members and friends, you should use the provided Softa Cloth CHX 2% cleaning wipe to disinfect the Bloom Lab and any other surface that may have come into contact with bodily fluids.

2.5. Use of second-hand Bloom products

The Bloom Lab and entire Bloom System should only be bought from authorized retailers and never second-hand.

2.6. Use of non-Bloom parts with the Bloom System

The Bloom System should only be used with products specified by Bloom Diagnostics.

2.7. Damaged products or packaging

Do not use Bloom products, including all element of the Bloom System, if the products or their packaging appears damaged or shows evidence of tampering.

Do not attempt to deconstruct, damage, or tamper with any elements of the Bloom System as this can result in malfunction and will void the product's warranty.

3. Terms and definitions

Safety notes: Safety-related information that warns users about hazards and instructs on how to avoid them:

- **Warning:** A warning points out a situation where a potentially hazardous condition exists that could result in a critical or catastrophic situation such as death or serious long-term injury if not avoided.
- **Caution:** A caution points out a situation where a potentially hazardous condition exists that could result in a serious situation such as short-term injury or impairment requiring medical intervention to correct if not avoided.
- **Information:** Information points out a situation where a potentially hazardous condition exists that could result in a non-serious situation such as a slight inconvenience if not avoided.

4. Product information

4.1. Overview of the Bloom System

The Bloom System is a multi-part, testing system. It provides quantitative screening and monitoring of target analytes and health guidance. For this purpose, the Bloom System combines lateral flow assay-based bodily fluid testing using an established diagnostic technology (Bloom Test strips) with an electronic reader for quantification of the test line (Bloom Lab) and user-friendly software for displaying and interpreting results (Bloom App). See visual N°1.


Warning

The Bloom System does not function as a diagnostic device.

Using the information provided by the Bloom System for diagnosis can lead to mistreatment and severe consequences.

Do not use the Bloom System to diagnose or treat medical conditions.


Visual N°1: Bloom System

4.2. Bloom Lab

The Bloom Lab is an electronic device, which contains an optical measurement system to quantify the test line intensities produced by the Bloom Test strip during use. The Bloom Lab is intended to be used in combination with the Bloom App and Bloom Tests.

4.3. Bloom Test

The Bloom Test (Visual N°2) contains disposables which can be used for single-use blood sampling.

Visual N°2: Bloom Test

The Bloom Test strip (see Visual N°3) is sold as one integral part of the Bloom Test. The Bloom Test strip is a handheld, disposable single-use assay, which is designed to detect the concentration of a target analyte in a finger pricked blood sample and produces qualitative results in the form of a test line intensity. An additional control line indicates the proper analytical function of the Bloom Test strip.

Visual N°3: Bloom Test strip

4.4. Bloom App

The Bloom App is a smartphone app, which connects to the Bloom Lab via Bluetooth, interprets the quantified value and displays the measurement result to the user. The Bloom App is connected to an online database, owned and maintained by Bloom Diagnostics, which administers user information and uses data analysis to generate individual user insights.

5. Bloom Lab components

Caution: Only use equipment specified by Bloom Diagnostics to guarantee proper function of the device.


The Bloom Lab package includes (required):

- 1 x Bloom Lab
- 1 x Power supply (model no RYS518050150UE)
- 1 x USB-A to USB-C cable (150cm)
- 1 x Instructions for use


See Visuals N°4, 5, 6 and 7.

Not included (required):

- Bloom Test
- Bloom App (for a list of compatible smartphones see www.bloomdiagnostics.com)

Visual N°4: Bloom Lab

Visual N°5: USB-C cable

Visual N°6: Power supply

Visual N°7: Instructions for use

6. Using the Bloom System for the first time

Caution: Exceeding or not meeting the specified conditions can cause damage to the device.

The Bloom Lab works only in conjunction with the Bloom App and Bloom Tests, collectively known as the “Bloom System”. To get started, download and install the Bloom App on your smartphone and ensure you have access to a Bloom Test.

Visit www.bloomdiagnostics.com to learn more.

7. General operating steps

Caution: Do not use the Bloom Lab in rooms with high humidity such as the bathroom as this could lead to failures regarding the evaluation of the retrieved sample.

Caution: Do not use the Bloom Lab for self test without proper guidance from a healthcare professional.

In order to use the Bloom System, you must have a smartphone with a working internet connection and the installed Bloom App ready, as well as a Bloom Test. Check the Bloom Lab package components for completeness and any damage to the Bloom Lab.

Step 1: Set up Bloom Lab

Caution: Do not use the Bloom Lab if the packaging or the device are obviously damaged, otherwise correct functioning cannot be guaranteed.

Do not use the Bloom Lab when it is exposed to environmental conditions such as magnetic or electric fields, external electrical influences, electrostatic discharge, or thermal ignition sources.

Do not use the Bloom Lab with power supplies other than those provided by Bloom Diagnostics in order to avoid electronic malfunction of the device or electric shock.

1. Unpack the contents of the Bloom Lab package and check that all components have been received (see section 5.).
2. Plug in the Bloom Lab using the provided power supply. Wait until you see the white LED ring on the front lighting up and turning clockwise (see section 13.1.).
3. After less than 1 minute, the Bloom Lab will show a ready-state, indicated by the white LED ring being permanently lit (see section 13.1.).

Step 2: Pair your Bloom Lab with your Bloom App

1. Open your Bloom App on your smartphone.
2. Log into the Bloom App.
3. Click "Add Test".
4. Tap the upper right corner of your smartphone to the Bloom Lab (see visual N°8) or follow the instructions demonstrated in the Bloom App.

Visual N°8: Connect your Bloom App and Lab

Information: Successful pairing is shown by the white heart beat pattern on the Bloom Lab.

Step 3: Prepare to perform the test and collect your sample

Follow the steps according to the Bloom Test instructions for use.


Step 4: Analyze your retrieved sample

Warning
Using the information provided by the Bloom System for diagnosis can lead to mistreatment and severe consequences.
Do not use the Bloom System to diagnose or treat medical conditions.

Caution: Place the Bloom Lab on a stable, flat surface and do not move the device during measurement.
Please do not expose the Bloom Lab to direct sunlight during measurement, as this could influence your test results.

1. Insert the Bloom Test strip into the Bloom Lab (Visual N°9). You should feel the Bloom Test strip lock into place. The Bloom Lab will show a white heartbeat pattern on the LED ring when the strip is inserted correctly.

Visual N°9: Insert Bloom Test strip into the Bloom Lab

2. While the Bloom Lab analyses your sample, the Bloom Lab shows a spinning white loader and the Bloom App shows a countdown-timer.
3. After approximately 10 minutes the measurement will be completed. The LED ring on the Bloom Lab will pulse blue light if the measurement was successful. The Bloom App displays the measured value together with a personalized report.

Information: For any error cases see section 13.2..

Information: If the Bloom App disconnects from the Bloom Lab during the measurement process, the Bloom Lab completes the measurement process and saves the measurement value locally for a 3 hour period. If the same Bloom App reconnects with the Bloom Lab during this period, the saved measurement value will be sent to the Bloom App and the Bloom Lab will subsequently delete the measurement value and any unique identifier of the Bloom App. After a 10 minute period other Bloom Apps will be able to connect to the Bloom Lab and perform new tests, which will automatically delete your results.

Step 5: Switch off your Bloom Lab

1. Unplug the power supply from the device.

Step 6: Clean your Bloom Lab

Caution: Not cleaning the device as instructed can result in contamination of the device.

1. Use the cleaning wipe provided with each Bloom Test to wipe the Bloom Lab surface after each use.

8. Store-use operating steps

For the use of in-store Bloom Labs (device usable by multiple users), please follow the instructions below.

The following steps should be carried out by a store employee:

Step 1: Set up Bloom Lab

Caution: Do not use the Bloom Lab if the packaging or the device is obviously damaged. If you do so, the correct functioning cannot be guaranteed.

Caution: do not use the Bloom Lab when it is exposed to environmental conditions such as magnetic or electric fields, external electrical influences, electrostatic discharge, or thermal ignition sources.

Caution: Do not use the Bloom Lab with power supplies other than those provided by Bloom Diagnostics in order to avoid electronic malfunction of the device or electric shock.

Caution: Using the Bloom System in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets, etc.) may cause damaging electrostatic discharges that may produce erroneous results.

1. Unpack the contents of the Bloom Lab packaging.
2. Plug in the Bloom Lab using the provided power supply. You will see the white LED ring on the front lighting up and turning clockwise (see section 13.1.). This will take less than 1 minute.
3. Once the Bloom Lab is ready this will be indicated by the white LED ring lighting up permanently (for pictures see section 13.1.).

Step 2: Setting up an administrator device

1. Open the Bloom App on the smartphone you want to use as the administrator device.
2. Log into the Bloom App.
3. Click "Add Test".
4. Tap the upper right corner of your smartphone to the Bloom Lab. When pairing with the Bloom Lab is successful, you will see a white heartbeat pattern on the LED ring (see section 13.1.).
5. The first smartphone connected to the Bloom Lab is set as the administrator device. The administrator device can be used to update the Bloom Lab.

Optional step: Changing the administrator device

1. To reset the administrator device please press the reset button on the bottom of the device (Visual N°10).
2. After resetting the device, the first smartphone connected with the Bloom App will be defined as the administrator device.

Step 3: Renaming the Bloom Lab

1. After connecting the Bloom Lab to the Bloom App, please rename the Bloom Lab to give it an individual identifier.

To use the Bloom Lab in-store users should now perform steps 2–4 of section 7.

Step 4: Clean the Bloom Lab

Caution: Please use protective equipment, such as gloves, when cleaning to avoid contamination.

Caution: Not cleaning the device as instructed can result in contamination of the device.

1. Use the cleaning wipe provided with each Bloom Test to thoroughly wipe all surfaces of the Bloom Lab after each use.
2. Check gaps of the Bloom Lab casing for stains and clean again when needed.

Step 5: Switch off your Bloom Lab

1. Unplug the power supply from the device after the completion of the last measurement.

9. Cleaning the Bloom Lab

Caution: Please use protective equipment, such as gloves, when cleaning the surfaces of the Bloom Lab to avoid contamination.

Thoroughly wipe all surfaces of the Bloom Lab after every use using a disinfection cleaning wipe, Braun Softa Cloth CHX 2% (Ref.: 19581) or similar. Check gaps for stains and clean thoroughly when needed.

10. Maintenance

10.1. Updating a Bloom Lab


The Bloom App informs you about any existing updates and guides you through installation. The installation can only be done by the paired administrator device of the Bloom Lab using a Bluetooth connection.

10.2. Setting up an administrator for a Bloom Lab

The first Bloom App connected with the Bloom Lab is defined as the administrator device, allowing it to update your Bloom Lab. This device is also allowed to rename the Bloom Lab to give it a unique identifier.

10.3. Changing the administrator of a Bloom Lab

To reset the Bloom Lab, please press the reset button on the bottom of the Bloom Lab for more than 5 seconds (Visual N°10).

Visual N°10: Reset button on the Bloom Lab

After resetting the device, the first Bloom App to be connected to it is defined as the administrator device.

11. Storage and transport

In order to avoid hazardous situations, the Bloom Lab should be stored with care at all times, this includes:

- The recommended storage temperature of the Bloom

Lab is between -20°C and $+60^{\circ}\text{C}$. Please store the Bloom Lab in a dry place, away from direct or excessive heat. Make sure to store the Bloom Lab in an environment with humidity between 30 and 85%.

- The Bloom Lab should not be stored outside: it is intended for indoor use only.
- The Bloom Lab should not be stored in: brightly sunlit, moist, dirty, dusty, physically unstable, compressed, or humid areas, or areas that experience rapid temperature changes, areas that have strong magnetic fields (e.g. created by multiple mobile phones, batteries, magnets, high voltage power transmitters, radioactive materials, microwaves, TVs), or areas that experience strong power surges.

12. Disposal


Dispose of the Bloom Lab in accordance with your national regulations regarding electronic waste. Please contact your local authority to determine the proper method of disposal of potentially biohazardous parts.

Before disposal or transportation the unit shall be cleaned and disinfected. Use Braun Softa Cloth CHX 2% (Ref.: 19581) or similar.

Do not dispose the Bloom Lab as unsorted municipal waste. Every person should reuse, recycle and recover electronic waste responsibly to avoid hazardous substances affecting health and the environment.

13. LED indications and error handling

The Bloom Lab performs a self test during each start-up of the device. Additional internal tests are performed during each measurement to check for the correct function of the device. The proper function of the Bloom Lab is indicated via the LED indicator at the front of the device. See section 13.1. for more information. In case of malfunction, the user is informed by the Bloom App and the LED indicator (Visual N°11).

Visual N° 11:
LED indications on
the Bloom Lab

13.1. Proper function

The correct function and readiness of the Bloom Lab are indicated by the LED ring. Please see Table N° 1 for the meaning of the various indications:

Bloom Lab is turning on

A white light circle is closing representing the start-up progress of the Bloom Lab.
Duration: <1 minute.

Bloom Lab is ready

A dimmed but stable ring of white light awaits user interaction.

1.

Pairing success

A white light merging together from both ends to form a circle indicates that the Bloom App and the Bloom Lab are pairing. Successful pairing is notified by a white heart beat pattern. Duration: <10 seconds.

2.

Waiting for sample collection

A white, soft, heartbeat pulse indicates that the Bloom Lab is waiting for a sample to be collected.

Confirmation of sample transfer

The LED ring displays a white half moon to signal that the blood sample has been transferred to the Bloom Test strip.

1.

Strip inserted

A white light merging together from both ends to form a circle indicates that the Bloom Lab recognizes that a Bloom Test strip has been inserted. Successful insertion is notified by a white heartbeat pattern. Duration: <10 seconds.

2.

Measurement in progress

A spinning, white loader is shown.

New results available

A blue, circular light will pulse until the new results have been read in the Bloom App.

Table N°1: LED ring indications: proper function

13.2. Error case

Errors are indicated via the Bloom Lab's LED ring in red. Error messages and their meanings are described in Table N° 2. Further information about the error is displayed in the connected Bloom App.

Recoverable errors

A blinking red and white ring means that something went wrong. The user is invited to check their Bloom App to learn how to solve the problem.

Fatal error

A steady, bold red LED circle indicates that an error has occurred which will require servicing of the Bloom Lab to rectify (for contact details see section 18.).

Table N°2: LED ring indications: error case

In the case that the LED indication differs from the error message demonstrated in the Bloom App, please contact our technical support (see section 18.)

14. Troubleshooting

Report fatal errors and malfunctions of the Bloom Lab to Bloom by contacting our technical support (see section 18.).

The following table (Table N°3) highlights possible errors you might face when performing Bloom Tests, determines possible malfunction messages, and provides instructions on how to correct those problems.

Bloom Lab test processing		
Error	Possible cause	Action to correct the error
Caution: The result is lost if the connection is re-established with another device. Bluetooth connection lost or cannot be established.	Too large distance or obstructions between Bloom Lab and smartphone.	Bring your smartphone close to the Bloom Lab and ensure there are no obstructions in between. Re-establish the Bluetooth connection.
	Bloom Lab or smartphone switched off.	Check to make sure that your Bloom Lab is connected to a power supply and your smartphone is switched on. Re-establish the Bluetooth connection.
	Bluetooth not enabled on smartphone.	Check that your smartphone's Bluetooth is enabled and ready to pair. Re-establish the Bluetooth connection.
	Bloom Lab not responding.	Disconnect and reconnect the Bloom Lab from the power supply and wait for the boot process to be finished. Establish the Bluetooth connection.

Bloom Lab shows recoverable error.	Bloom Test strip was inserted earlier than expected.	Remove the Bloom Test strip from the Bloom Lab and follow the steps according to the Bloom App.
	QR code on the Bloom Test strip cannot be read.	Remove the Bloom Test strip from the Bloom Lab and follow the steps according to the Bloom App. Try inserting the Bloom Test slower or faster.
	Bloom Test strip has been used already or does not match the information provided on the registered test package.	Remove the Bloom Test strip from the Bloom Lab and follow the steps according to the Bloom App.
	Bloom Test strip was not inserted properly.	Remove the Bloom Test strip from the Bloom Lab and follow the steps according to the Bloom App.
	Bloom Test strip was removed during measurement.	Re-insert the Bloom Test strip into the Bloom Lab and follow the steps according to the Bloom App.
	Bloom Test strip inserted too late for measurement.	Please repeat the testing procedure by using a new test. Information: When a Bloom Test strip is inserted later than expected, the sample can not be evaluated.
	The temperature that the Bloom Lab is exposed to is above 40°C.	Ensure an ambient temperature below +40°C. Take care that the Bloom Lab is not exposed to direct sunlight.
	The temperature the Bloom Lab is exposed to is below 15°C.	Ensure an ambient temperature above 15°C.
Bloom Lab shows fatal error.	Bloom Lab internal error, not recoverable.	Contact technical support for assistance.
Bloom Lab does not work.	Bloom Lab has no power supply.	Make sure that the Bloom Lab is plugged in as instructed.
	Bloom Lab internal error.	Contact technical support for assistance.
	Bloom Lab not responding.	Disconnect and reconnect the Bloom Lab from the power supply and wait for the boot process to be finished.

Table N°3: Troubleshooting

15. Technical details

15.1. Technical specifications

Power supply	100-240VAC (50/60 Hz) / 5VDC, 1.5 A via USB-C
Size	100x115x110mm (WxLxH)
Weight	400g
IP class	IP51 (dust and dripping water protected)
Environmental	+15°C to +40°C (Measurement) -20°C to +60°C (Storage) <80% Relative Humidity
Communication	Bluetooth Low Energy 4.2 Module Type Class 2 NFC , ISO/IEC 14443A compliant, 106 kbps transmission rate
Firmware Version	1.1
End of life	5 years

Table N°4: Technical specification

The Bloom Lab complies with the electromagnetic interference requirements of EN 61326-2-6:2013 and EN 61326-1:2013. Further, it complies with the particular requirements for electromagnetic compatibility for medical devices according to EN 60601-1-2 (see Table N°5 for immunity and emission levels).

Caution: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Bloom Lab, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

The Bloom System should be used in a home or professional healthcare environment.

Emission	
Test	Limits
Conducted emission	CISPR 11, Group 1, Class B
Radiated emission	CISPR 11, Group 1, Class B

Harmonic current emission	IEC 61000-3-2, Class A, power < 75W
Voltage fluctuations and flicker	IEC 61000-3-3
Immunity test levels	
Test	Test levels
Electrostatic discharge (IEC 61000-4-2)	Contact Discharge: ± 8 kV Air Discharge: ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV
Radiated RF EM fields (IEC 61000-4-3)	80-2700 MHz; 1kHz AM 80 %; 10 V/m
Proximity fields from RF wireless communications equipment (IEC 61000-4-3)	385 MHz; Pulse Modulation: 18 Hz; 27 V/m 450 MHz; FM: ± 5 kHz Hub, 1 kHz Sinus; 28 V/m 710, 745, 780 MHz; Pulse Modulation: 217 Hz; 9 V/m 810, 870, 930 MHz; Pulse Modulation: 18 Hz; 28 V/m 1720, 1845, 1970 MHz; Pulse Modulation: 217 Hz; 28 V/m 2450 MHz; Pulse Modulation: 217 Hz; 28 V/m; 5240, 5500, 5785 MHz; Pulse Modulation: 217 Hz; 9 V/m
Electrical fast transients and bursts (IEC 61000-4-4)	Power supply lines: ± 2 kV, 100 kHz repetition freq. Signal / Data lines: ± 1 kV, 100 kHz repetition freq.
Surge Immunity (IEC 61000-4-5)	0.5, 1 kV; Differential mode Surge Waveform: Combination Wave ($2 \mu\text{s} \times 50 \mu\text{s}$ Voltage, $8 \mu\text{s} \times 20 \mu\text{s}$ Current)
Conducted disturbances induced by RF fields (IEC 61000-4-6)	0.15-80 MHz; 1 kHz AM 80 %; 3 Vrms, 6 Vrms ISM band
Rated power frequency magnetic fields (IEC 61000-4-8)	30 A/m, 50 Hz / 60 Hz

Table N°5: EMI test levels

15.2. Operating principle

The Bloom Lab is an electronic device, which contains an optical measurement system (optical reader) to quantify the test line intensities produced by the Bloom Test strip during use. It captures the Bloom Test strip read-out window, automatically detects the test and control line and calculates the mean intensity value for each line by image analysis. This mean value, in connection with further image analysis, corresponds to the concentration of the captured target analyte. The actual concentration of the target analyte is calculated differently for each batch of Bloom Tests using the Bloom App, thus enabling measurements with the Bloom Lab for a wide variety of analytes.

15.3. Performance

For the analyte specific performance please see the instructions for use of each Bloom Test.

15.4. Trademarks



The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Bloom Diagnostics is under license. Other trademarks and trade names are those of their respective owners.

16. Warranty

For consumers, who are covered by consumer protection laws or regulations in their country of purchase or, if different, their country of residence, the benefits conferred by Bloom Diagnostics' Limited Warranty are in addition to all rights and remedies conveyed by such consumer protection laws and regulations, including but not limited to these additional rights.

For further information on the rights provided by consumer law visit: www.bloomdiagnostics.com.

17. Symbol information

	Manufacturer		CE Mark
	Do not dispose of waste of electrical and electronic equipment as unsorted municipal waste		Consult instructions for use
	Fragile, handle with care		Keep dry
	Keep away from sunlight		Do not use if package is damaged
	Humidity exposure limits between 30% and 80% for device operation		Temperature exposure limits between +15°C and +40°C for device operation
	Atmospheric pressure exposure limits between 700 hPa to 1060 hPa for device operation		Humidity exposure limits between 30% and 85% for transport and storage
	Temperature exposure limits between -20°C and +60°C for transport and storage		Atmospheric pressure exposure limits between 700 hPa to 1060 hPa for transport and storage
	In vitro diagnostic medical device		IP51 (dust and dripping water protected)
	Serial Number		Catalog number

18. Contact information

 Bloom Diagnostics AG
Susenbergstrasse 185
8044 Zürich, Switzerland

For assistance or technical support contact:
Email - support@bloomdiagnostics.com
Web - www.bloomdiagnostics.com/contact

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