

Ziegelei 1 72336 Balingen-Frommern Germany

KERN & Sohn GmbH www.kern-sohn.com

+0049-[0]7433-9933-0

+0049-[0]7433-9933-149

info@kern-sohn.com

Operating instructions Personal balances

KERN MPN

TMPN 200K-1HM-A TMPN 200K-1M-A **TMPN 200K-1PM-A TMPN 300K-1LM-A**

Version 1.3 2022-06 GB



D	Weitere Sprachversionen finden Sie online unter www.kern-sohn.com/manuals
BG	Други езикови версии ще намерите в сайта <u>www.kern-sohn.com/manuals</u>
DK	Flere sprogudgaver findes på websiden www.kern-sohn.com/manuals
EST	Muud keeleversioonid leiate Te leheküljel www.kern-sohn.com/manuals
E	Más versiones de idiomas se encuentran online bajo www.kern-sohn.com/manuals
GR	Άλλες γλωσσικές αποδόσεις θα βρείτε στην ιστοσελίδα <u>www.kern-sohn.com/manuals</u> s
F	Vous trouverez d'autres versions de langue online sous www.kern-sohn.com/manuals
LV	Citas valodu versijas atradīsiet vietnē www.kern-sohn.com/manuals
FIN	Muut kieliversiot löytyvät osoitteesta www.kern-sohn.com/manuals
LT	Kitas kalbines versijas rasite svetainėje <u>www.kern-sohn.com/manuals</u>
GB	Further language versions you will find online under www.kern-sohn.com/manuals
RO	Alte versiuni lingvistice veţi găţi pe site-ul www.kern-sohn.com/manuals
ı	Trovate altre versioni di lingue online in www.kern-sohn.com/manuals
SK	Iné jazykové verzie nájdete na stránke <u>www.kern-sohn.com/manuals</u>
NL	Bijkomende taalversies vindt u online op www.kern-sohn.com/manuals
SLO	Druge jezikovne različice na voljo na spletni strani <u>www.kern-sohn.com/manuals</u>
Р	Encontram-se online mais versões de línguas em www.kern-sohn.com/manuals
CZ	Jiné jazykové verze najdete na stránkách www.kern-sohn.com/manuals
PL	Inne wersje językowe znajdą Państwo na stronie <u>www.kern-sohn.com/manuals</u>
SE	Övriga språkversioner finns här: www.kern-sohn.com/manuals
Н	A további nyelvi változatok a következő oldalon találhatók: www.kern-sohn.com/manuals
HR	Druge jezične verzije su dostupne na stranici: www.kern-sohn.com/manuals
NO	Andre språkversjoner finnes det på www.kern-sohn.com/manuals



KERN MPN

Version 1.3 2022-06

Operating instructions Personal balances with BMI function

Contents

1 1.1	Tolerances height meter	
2 2.1	Declaration of conformity Explanation of the graphic symbols for medical products	
3 3.1 3.2	Appliance overview Overview of displays Keyboard overview	11 14
4.1 4.1.1 4.1.2 4.2 4.3 4.4 4.5 4.6 4.7	Basic Information (General) Specific function Indication Contraindication Proper use Non-intended product use / contraindications Warranty Monitoring of Test Resources Plausibility check. Reporting serious incidents	16 16 16 17 18 18
5 5.1 5.2 5.3 5.4	Basic Safety Precautions Pay attention to the instructions in the Operation Manual. Personnel training Preventing contamination Preparation for use	20 20 20
6 6.1 6.2 6.2.1 6.3 6.3.1 6.4	Electromagnetic compatibility (EMC) General hints Electromagnetic emission of interferences Drop in performance Electromagnetic interference immunity Crucial features of performance Minimum distances	21 23 23 24
7 7.1 7.1.1	Transport and storage Testing upon acceptance Packaging / return transport	28
8 8.1 8.2 8.3 8.4 8.5 8.6 8.7	Unpacking, Installation and Commissioning Installation Site, Location of Use Unpacking Scope of delivery Balance assembly and installation Fix the body height measuring rod Battery operation Rechargeable battery operation using an optional battery power pack	29 29 30 31
5.7	resolar godolo battory operation doing an optional battory power pack	

8.8	Mains connection	34
8.9	Initial Commissioning	34
9	Operation	35
9.1	Weighing	
9.2	Taring	35
9.2.1	Subsequent tare weight	
9.2.2	Pretare	
9.3	HOLD function	
9.4	Calculation of the Body Mass Index	
9.4.1 9.4.2	Determing the body height (MPN-HM-A only)	
9.4.2 9.4.3	Calculating Body Mass IndexClassification of BMI values	
9. 4 .5 9.5	Automatic switch-off function "AUTO OFF"	
9.6	Signal tone when pressing button	
9.7	Setting time and date	
10	Menu	
10.1	Navigation in the menu	
10.2	Menu overview	46
11	Communication with peripheral devices via KUP connection	48
11.1	KERN Communications Protocol (KERN Interface Protocol)	49
11.2	Data output after pressing the PRINT button < manual >	
11.3	Automatic data output <auto></auto>	50
11.4	Continuous data output < cont >	51
11.5	Data format	52
11.6	WLAN	53
11.7	Print function	54
12	Bluetooth	55
13	Alibi memory optional	55
14	Error messages	56
15	Servicing, maintenance, disposal	
15.1	Cleaning	
15.2	Cleaning / disinfecting	
15.3	Sterilisation	
15.4	Servicing, maintenance	
15.5	Disposal	57
16	Instant help for troubleshooting	58
17	Verification	59
17.1	Verification validity period (current status in D)	
	,	
18	Adjustment	
18.1	Adjustment switch and seals	62

1 Technical data

KERN	MPN 200K-1HM	MPN 200K-1PM		
Item no./ Type	TMPN	TMPN		
Display	200K-1HM-A 200K-1PM-A 6-digit			
Weighing range (max)	250 kg			
Minimum load (Min)	2 kg			
Readability (d)	0.1			
Verification value (e)	100	g		
Accuracy for first verification	≤50 kg = >50 kg-200 >200 kg – 25) kg = 1 e		
Linearity ±	0.1			
Display	LCD with 25m	m high digits		
Recommended adjustment weight, not added, (Category)	200 (M:	•		
Stabilization time (typical)	3 se	ec.		
Warm-up time	10 min			
Operating temperature	10° C + 40° C			
Storage and transportation environment	-20 to +60°C, and 30% to 90% relative humidity			
Humidity of air	max. 80 % (not condensing)			
Atmospheric pressure (kPa)	70kpa-106kpa 6 V / 1 A			
Input Voltage	6 V / 6 x 1.5			
Battery use	Battery working range: 48 hours backgroud illumination off 24 hours background illumination on Loading time: 8 hours			
Auto Off	off, after 30 s / 1, 2, 5, 30, 60 min			
Dimensions fully mounted (W x D x H) mm	365 x 570 x 2134	365 x 570 x 1030		
Weighing plate (W x D x H) mm	365 x 36	60 x 80		
Weight kg (net)	11.5	10.8		
Verification in accordance with 2014/31/EU	Catego	ory III		
Medical product in accordance with 93/42/EEC	Category Im (with measuring function)			
Rechargeable battery operation (optional)	optional; 3.8 VDC – 4.2 VDC / 3700 mAh			
Height measuring rod in stand integrated, extendable (from 3 cm to 205 cm)	-			
Data interface	Intern: Wi-Fi Optional / Extern: KUP (RS232, Bluetooth, USB-D, Extension box)			

KERN	MPN 200K-1M	MPN 300K-1LM		
Item no./ Type	TMPN 200K-1M-A	TMPN 300K-1LM-A		
Display	6-digit			
Weighing range (max)	250 kg	300 kg		
Minimum load (Min)	2 kg	2 kg		
Readability (d)	0.1 kg	0.1 kg		
Verification value (e)	0.1 kg	0.1 kg		
Accuracy for first verification	>50 kg-20	= 0.5 e 0 kg = 1 e 50 kg = 1.5 e		
Linearity ±	0.1 kg	0.1 kg		
Display	LCD with 25n	nm high digits		
Recommended adjustment weight, (Category)	200 kg (M1)	200 kg (M1)		
Stabilization time (typical)	3 sec.	3 sec.		
Warm-up time	10 min	10 min		
Operating temperature	10° C + 40° C			
Storage and transportation environment	-20 to +60°C, and 30% t	to 90% relative humidity		
Humidity of air	max. 80 % (not condensing)			
Atmospheric pressure (kPa)	70kpa-	106kpa		
Input Voltage	6V	1A		
Battery use	6 x 1.5 V AA Battery working range: 48 hours background illumination off 24 hours background illumination on Loading time: 8 hours			
Auto Off	off, after 30 s / 1, 2, 5, 30, 6 (adjus	60 min without load change stable)		
Weighing plate (W x D x H) mm	365 x 360 x 80	400 x 500 x 120		
Weight kg (net)	8,4	10,0		
Verification in accordance with 2014/31/EU	Category III			
Medical product in accordance with 93/42/EEC				
Rechargeable battery operation (optional)	optional; 3.8 VDC - 4.2 VDC / 3700 mAh			
Data interface Optional / Extern: KUP (RS232, Bluetooth, USB-D, Extension box)				

1.1 Tolerances height meter

Measured value (cm)	Tolerance (cm)
< 90	± 0.5
100	± 1.0
150	± 1.0
200	± 1.0

2 Declaration of conformity

The current EC/EU Conformity declaration can be found online in:


www.kern-sohn.com/ce



For verified weighing scales (= weighing scales assessed for conformity) the declaration of conformity is included in the scope of delivery.

Only these balances are medical products.

2.1 Explanation of the graphic symbols for medical products



All medical balances with this mark fulfill the following guidelines:

- 1. 2014/31/EU: Guideline for non-automatic balances
- 2. 93/42/EEC: Guideline for medical products



Balances which carry this mark, are conformity-evaluated as per accuracy class III of the EC-guideline 2014/31/EU. Accuracy of the balance see chapter 1 Technical data.

WF 170012

Serial number of the device, applied at the device and on the packaging

Number here as example



Identification of the manufacturing date of the medical product.

Year and month here as example



Please observe the attached documents" or "Please observe the operating instructions"



Observe operating instructions



Observe operating instructions



Identification of manufacturer of medical product including address

Kern & Sohn GmbH D-72336 Balingen,Germany www.kern-sohn.com



Electro-medical device with supplement for type B



Device protection class II



Dispose of old appliances separately from your household waste!

Instead, take them to communal collection points.




Display of supply voltage for scales with polarity display



Supply voltage direct current



Information



Level balance before use



Seal mark KERN SEAL




Mains connection



Electrostatically endangered structural components

10 TMPN_A-BA-e-2213

3 Appliance overview






MPN-M-A



MPN-LM-A

Secondary display at rear (MPN-HM-A and MPN-PM-A)





3.1 Overview of displays

Display	Designation	Description
	Stability display	Scales are in a steady state
→0 ←	Zeroing display	Should the balance not display exactly zero despite empty weighing plate, press the button. Your balance will be set to zero after a short standby time.
NET	Net weight display	Illuminated when net weight is displayed Illuminated after weighing scale was tared
GROSS	Gross weight display	Illuminated when gross weight is displayed
HOLD	HOLD function	HOLD function active
ВМІ	BMI function	Illuminated while BMI function is enabled
	Battery symbol	Shows the charging capacity of the batteries
kg	Weighing unit	Displays the weighing unit
*	WiFi-interface	Shows the connection to the cableless network

3.2 Keyboard overview



Button	Designation	Function
ON OFF	ON/OFF-button	Turn on/off
HOLD	HOLD button	Hold function/Calculation of a stable weight value In menu: • Select menu items For numeric entry: • Lowering numeric value
ВМІ	BMI button	Calculation of the Body Mass Index In menu: Return to weighing mode
PRINT	Print button	Data transfer via interface In menu: • Select menu items For numeric entry: • Increase numerical value
F	Function key	Function quick key Quick call-up of a function stored before
TARE →0←	Zeroing key	Weighing scale will be reset to "0.0" In menu: Confirm selection For numeric entry: Change decimal digit Confirm entry
	TARE key	Tare balance

4 Basic Information (General)



Balances have to be verified for the purposes stated below in accordance with Directive 2014/31/EU. Article 1, paragraph 4. "Determination of mass in the practice of medicine that is, weighing patients for reasons of medical supervision during medical surveillance, examination and treatment."

4.1 Specific function

4.1.1 Indication

- Determining the body weight in the medical practice area
- Using as "non-automatic balance"
- The person steps carefully on the center of the weighing platform.

Once a steady display value is shown, you can read the weight result.

4.1.2 Contraindication

No contraindication known.

4.2 Proper use

This weighing scale is designed for determining the weight of a person whilst standing, such as in medical treatment rooms. The regularly used function of the balance consists of recognizing, prevention and treatment of illnesses.

• On personal weighing scales, the person should step onto the center of the weighing platform and remain standing without moving.

Once a steady display value is shown, you can read the weight result. The weighing scale is designed for continuous duty.



The weighing platform may only be stepped on by persons capable of standing on both feet on the weighing platform.

- The weighing platforms are fitted with an anti-slip surface that must not be covered during weighing a person.
- The balance should be checked for correct condition prior to each utilization by a person familiar with proper operation of the balance.
- When using balances with mounted body height measuring rod, ensure that the top flap is turned downwards immediately after use in order to avoid risk of injury.

The WIFI interface allows a wireless transfer of the measurement results to a PC.



Scales fitted with a serial interface may only be connected to appliances in compliance with Directive EN60601-1.

If the balance doesn't have any contact with the transfer cable, do not touch the transfer port in order to avoid an ESD-failure.





4.3 Non-intended product use / contraindications



- Do not leave permanent load on the weighing pan. This may damage the measuring system.
- Impacts and overloading exceeding the stated maximum load (max) of the weighing plate, minus a possibly existing tare load, must be strictly avoided. This could cause damage to the balance.
- Never operate the balance in explosive environment. The serial version is not explosion protected. It should be noted that a flammable mixture of anaesthetics and oxygen or laughing gas may occur.
- The structure of the balance may not be modified. This may lead to incorrect weighing results, safety-related faults and destruction of the balance.
- The balance may only be used according to the described conditions. Other areas of use must be released by KERN in writing.
- If the balance is not used for a longer time, take out the batteries and store them separately. Leaking battery liquid could damage the balance.
- The balance may only be used for weighing persons.
 Persons who are heavier than the indicated maximum load, may not step onto the balance.



Non-intended use of the optional body height measuring rod



- The body height measuring rod may only be assembled as specified in the operating instructions.
- The structure of the body height measuring rod may not be modified. This may result in incorrect measuring results, safety-related defects as well as destruction.
- The body height measuring rod may only be used according to the described conditions. Other areas of use must be released by KERN in writing. For more details please see the user manuals of the body height measuring rod.

4.4 Warranty

Warranty claims shall be voided in case:

- Our conditions in the operation manual are ignored
- The appliance is used beyond the described uses
- The appliance is modified or opened
- Mechanical damage and damage caused by media, liquids,
- Natural wear and tear
- Improper installation or faulty electrical connection
- The measuring system is overloaded
- Dropping the balance

4.5 Monitoring of Test Resources

In the framework of quality assurance the measuring-related weighing properties of the balance and, if applicable, the testing weight, must be checked regularly. The responsible user must define a suitable interval as well as type and scope of this test. Information is available on KERN's home page (www.kern-sohn.com) with regard to the monitoring of balance test substances and the test weights required for this. In KERN's accredited DKD calibration laboratory test weights and balances may be calibrated (return to the national standard) fast and at moderate cost.

For personal balances with body height measuring rods, we recommend a metrological examination of the accuracy of the body height measuring rod, however, this is not mandatory as the determination of human body height involves rather large, intrinsic inaccuracies.

4.6 Plausibility check

Please make sure that the measurement values computed by the appliance are plausible and are allocated to the respective patient, before storing and using the values for further purposes. This applies especially also for values transferred via interface.

4.7 Reporting serious incidents

All serious incidents appeared related to this product must be reported to the manufacturer and the responsible authority of the member state where the user and/or the patient are residents.

"Serious incident" that means an incident which directly or indirectly had, could have or could have had one of the following consequences:

- > the death of a patient, a user or another person,
- the temporary or permanent fatal deterioration of the health status of a patient, a user or other persons,
- a serious danger for public health.

5 Basic Safety Precautions

5.1 Pay attention to the instructions in the Operation Manual



 □ Carefully read this operation manual before setup and commissioning, even if you are already familiar with KERN balances.



5.2 Personnel training

The medical staff must apply and follow the operating instructions for proper use and care of the product.

The balance must be set up via the interfaces and integrated into a network only by experienced administrators or hospital technicians.

5.3 Preventing contamination

The prevention of cross-contamination (fungal skin infections,.....) requires regular cleaning of the weighing platform. Recommendation: After any weighing procedure that could potentially result in contamination (e. g. after weighing that involves direct skin contact).

5.4 Preparation for use

- Check the personal balance for damage before any use
- Maintenance and reverification (in Germany MTK): The personal balance must be serviced and reverified at regular intervals.
- Do not use the appliance on slippery surfaces or in facilities with risk of vibration
- During installation the personal balance must be levelled
- If possible, the product must remain in its original packaging for transportation purpose. Should this not be possible, make sure that the product is protected against damage
- Step onto and leave the personal balance only when a qualified person is present

6 Electromagnetic compatibility (EMC)

6.1 General hints

This device complies with the limits set for medical electrical devices of group 1, class B (as per EN 60601-1-2). The device is suitable for home healthcare and professional health institutes environments.



The installation and use of this electrical medical device requires special precautionary measures as outlined in the EMC information below.



Do not install the appliance near active surgical high-frequency devices and in radio-frequency-screened rooms of a ME system for magnetic resonance reproduction where high intensity of electro-magnetic interferences occurs.



Please avoid to operate the appliance beside or stacked on other devices, as this may cause incorrect measuring results. If such use should be required, this appliance and the other devices should be observed, in order to ensure that they work normally.



Using accessories, transformers and other cables than the specified ones or delivered by the manufacturer together with the appliance, could have as consequence reinforced electromagnetic radiation or reduced electromagnetic immunity to interference and by that way reduced functionability.



Portable radio-frequency communication equipment (including periphery as well as antenna cable and external antennae) should be placed from any part of the MPN (including the cables authorised by manufacturer) at a distance of at least 30 cm (12 inches). Otherwise the performance of the appliance could drop down.

Note: The emission characteristics of this appliance allow its use in industrial areas and in hospitals (CISPR 11 category A). If it is used in residential areas (where CISPR 11 category B normal is required), this appliance cannot ensure sufficient protection against radio-frequency-communication services. In order to obtain reasonable results, the user should apply measures for weakening the intensity, e.g. install the device in another site or re-align it.

Electromagnetic compatibility (EMC) describes a device's ability to perform reliably within an electromagnetic environment without causing there inadmissible electromagnetic interference at the same time. Amongst other things, such disturbances may be transmitted by connecting cables or by the air.

Inadmissible interferences from the environment may result in wrong displays, inaccurate measured values or incorrect behaviour of the medical device. The

performance regulation is less than ±1kg instable reading when measuring with evaluated weight capacity.

In the same way, the personal balance MPN in certain cases may cause such disturbances in other devices. To eliminate problems of that kind, we recommend you to take one or several of the measures listed below:

- Change the alignment or distance of the device to the source of EMI.
- Install or use the personal balance MPN at another location.
- Connect the personal balance MPN to another power source.
- For further questions please contact our customer services.

Disturbances may be caused by improper modification or add-ons to the device or not recommended accessories (such as power supply units or connecting cables). The manufacturer will not be responsible for these. Modifications may also result in a loss of authorisation relating to the use of the device.



Devices emitting high frequency signals (mobile telephones, radio transmitters, radio receivers) may cause interference in the medical device. For that reason do not use them near the medical device. Chapter 6.4 contains details about recommended minimum distances.

6.2 Electromagnetic emission of interferences

All required instructions to preserve the BASIC SAFETY and REQUIRED PERFORMANCE considering the electromagnetic interferences for the expected service life.

The tables below refer to the product operated by mains current

Guidelines and manufacturer's declaration - Electromagnetic emissions

The personal balance MPN is designed for use in an electromagnetic environment that meets the requirements stated below. The customer or user of the personal balance MPN should ensure that operation takes place in such an environment.

Emission check	Fulfillment
Radio-frequency emissions	Assembly 1
CISPR 11	
Radio-frequency emissions	Category [A]
CISPR 11	
Harmonic emissions	Category A
IEC 61000-3-2	
Voltage fluctuations/Flickering	Keeping
IEC 61000-3-3	

Do not operate the personal balance MPN directly next to other appliances or stacked with other devices. If this type of operation cannot be avoided, observe the personal balance MPN to ensure normal operation in such an arrangement.

6.2.1 Drop in performance



Powerful fields with electromagnetic interferences which e.g. are arising from electric motors or inductive charging devices, may cause drop down of performance if they are located near the personal balance MPN. Drop down of performance may cause instably displayed weighing values.

6.3 Electromagnetic interference immunity

Guidelines and manufacturer's declaration - electromagnetic interference immunity

The personal balance MPN is designed for use in an electromagnetic environment that meets the requirements stated below. The customer or user of the personal balance MPN should ensure that operation takes place in such an environment.

Immunity to interference tests	IEC 60601-1-2 Test level	Fulfillment level		
Electrostatic discharge (ESD)	± 8 kV contact	± 8 kV contact		
IEC 61000-4-2	± 2 kV, ± 4 kV, ± 8 kV, ± 15kV air	±2 kV,± 4 kV, ±8kV, ±15kV air		
Electrical quick compensation / Burst	± 2 kV for power lines	±2 kV for power supply lines		
IEC 61000-4-4	± 1 kV Signal input/output	±1 kV signal input/output		
1EC 61000-4-4	100kHz repetition frequency	100 kHz repetition frequency		
Excess voltage	±0.5 kV, ±1 kV differential mode	±0.5 kV, ±1 kV differential mode		
IEC 61000-4-5	±0.5 kV, ±1 kV, ±2 kV General mode	±0.5 kV, ±1 kV, ±2 kV common mode		
Voltage low points,	0 % UT; 0.5 cycle. At 0°, 45°, 90°,	0 % UT; 0.5 cycle. At 0°, 45°, 90°,		
short interruption and	135°, 180°, 225°, 270° and 315°.	135°, 180°, 225°, 270° and 315°.		
voltage fluctuations in				
the power supply cables	0 % UT; 1 cycle and 70 % UT;	0 % UT; 1 cycle and 70 % UT;		
IEC 61000-4-11	25/30 cycles; single phase: at 0°.	25/30 cycles; single phase: at 0°.		
	0 % UT; 250/300 cycle	0 % UT; 250/300 cycle		
Immunity to	30A/m	30A/m		
interferences against magnetic fields with energy-technical frequencies	50Hz/60Hz	50Hz/60Hz		
IEC 61000-4-8				

Interference against wire-guided interference variables, inducted by high-frequency fields IEC 61000-4-6	3V 0.15MHz - 80MHz 6V in the ISM and in the amateur radio strips between 0.15MHz and 80MHz 80% AM at 1kHZ	3V 0.15MHz - 80MHz 6V in the ISM and in the amateur radio strips between 0.15MHz and 80MHz 80% AM at 1kHZ
Testing the immunity to interferences against high-frequency electromagnetic fields	10V/m 80MHz – 2.7GHz 80% AM at 1kHZ	10V/m 80MHz – 2.7GHz 80% AM at 1kHZ
IEC 61000-4-3		

 $\textbf{Note} \hbox{:}\ U_T\hbox{ is the mains AC voltage before applying the test level}.$

iated	Test fre- quency (MHz)	Strip (MHz)	Service	Modulation	Maximum Power (W)	Dis- tance (m)	IEC 60601-1-2 Test Level (V/m)	IMMUNITY TO INTER- FERENCE TEST LEVEL (V/m)
cy 00-	385	380 – 390	TETRA 400	Impulse Modulation 18 Hz	1,8	0.3	27	27
i -	450	430 – 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28	28
	710 745 780	704 – 787	LTE strip 13, 17	Impulse Modulation 217 Hz	0,2	0.3	9	9
, _	810 870 930	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE strip 5	Impulse Modulation 18 Hz	2	0.3	28	28
	1720 1845 1970	1 700 - 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE strip 1, 3, 4, 25; UMTS	Impulse Modulation 217 Hz	2	0.3	28	28
	2450	2 400 - 2 570	Bluetooth WLAN 802.11 b/g/n, RFID 2450, LTE strip 7	Impulse Modulation 217 Hz	2	0.3	28	28
	5240 5500 5785	5 100 - 5 800	WLAN 802.11 a/n	Impulse Modulation 217 Hz	0,2	0.3	9	9

The field intensity of stationary radio transmitters such as base stations of wireless telephones and mobile radio sets, amateur radio stations, AM and FM radio and television stations cannot be reliably predicted in advance. To determine the electromagnetic environment in respect of stationary transmitters, you should consider a study of electromagnetic phenomena at the location. If the measured field strength of the location where the personal balance MPN is operated, is exceeding the above mentioned conformity level, the personal balance MPN should be observed in order to make sure the intended function. If you observe unusual features of performance you may have to take additional measures such as a change in alignment or a different location for the medical device.

Note: The EMISSION characteristics of this appliance allow its use in industrial areas and in hospitals (CISPR 11 category A). If it is used in residential areas (where CISPR 11 category B normal is required), this appliance cannot ensure sufficient protection against radio-frequency-communication services. In order to obtain reasonable results, the user should apply measures for weakening the intensity, e.g. install the device in another site or re-align it.

6.3.1 Crucial features of performance



The personal balance MPN does not have any crucial features of performance as per IEC 60601-1. The system may be subject to interference by other devices even if these devices conform to current emission requirements as per CISPR.

6.4 Minimum distances

Recommended safety distances between portable and mobile HF telecommunication devices and the medical device

The personal balance MPN is designed for use in an electromagnetic environment in which high-frequency disturbance variables are controlled. The customer or user of the personal balance MPN can help avoiding electromagnetic disturbances by keeping the minimum distance between portable and mobile high-frequency telecommunication devices (transmitters) and the medical device – depending on the output performance of the communication equipment, as stated below.

Rated capacity of transmitter %W	The safety distance depends on the transmission frequency %m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.20	1.20	2.30
10	3.80	3.80	7.30
100	12.00	12.00	23.00

For transmitters with a maximum rated capacity not stated in the table above you can calculate the recommended safety distance in meters (m) yourself by using the equation belonging to each column, whereby P equals the maximum rated capacity of the transmitter in Watt (W) as per details provided by the transmitter manufacturer.

NOTE 1 Higher frequency range applies to 80 MHz and 800 MHz.

NOTE 2 These guidelines may not be applicable in all cases.

The spreading of electromagnetic variables is influenced by absorption and reflections from buildings, objects and persons.

7 Transport and storage

7.1 Testing upon acceptance

When receiving the appliance, please check packaging immediately, and the appliance itself when unpacking for possible visible damage.

7.1.1 Packaging / return transport



- ⇒ Keep all parts of the original packaging for a possibly required return.
- ⇒ Only use original packaging for returning.
- ⇒ Prior to dispatch disconnect all cables and remove loose/mobile parts.
- ⇒ Reattach possibly supplied transport securing devices.
- ⇒ Secure all parts such as the weighing plate, power supply unit etc. against shifting and damage.

8 Unpacking, Installation and Commissioning

8.1 Installation Site, Location of Use

The balances are designed in a way that reliable weighing results are achieved in common conditions of use.

You will work accurately and fast, if you select the right location for your balance.

On the installation site observe the following:

- Place the balance on a stable, even surface
- Avoid extreme heat as well as temperature fluctuation caused by installing next to a radiator or in the direct sunlight
- Protect the balance against direct draughts due to open windows and doors
- Avoid jarring during weighing
- Protect the balance against high humidity, vapors and dust
- Do not expose the device to extreme dampness for longer periods of time. Nonpermitted condensation (condensation of air humidity on the appliance) may occur if a cold appliance is taken to a considerably warmer environment. In this case, disconnect the appliance from mains and acclimatize it for approx. 2 hours at room temperature.
- Avoid static charge of the balance and of the person to be weighed.
- Avoid contact with water.

Major display deviations (incorrect weighing results) may be experienced, should electromagnetic fields (e.g. due to mobile phones or radio equipment), static electricity accumulations or instable power supply occur. Change location or remove source of interference.

8.2 Unpacking

Take the balance out of its packaging and place it at the intended position. When using the power supply unit, ensure that the power cable does not produce any risk of stumbling.


8.3 Scope of delivery

- Balance
- Mains adapter (in conformity with EN 60601-1)
- Protective hood
- Wall fixture (only for models TMPN-1M-A and TMPN-1LM-A)
- Operating instructions


8.4 Balance assembly and installation


Assembly:

Ensure that the cable (3) does not get clamped!



⇒ Fasten stand, using 4 screws





- ⇒ Level balance with foot screws until the air bubble of the water balance is in the prescribed circle.
- ⇒ Check levelling regularly



Once installation is complete, check all screws for tight fit. Otherwise the person to be weighed may suffer injury.

General instructions for installing the aforementioned balances:

Install the personal balance at the intended location and align it using the screwed height-adjustable rubber feet, until the air bubble of the bubble level (on the weighing plate) is in the center.

Wall fixture for models TMPN-1M-A and TMPN-1LM-A:




1	Screws for affixing the wall fixture on the display unit
2	Position of the screws to affix the display unit to the wall

8.5 Fix the body height measuring rod

The force necessary for extending the telescopic height meter can be adapted with the help of two adjustment screws on the stand (see fig.).

For this purpose proceed as follows:



- Using a suitable slotted screw driver turn the two adjustment screws to adjust the desired effort
- (if frequently used, this procedure possibly has to be repeated after a certain time)







Once installation is complete, check all screws for tight fit. Otherwise the person to be weighed may suffer injury.

8.6 Battery operation

As an alternative for the rechargeable battery operation, the balance offers also the possibility to be operated with 6x AA-batteries.

Open the battery cover (1) at the lower side of the display unit and insert the batteries according to the example shown below. Lock again the battery compartment cover.

If the batteries are exhausted, the symbol and "Lo bAt" is displayed on the balance; replace the batteries. To save the battery, the balance switches automatically off (see chap. 9.5).

	Capacity of batteries exhausted.	
+		
Lo bAt		
	Capacity of batteries will soon be exhausted.	
	Batteries completely charged	

Disclaimer:

Only use KERN type rechargeable battery YMR-01 (RC 193650) or batteries AA 1.5 V (6x). Other may cause damages to the product, which can lead to injuries of persons.



If the balance is not used for a longer time, take out the batteries and store them separately. Leaking battery liquid could damage the balance.



Replacement for rechargeable battery by inadequately trained personnel could result in a hazard.

Insert batteries:

Remove the battery compartment lid at the side of the display unit (1)



Remove the battery holder (2)



Insert the batteries into the battery holder



Insert batteries with the battery holder into the battery compartment and lock it with battery compartment cover.



8.7 Rechargeable battery operation using an optional battery power pack



Open the battery compartment lid (1) at the base of the display unit and insert the rechargeable battery. Charge the rechargeable battery for at least 12 hours before initial use.

The appearance of the symbol in the weight display indicates that the rechargeable battery is almost exhausted. The weighing scale will remain ready for operation for a few more minutes before switching off automatically in order to save the rechargeable battery (see chap. 9.5). Charge rechargeable battery.

Voltage has dropped below prescribed minimum.
Rechargeable battery very low.
Rechargeable battery completely recharged



• If the rechargeable battery is exhausted, "LoBAt" is displayed. The rechargeable battery is charged via the provided plug-in power supply unit (charging time 12 h for a complete charging).

- If the balance is not used for a longer time, take out the rechargeable battery and store it separately. Leaking liquid could damage the balance.
- When the optional WIFI-interface is used, power consumption increases

8.8 Mains connection

Power is supplied by the external mains adapter which also serves to isolate the mains supply from the scale. The stated voltage value must be the same as the local voltage. Always use genuine approved KERN power pack units as per directive EN 60601-1.

If the symbol appears in the display, the rechargeable battery will soon be exhausted. Connect the power supply unit and charge the rechargeable battery.

The flashing symbol informs you during charging about the charging status of the rechargeable battery.


8.9 Initial Commissioning

In order to obtain exact results with the electronic balances, your balance must have reached the operating temperature (see warming up time chap.1). During this warm-up time the balance must be connected to the power supply (mains, rechargeable battery) and be switched on.

The accuracy of the balance depends on the local acceleration of gravity. The value of gravity acceleration is shown on the type plate.

9 **Operation**

Weighing





- At any moment you can reset the balance to zero by pressing TARE **→**0← button.
- ⇒ Have person stand in the center of the scales. Wait until the stability display \(\simeq \) appears, then read the weighing result.




(=overload) will appear in the display.


If the person is heavier than the weighing range,

9.2 Taring

The tare weight of any preloads can be deducted by pressing a button so that the actual weight of the person is displayed in subsequent weighings.




Place the object on the weighing plate.



Press I the zero display appears. "NET" is shown at the bottom on the left.

TARE

→0←



 \Rightarrow Allow the person to step onto the center of the weighing platform.

Wait until the stability display **appears**, then read the weighing result.



- When the balance is unloaded the saved tare value is displayed with negative sign.
- To delete the stored tare value, unload the balance and press


9.2.1 Subsequent tare weight

The balance can be tared several times successively.


9.2.2 Pretare

The possibility exists to enter either a known pre-tare value via the keyboard or to store the weight of an object placed on the weighing plate as a pre-tare value.

Store the weight of an object placed on the weighing plate:



Enter the pre-tare value manually via the keyboard:



Delete pre-tare value:


Unload the weighing plate and press, the balance goes to zero display.

TARE

HOLD

9.3 HOLD function

The balance has an integrated standstill function (mean value calculation). With this function it is possible to weigh people accurately even if they do not stand still on the weighing plate.





You need to know a person's body height before you can calculate the BMI for that person. It should either be known or can be determined directly with the MPN-HM-A model.

9.4.1 Determing the body height (MPN-HM-A only)



- ⇒ Push measuring rod upwards and set the headpiece horizontally.
- ⇒ Push measuring rod carefully down until the headpiece touches the person's head (measuring without shoes).



A fixed headpiece pointing outwards poses a risk of injury.




⇒ Read body height on measuring stick.



A correctly exercised height measurement will achieve an accuracy of up to 5 mm.

9.4.2 Calculating Body Mass Index





- Reliable calculation of BMI is restricted to a body height of 100 cm to 200 cm and a weight of >10 kg.
- If weighing has to take place under unsteady conditions, the display can be stabilised via the Hold function.

9.4.3 Classification of BMI values


Weight classification for adults over 18 years of age using the BMI in accordance with WHO, 2000 EK IV and WHO 2004.


Category	BMI (kg/m²)	Risk of diseases associated with overweight
Underweight	< 18.5	Low
Normal weight	18.5 – 24.9	Average
Overweight	<u>></u> 25.0	
Pre-adipose	25.0 – 29.9	A bit increased
Adipose degree I	30.0 – 34.9	Increased
Adipose degree II	35.0 – 39.9	High
Adipose degree III	<u>≥</u> 40	Very high

9.5 Automatic switch-off function "AUTO OFF"


The weighing scale will switch off automatically after the allotted time as long as neither the display unit nor the weighing plate is operated.








To set a determined switch-off time, proceed as follows:




[2 min]	Weighing system will be turned off after 2 min.
[5 min]	Weighing system will be turned off after 5 min.
[30 min]	Weighing system will be turned off after 30 min.
[60 min]	Weighing system will be turned off after 60 min.
[30 S]	Weighing system will be turned off after 30 sec.
[1 min]	Weighing system will be turned off after 1 min.



9.6 Signal tone when pressing button

The balance offers the possibility to switch off or on an acoustic signal by pressing the keys.






9.7 Setting time and date

(only available with Real Time Clock)



Menu settings:

[dAtiME]



Setting date:

TARE

- ⇒ In menu call-up the menu item "dAtiME"
- Press repeatedly until the prompt for entering the year appears "YY20xx". Enter the current year with or TARE and confirm by
- ⇒ The display changes automatically to the prompt for entering the month: "MM xx"
- Enter the current month with or and confirm by
 - The display changes to the prompt for entering the day:

 "dd xx". Enter the current day with

 TARE

 confirm by

Setting time:

- Now the display changes automatically to the prompt for entering the time, first the hour: "hh xx".
- ⇒ Enter the hour with or PRINT and confirm by
- The minutes appear: "MM xx". Enter the current minutes and confirm on
- The seconds appear: "SS xx". Enter the seconds and confirm on
- □ Date and time have now been entered; the display changes to the weighing mode.

10 Menu



Access to service menu "x10" is locked in verified balances.

To disable the access lock, destroy the seal mark and actuate the adjustment switch. Position of the adjustment switch see chap. 18.1

Attention:

After destruction of the seal mark the weighing system must be reverified by an authorized agency and a new seal mark fitted before it can be reused for applications subject to verification.

10.1 Navigation in the menu

Call up menu	⇒ In weighing mode press simultaneously, "SetuP" followed by "coM" will appear.
Select function	⇒ With help of or PRINT, the individual functions can be selected one after the other.
Change settings	 Confirm selected function by will be displayed. ⇒ Select desired setting with PRINT and confirm by TARE, the balance returns to the menu.
Exit menu/ Return to weighing mode	⇒ Press repeatedly until the zero display appears. The balance is now in weighing mode.

10.2 Menu overview

Menu block Main menu	Menu item Submenu	Available settings / explanation	
SEtuP			
соМ	r\$232	RS232 mode	
Interface parameter	bAUd	Baud rate 9600, 14400, 19200, 38400, 57600, 115200, 128000, 256000, 600, 1200, 2400, 4800	
	dAtA	Data bits: 8dbitS, 7dbitS	
	PAritY	Parity nonE, odd, EVEn	
	StoP	Stop bits: 1Sbit, 2SbitS	
	HAndSh	Handshake: nonE	
	Protoc	Communication Protocol: KCP,	
	uSb-d	USB interface	
	bAud	Baud rate 9600, 14400, 19200, 38400, 57600, 115200, 128000, 256000, 600, 1200, 2400, 4800	
	dAtA	Data bits: 8dbitS, 7dbitS	
	PAritY	Parity nonE, odd, EVEn	
	StoP	Stop bits: 1Sbit, 2SbitS	
	HAndSh	Handshake: nonE	
	Protoc	Communication Protocol: KCP,	
	WLAn	WLAN interface: on, oFF	
Print	intFcE	Interface	
Printer settings	rS232	RS232 interface	
	uSb-d	USB interface	
	PrModE	Printer settings	
	MAnUAL	The weighing value is printed out when pressing	
	Auto	The weighing value is automatically printed out	
	cont	Continuous data output	
	ForMAt	long , short	
	LAYout	USEr	

bEEPEr	KEVE	Signal tang when proceing button
Signal tone	KEYS	Signal tone when pressing button
	on	Acoustic signal on
		Acoustic signal off
	oFF	
A	Mada	Settings autom. switch-off function
AutoFF Automatic switch-	Mode	
on function	oFF	Automatic switch-off function switched off
	Auto	Automatic switch-off independent on the load on the balance
	onLY0	Automatic switch-off only with zero display
	timE	Setting the switch-off time
	30 S, xMin	Automatic switch off after: 30 sec, 1 min, 2 min, 5 min, 30 min, 60 min
dAt iME		Setting date and time
Date and time	YY.2021	Setting date: Year
	MM 06	Month
	dd 22	Date
	hh 11	Setting time: Hour
	MM 53	Minutes
	SS 33	Seconds
	· ·	
rESEt back to default setting		Balance is reset to default setting

11 Communication with peripheral devices via KUP connection

Via the interfaces weighing data may be exchanged with connected peripheral devices. Issue may be made to a printer, PC or check displays. In reverse order, control orders and data inputs may be made via the connected devices.

The balances of the TMPN series are equipped with a KUP connection (KERN Universal Port) as per standard.


The following options are available as interfaces:

	Interface adapter with cable	
	Model Example for application	
RS232	YKUP-01	Serial printer
USB	YKUP-03	PC
Ethernet	YKUP-04	PC
Bluetooth	YKUP-06	Android terminal unit or PC
Kern Extension Box	YKUP-13	Several interfaces in parallel

The available interfaces can be used via the KUP (YKUP-13) in parallel manner.

Additional equipment connected to medical electrical equipment must comply with the respective IEC or ISO standards (e.g. IEC 60950 for data processing equipment). Furthermore all configurations shall comply with the requirements for medical electrical systems (See IEC 60601-1 or clause 16 of the 3Ed. of IEC 60601-1, respectively). Anybody connecting additional equipment to medical electrical equipment configures a medical system complies with the requirements for medical electrical systems. Attention is drawn to the fact that local laws take priority over the above mentioned requirements. If in doubt, consult your local representative or the technical service department.

Connector assignment balance:



Warning: Only for use with KUP interfaces



Externally attached third party interface cables, which are connected to the KUP interfaces, may not exceed 10 m in length.

11.1 KERN Communications Protocol (KERN Interface Protocol)

KCP is a standardized set of interface orders for KERN balances, which allows many parameters and device functions to be called up and controlled. KERN devices that have KCP can use it to connect easily to computers, industrial control systems and other digital systems. A detailed description you will find in the "KERN Communications Protocol" manual, available in the download area on our KERN homepage (www.kern-sohn.com).

To activate KCP please observe the menu overview of your balance's operating instructions.

KCP is based on simple ASCII orders and replies. Every interaction consists of an order, possibly with arguments separated by spaces and finished by <CR>< LF>.

The KCP orders supported by your balance may be queried emitting the order "I0" followed by CR LF.

Extract of the mostly used KCP orders:

10	Shows all implemented KCP orders
S	Sending stable value
SI	Sending current value (also instable)
SIR	Sending current value (also instable) and repeating
Т	Taring
Z	Zeroing

Example:

Order	S	
Possible replies	S_S100.00_g S_I S_+ or S	Order accepted, execution of the order started, currently another order is executed, timeout reached, over- or underload

11.2 Data output after pressing the PRINT button < manual > Activate function:

- ⇒ In Setup menu invoke the menu setting <Pr → PrnodE> and confirm with → button.
- \Rightarrow Use the navigation keys $\downarrow\uparrow$ to select the setting $< \Box \Box >$ and confirm on \Rightarrow button.
- ⇒ To exit the menu press the navigation key ← repeatedly.

Place goods to be weighed on balance:

- ⇒ If required, place empty container on scale and tare.
- ⇒ Place goods to be weighed. The weighing value is edited by pressing the PRINT-button.

11.3 Automatic data output <auto>

Data output happens automatically without having to press the **PRINT**-key as soon as the corresponding output condition has been met, dependent on the setting in the menu.

Enable function and set the output condition:

- ⇒ In Setup menu invoke the menu setting <Pr → Pr□□□E> and confirm with → button.
- ⇒ For an automatic data output select the menu setting < ☐□□□ > using the navigation keys ↓↑ and confirm by the → button.
- ⇒ Acknowledge by →-button and set the required output condition with the navigation keys ↓↑.
- ⇒ Acknowledge by →-button.
- ⇒ To exit the menu press the navigation key ← repeatedly.

Place goods to be weighed on balance:

- ⇒ If required, place empty container on scale and tare.
- ⇒ Place weighed goods and wait until the stability display (► △) appears. The weighing value is issued automatically.

11.4 Continuous data output < cont >

Enable function and set the output interval:

- ⇒ In Setup menu invoke the menu setting <Pr → Pr∏□dE> and confirm on button.
- Use the navigation keys ↓↑ to select the setting < □□ > and confirm on → button.
- ⇒ <5PEEd> is displayed.
- Acknowledge by →-button and set the required time interval using the navigation keys ↓↑
- ⇒ To exit the menu press the navigation key ← repeatedly.

Place goods to be weighed on balance

- ⇒ If required, place empty container on scale and tare.
- ⇒ Place goods to be weighed.
- ⇒ The weighing values are issued according to the defined interval.

Sample log (KERN YKB-01N):

```
S D 1.9997 kg
S D 1.9999 kg
S D 1.9999 kg
S D 1.9999 kg
S D 1.9999 kg
S S 2.0000 kg
S S 2.0000 kg
S S 2.0000 kg
S D 1.9998 kg
S D 1.9998 kg
S D 1.9998 kg
S D 2.0002 kg
S D 2.4189 kg
S D 2.9996 kg
S D 2.9996 kg
S D 2.9996 kg
S D 2.9997 kg
S D 2.9996 kg
S D 2.9997 kg
S D 2.9996 kg
```



The balances which have a serial interface, may only be connected to electrical office machines.

11.5 Data format

- ⇒ In Setup menu invoke the menu setting < Pr int → Pr∏odE > and confirm with button →.
- ⇒ Use the navigation keys \$1\$ to select the menu setting < F□□□□□□□ > and confirm on → button.
- - < 与hort > Standard measuring protocol
 - < L□□□ > Detailed measuring protocol
- ⇒ Confirm setting with →-button.
- ⇒ To exit the menu press the navigation key ← repeatedly.

Sample log (KERN YKB-01N):


For	1AE → Shor	-E	ForNAL → Lon	<u> </u>
N: T: G:	5 5	2.0000 kg 0.5000 kg 2.5000 kg	N: 5 D Tara weight after x Gross weight:	2.0000 kg 0.5000 kg 2.5000 kg

11.6 WLAN


- > W-LAN standard: IEEE 802.11 b/g/n (Wi-Fi)
- Network log: TCP/IP with DHCP
- > Supported encryption methods: WPA, WPA2
- > Transmission frequency: 2412 2472 MHz
- > Maximum transmission performance: < 20dBm
- > Application log: KCP (KERN Communications Protocol):

Set up WIFI connections:


- 1. The balance creates a WIFI access point as soon as it has started up (WLAN symbol in the balance display appears).
 - Use your computer to connect to this access point.
 - The SSID (name of the balance's access point) is "Al THINKER xxxxxxx"
- 2. Using a web browser visit the website http://192.168.4.1/. In the website:
 - A. Set the mode "Mode" to "apsta"
 - B. Enter information about the network which you want to integrate the balance into (network "AP Name" and password "AP Password")
 - C. Save the settings "Save" and update the site



- 3. Separate the access point from the computer
- 4. Cut the power supply to the balance for a short moment
- 5. Connect the computer again to the access point of the balance and update the website
 - D. Now the IP-address "IP address" is displayed



- 6. Close the website
- 7. Connect the computer to the selected network
- 8. Enter the IP / Port into the target software: 23



11.7 Print function

After a correct connection of the software and the balance, the weighing data can be transferred using the button on the balance.

12 Bluetooth

The balance is optionally available with **Bluetooth Low Energy** (BLE) and will then be visible to Bluetooth Master devices using its serial number.

To access this, please use an appropriate software program / app which supports Bluetooth Low Energy (BLE). Applications exclusively using Bluetooth Classic (BLC) will not work.

The following profile must be adjusted:

Service UUID

0000fff0-0000-1000-8000-00805f9b34fb

Read characteristic UUID

0000fff1-0000-1000-8000-00805f9b34fb

Write characteristic UUID

0000fff2-0000-1000-8000-00805f9b34fb

- > Transmission frequency: 2402 2480 MHz
- > Maximum transmission performance: < 20dBm

13 Alibi memory optional

The balance stores weighing data after pressing in the alibi memory. They can be retrieved from there at any time.

PRINT

14 Error messages

Display	Description
	Battery capacity exhausted
LobAt	Battery capacity soon exhausted
[]Eroh (Zero setting range exceeded
ZEroLo	Zero setting range not achieved
AronG	Adjustment error
(uPFHP)	Load instable
	Underload
[]	Overload

Should other error messages occur, switch balance off and then on again. If the error message remains inform manufacturer.

15 Servicing, maintenance, disposal

15.1 Cleaning



Before any maintenance, cleaning and repair work disconnect the appliance from the operating voltage.

15.2 Cleaning / disinfecting

Clean weighing platform (such as seat pan) as well as casing with household detergents or commercially available disinfectants, e.g. 70% isopropanol. We recommend a disinfectant suitable for wiping disinfection. Please follow manufacturer's instructions.

Do not use abrasive or aggressive cleaners such as spirits or alcohol or similar as they might damage the high-quality surface.

To prevent cross-contamination (fungal skin infection) please observe the following time intervals for disinfection:

- Weighing plate before and after any measurement with direct skin contact
- When required:
 - Display
 - Touch-sensitive keyboard



Do not spray the device with disinfectant, just wipe it.

Make sure that disinfectant does not penetrate the interior of the balance.

Remove dirt immediately.

15.3 Sterilisation

Sterilisation of the appliance not allowed.

15.4 Servicing, maintenance

The appliance may only be opened by trained service technicians who are authorized by KERN.

We recommend a regular safety-related technical check (STK) by a qualified expert. Disconnect scales from mains before opening.

15.5 Disposal

Disposal of packaging and appliance must be carried out by operator according to valid national or regional law of the location where the appliance is used.

16 Instant help for troubleshooting

In case of an error in the program process, briefly turn off the balance and disconnect from power supply. The weighing process must then be restarted from the beginning.

Fault

Possible cause

The weight display does not • The balance is not switched on glow

- The mains supply connection has been interrupted (mains cable not plugged in/faulty).
- Power supply interrupted.
- The rechargeable battery / the batteries is/ are inserted incorrectly or empty
- No rechargeable battery / no battery is/ are inserted

The displayed weight permanently changing

- is Draught/air movement
 - Table/floor vibrations
 - The weighing plate is in contact with foreign bodies or is not correctly positioned.
 - Electromagnetic fields / static charging (choose different location/switch off interfering device if possible)

The weighing result obviously incorrect

- is The display of the balance is not at zero
 - Adjustment is no longer correct
 - Great fluctuations in temperature.
 - Warm-up time was ignored.
 - Electromagnetic fields / static charging (choose different location/switch off interfering device if possible)

the WIFI interface

- No data can be transferred to The mains signal is not stable or too weak
 - Wrong interface configuration

Should other error messages occur, switch balance off and then on again. If the error message remains inform manufacturer.

17 Verification

General:

According to EU directive 2014/31/EU balances must be officially verified if they are used as follows (legally controlled area):

- a) For commercial transactions if the price of goods is determined by weighing.
- b) For the production of medicines in pharmacies as well as for analyses in the medical and pharmaceutical laboratory.
- c) For official purposes
- d) For manufacturing final packages
- e) Determination of mass in the practice of medicine that is, weighing patients for reasons of medical supervision during medical surveillance, examination and treatment,

In cases of doubt, please contact your local trade in standard.

Verification notes:

An EU type approval exists for balances described in their technical data as verifiable. If the balance is used where obligation to verify exists as described above, it must be verified and re-verified at regular intervals.

Re-verification of a balance is carried out according to the respective national regulations. For verification validity period, see chap. 17.1.

The legal regulation of the country where the balance is used must be observed!



Verification of the balance is invalid without the seal.

The seal marks attached on balances with type approval point out that the balance may only be opened and serviced by trained and authorized specialist staff. If the seal mark is destroyed, verification looses its validity. Please observe all national laws and legal regulations. In Germany a re-verification will be necessary.

Balances with obligation to verify must be taken out of operation if:

- The **weighing result** of the balance is outside **the error limit.** Therefore, in regular intervals load balance with known test weight (ca. 1/3 of the max. load) and compare with displayed value.
- The reverification deadline has been exceeded.

17.1 Verification validity period (current status in D)

Personal scales (including chair and wheelchair scales) in hospitals	4 years
Personal balances, when not located in hospitals (for example, doctor's offices and nursing homes)	unlimited
Baby weighing scales and mechanical birth weight scales	4 years
Bed scales	2 years
Scales in dialysis stations	unlimited

Notes:

- Also rehab clinics and health authorities are treated as hospitals
- Not treated as hospitals (verification validity not limited) are dialysis stations, nursing homes and doctor's surgeries.

(Data source : "Bureau of Standards News, Weighing Instruments in Medicine" Original title "Bureau of Standards News, Weighing Instruments in Medicine"}

18 Adjustment

As the acceleration value due to gravity is not the same at every location on earth, each display unit with connected weighing plate must be coordinated - in compliance with the underlying physical weighing principle - to the existing acceleration due to gravity at its place of location (only if the weighing system has not already been adjusted to the location in the factory). This adjustment process must be carried out for the first commissioning, after each change of location as well as in case of fluctuating environment temperature. To receive accurate measuring values it is also recommended to adjust the display unit periodically in weighing operation.



- Prepare the required adjustment weight. The adjustment weight to be applied depends on the capacity of a weighing scale, see chap. 1. Carry out adjustment as closely as possible to admissible maximum load of weighing scale. Info about test weights can be found on the Internet at: http://www.kern-sohn.com.
- Observe stable environmental conditions.
 For warm-up time required for stabilisation see chap.1.




Access to service menu "x10" is locked in verified balances.


To disable the access lock, destroy the seal mark and actuate the adjustment switch. Position of the adjustment switch see chap. 18.1.

Attention:

After destruction of the seal the weighing system must be re-verified by an authorized agency and a new verification wire/seal mark fitted before it can be reused for applications subject to verification.

Procedure:






If errors during adjustment occur or an incorrect adjustment weight is used, an error message appears ("WronG"), Repeat adjustment process.

18.1 Adjustment switch and seals

Position adjustment switch and seal marks:



- 1. Adjustment switch
- 2. Self-destroying seal mark

TMPN A-BA-e-2213