

Abbott

Clearview™
hCG Cassette (25mIU/mL)

REF CV506788C

Clearview™ hCG Cassette (25 mIU/mL)

Intended Use

Clearview™ hCG Cassette (25 mIU/mL) is intended for the qualitative detection of human chorionic gonadotropin (hCG) in urine as an aid in the diagnosis of pregnancy. For professional *in vitro* diagnostic use only.

Introduction

hCG is a glycoprotein hormone produced by the blastocyst.^{1,2} hCG normally begins to be detected in the urine from 7 days after conception. The sudden rapid rise in concentration of hCG in urine following conception makes it an excellent marker for pregnancy.^{3,4}

Test Principle

Clearview™ hCG Cassette (25 mIU/mL) is a rapid chromatographic immunoassay for the qualitative detection of human chorionic gonadotropin in urine to aid in the early detection of pregnancy. The test uses two lines to indicate results. The test line utilizes a combination of antibodies including a monoclonal hCG antibody to selectively detect elevated levels of hCG. The control line is composed of goat polyclonal antibodies and colloidal gold particles. The assay is conducted by adding a urine sample to the sample well of the test device and observing the formation of colored lines. The sample migrates via capillary action along the membrane to react with the colored conjugate.

Positive samples react with the specific antibody-hCG-colored conjugate to form a colored line at the test line region of the membrane. Absence of this colored line suggests a negative result. To serve as a procedural control, a colored line will always appear in the control line region indicating that the correct volume of sample has been added and membrane wicking has occurred.

Kit Contents And Storage

Materials Provided

Each **Clearview™ hCG Cassette (25 mIU/mL)** kit contains sufficient materials for 20 tests: The device contains anti-hCG particles and anti-hCG coated on the membrane.

- 20 Pouched devices: Each sealed foil pouch contains 1 device and 1 desiccant packet
- 20 Disposable pipettes
- 1 Package insert

Store at 2-30°C. Do not use after the expiry date.

Materials Required But Not Provided

- Sample collection container
- Clock, timer or stopwatch

Precautions

1. Do not open the foil pouch until ready to test.
2. Do not use devices that have become wet, or if the foil pouch has been damaged.
3. Do not use kit beyond expiration date printed on the outside of the kit carton.
4. Do not reuse device.
5. Do not eat, drink or smoke in the area where the specimens or kits are handled.
6. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing.
7. It is recommended that disposable gloves should be worn while handling specimens.

8. All specimens and specimen contaminated materials should be disposed of in accordance with local biohazard waste disposal protocol.
9. To obtain accurate results, the package insert instructions must be followed.
10. It is recommended to immediately use the device after it is taken from the open pouch. If the testing is not done immediately after opening the device pouch, at extreme temperature and relative humidity condition the performance of the assay may decay or lead to invalid results.

Sample Collection And Storage

A urine sample collected at any time of the day is suitable, but a first morning sample is recommended.⁵ Urine samples must be collected in a clean and dry container. Samples may be stored in the refrigerator (2-8°C) for up to 48 hours or frozen below -20°C. Samples must be allowed to reach room temperature (15-30°C) prior to testing. Urine samples exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear sample for testing. **Clearview™ hCG Cassette (25 mIU/mL)** is not validated for use with samples containing preservatives.

Assay Procedure

Allow the device and urine sample to equilibrate to room temperature (15-30°C) prior to testing.

1. Remove the device from the sealed pouch and place on a clean and level surface.
2. Using the plastic pipette supplied, insert the tip into the sample and squeeze the top bulb fully (A). DO NOT squeeze the bottom bulb (B) when using the pipette. Release the top bulb to draw up the liquid. The sample (100 µL) will be drawn into the lower part of the pipette. Excess sample will be drawn into the bottom bulb.
3. Remove the pipette from the sample. Squeeze the top bulb fully to dispense the sample **carefully** onto the sample well (S) and start the timer. Use a new pipette for each test performed, even if using the same urine sample.

4. **Read the result at 3 minutes after dispensing the sample. Do not interpret results after the 3 minute read time.** It is important that the background is clear before the result is read.

Interpretation Of Results

(Please refer to the illustration)

- **POSITIVE: Two distinct colored lines appear.** One line should be in the control line region (C) and another line should be in the test line region (T). The color intensities of lines may vary. Therefore, any shade of color in the test line region should be considered positive.
- **NEGATIVE: One colored line appears in the control line region (C).** No apparent colored line appears in the test line region (T).
- **INVALID: Control line fails to appear.** Insufficient sample volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new device. If the problem persists, discontinue using the kit immediately and contact your local distributor.

Limitations

1. Positive results from very early pregnancy may later prove negative due to natural termination of the pregnancy.⁶ It is therefore recommended that weak positive results are re-tested 48-72 hours later with a first morning urine sample.
2. A negative result may be obtained with too dilute urine sample, as indicated by a low specific gravity, may not contain representative levels of hCG. If it is inconsistent with clinical evidence, results should also be confirmed by an alternative hCG method prior to medical procedure; or if pregnancy is still suspected, it is recommended the patient should be retested 48-72 hours later with a first morning sample.

- Concentrations of hCG are generally lower in ectopic pregnancy than expected normal values for a given gestational age. Abnormal pregnancy cannot be distinguished from normal pregnancy by hCG levels alone.^{2,7}
- hCG remains elevated for a time after pregnancy.⁸ Pregnancy tests carried out less than 3 weeks after giving birth or 9 weeks after natural loss or termination may need further evaluation.
- A number of conditions other than pregnancy can cause elevated levels of urinary hCG e.g. menopause, trophoblastic disease and certain non-trophoblastic neoplasms.⁹
- Drugs containing hCG may interfere with **Clearview™ hCG Cassette (25 mIU/mL)**, and produce misleading results.
- False positive and false negative pregnancy tests may be observed in patients with abnormal bladder or kidney function e.g. enterocystoplasties¹⁰ and renal failure.
- If the test result with **Clearview™ hCG Cassette (25 mIU/mL)** is not consistent with clinical evidence, further evaluation may be required.
- Concentrations of hCG greater than 500,000 mIU/mL may elicit a prozone effect.
- Clearview™ hCG Cassette (25 mIU/mL)** is not validated for use with samples containing preservatives.
- This test provides a presumptive diagnosis for pregnancy. A confirmed pregnancy diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.

Expected Values

hCG is not normally detected in the urine and serum specimens of healthy men and non-pregnant women. Healthy pregnant women have hCG present in their urine and serum specimens. The amount of hCG will vary greatly with gestational age and between individuals.

Clearview™ hCG Cassette (25 mIU/mL) has a sensitivity of 25 mIU/mL, and is capable of detecting pregnancy as early as 1 day after the first missed menses.

Performance Characteristics

Accuracy

A multi-center clinical evaluation was conducted comparing the results obtained using **Clearview™ hCG Cassette (25 mIU/mL)** and another commercially available urine membrane hCG test. The urine study included 159 samples, and both assays identified 88 negative and 71 positive results. The results demonstrated a >99% overall accuracy of **Clearview™ hCG Cassette (25 mIU/mL)** when compared to the other urine membrane hCG test.

Method		Other hCG Rapid Test		Total Results
Clearview™ hCG Cassette (25 mIU/mL)	Results	Positive	Negative	
	Positive	71	0	71
	Negative	0	88	88
Total Results		71	88	159

Sensitivity: >99.9% (95%-100%)*

Specificity: >99.9% (96%-100%)*

Accuracy: >99.9% (98%-100%)*

* 95% Confidence Intervals

Sensitivity and Specificity

Clearview™ hCG Cassette (25 mIU/mL) detects hCG at a concentration of 25 mIU/mL or greater. The test has been standardized to the W.H.O 4th International hCG Standard. The addition of LH (300 mIU/mL), FSH (1,000 mIU/mL), and TSH (1,000 µIU/mL) to negative (0 mIU/mL hCG) and positive (25 mIU/mL hCG) samples showed no cross-reactivity.

Prozone Effect

Clearview™ hCG Cassette (25 mIU/mL) has been shown to produce positive results with samples containing up to and including 500,000 mIU/mL hCG, which is higher than the maximum level expected during a typical pregnancy.

Interfering Substances

The following potentially interfering substances were added to hCG negative and positive samples.

Acetaminophen	20 mg/dL	Caffeine	20 mg/dL
Acetylsalicylic Acid	20 mg/dL	Gentisic Acid	20 mg/dL
Ascorbic Acid	20 mg/dL	Glucose	2 g/dL
Atropine	20 mg/dL	Hemoglobin	1 mg/dL
Bilirubin	2 mg/dL		

None of the substances at the concentration tested interfered in the assay.

Product Support

Contact one of the following Product Support Care Centers or your local distributor if you have any questions regarding the use of your product. You may also contact us at www.globalpointofcare.abbott.

Region:

EUROPE

Phone: +44 161 483 9032 Option 3
Email: EME.TechSupport@abbott.com

MIDDLE EAST

Phone: +965-2202 2828
Email: EME.TechSupport@abbott.com

APAC (Asia-Pacific)

Phone: +61 7 3363 7711
Email: AP.TechSupport@abbott.com

AFRICA

Phone: +27 10 500 9700 Option 3
Email: arcis.techsupport@abbott.com

RCIS (Russia and Commonwealth of Independent States)

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1. Hsu M., Kolm, P., Leete, J., Dong, K., Muasher, S., Oehninger, S. (1998). Analysis of Implantation in Assisted Reproduction through the use of Serial Human Chorionic Gonadotropin Measurements. *J. Assist. Reprod. Genet.* 15 (8): 496-503.
2. Catt K., Dufau M. & Vaitukaitis, J. (1975). Appearance of hCG in pregnancy plasma following the initiation of implantation of the blastocyst. *J. Clin. Endocrinol. Metab.* 40: 537-540.
3. Chard T. (1992). Pregnancy tests: review. *Hum. Reprod.* 7 (5): 701-710.
4. Corson S., Batzer, F. Schlaff, S. (1981). A Comparison of serial quantitative serum and urine tests in early pregnancy. *J. Reprod. Med.* 26 (12): 611-614.
5. Kaplan L.A. and Pesce A.J. (1989). The C. V. Mosby Company.
6. Wilcox A., et al. (1988). Incidence of early loss of pregnancy. *N. Engl. J. Med.* 319 (4): 189-194.
7. Braunstein G., Karow W., Gentry W., Rasor J., Wade M. (1978). First trimester chorionic gonadotropin measurements as an aid in the diagnosis of early pregnancy disorders. *Am. J. Obst. Gynec.* 131 (1): 25-32
8. Steier J. & Myking, O. (1984). Human Chorionic Gonadotropin in Maternal Plasma after Induced Abortion, Spontaneous Abortion, and Removed Ectopic Pregnancy. *Obstet. Gynecol.* 64 (3): 391-394.
9. Braunstein G., Vaitukaitis, J., Carbone, P., Ross, G. (1973). Ectopic Production of Human Chorionic Gonadotrophin by Neoplasms. *Ann. Intern. Med.* 78: 39-45.
10. Nethercliffe J., Trewick, A., Samuell, C., Leaver, R., Woodhouse, C. (2001). False-positive pregnancy tests in patients with enterocystoplasties. *BJU International.* 87 (9) : 780-782.

**Index of Symbols/Индекс на символите/Oversigt over symboler/Verzeichnis der Symbole/Πίνακας συμβόλων/Índice de símbolos/
Symboliluetelo/Index des symboles/Indice dei simboli/기호 색인/Index van symbolen/Symbolforklaring/Índice de símbolos/
Förteckning över symboler/Semboller Dizini**

 Consult instructions for use Консултирайте се с инструкциите за употреба Se bruksanvisningen Gebrauchsanweisung beachten Συμβουλευτείτε τις οδηγίες χρήσης Consulte las instrucciones de uso Lue käyttöohjeet Consultez le mode d'emploi Consultare le istruzioni per l'uso 사용 지침 참조 Raadpleeg gebruiksaanwijzing Se bruksanvisningen Consultar as instruções de utilização Se bruksanvisningen Kullanım talimatlarına başvurun	 Contains sufficient for <n> tests Съдържа достатъчно за <n> теста Indeholder tilstrækkeligt til <n> tests Ausreichend für <n> Prüfungen Περιέχει επαρκή ποσότητα για <n> εξετάσεις Contiene suficiente material para <n> pruebas Sisältö riittää <n> testiin Contenu suffisant pour <n> tests Contento suficiente per <n> test <n> 번 검사할 수 있음 Valdoende inhoud voor <n> tests Inneholder nok til <n> tester Contém o suficiente para <n> testes Innehållet räcker till <n> tester <n> test için yeterli miktarda içerir	 Manufacturer Производител Producent Hersteller Κατασκευαστής Fabricante Valmistaja Fabricant Produttore 제조사 Fabrikant Produsent Fabricante Tilverkare Üretici	 <i>In vitro</i> diagnostic medical device Μεταγεννητικό ιατροτεχνολογικό προϊόν Dispositivo para diagnóstico médico <i>in vitro</i> Lääkinnällinen laite <i>in vitro</i> -diagnostiikkaan Dispositif médical de diagnostic <i>in vitro</i> Dispositivo medico-diagnostico <i>in vitro</i> 체외진단용 의료기기 Medisch hulpmiddel voor <i>in vitro</i> diagnostiek <i>In vitro</i>-diagnostisk mediskinsk utstyr Dispositivo médico para diagnóstico <i>in vitro</i> Medicinteknisk produkt för <i>in vitro</i>-diagnostik <i>In vitro</i> tani amali tibbi cihaz
 Use-by date Дата на срок на годност Sidste anvendelsesdato Verwendbar bis Ημερομηνία λήξης Fecha de caducidad Viimeinen käyttöpäivämäärä Date d'expiration Data di scadenza 사용기한 Te gebruiken tot Brukes innen Prazo de validade Använd före-datum Son kullanna tarihi	 Do not reuse Не използвайте повторно Må ikke genbruges Nicht wiederverwenden Μην επαναχρησιμοποιείτε No reutilizar Et saa käyttää uudelleen Ne pas réutiliser Non riutilizzare 재사용 금지 Niet hergebruiken Må ikke brukes på nytt Não reutilizar Får inte återanvändas Yeniden kullanmayın	 Store between 2-30°C Съхранявайте между 2 – 30°C Opbevares ved 2 - 30 °C. Temperaturbegrenzung 2-30 °C Φυλάσσεται σε θερμοκρασία 2-30 °C Almacenar a entre 2 y 30 °C Säilytä 2 - 30 °C:ssa Conservation entre 2 et 30 °C Conservare a 2-30 °C 2-30°C에서 보관 Bewaren bij 2-30 °C Oppbevares mellom 2-30°C Conservar entre 2 °C e 30 °C Förvaras vid 2-30 °C 2-30°C'de saklayın	 Batch code Код на партида Batchcode Chargenbezeichnung Αριθμός παρτίδας Código do lote Eräkoodi Code de lot Codice lotto 배치 코드 Batchcode Batch- code Código do lote Batchcod Parti kodu
 Catalogue number Καταλογικη номер Katalognummer Artikelnummer Αριθμός καταλόγου Número de catálogo Tuotenumero Référence catalogue Número de catalogo 카탈로그 번호 Catalogusnummer Katalognummer Katalognummer Katalog numarası	 Authorized representative in the European Community Упълномощен представител в Европейската общност Autoriseret representant i EU Bevollmächtigter in der Europäischen Gemeinschaft Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα Representante autorizado en la Unión Europea Valtuutettu edustaja Euroopan yhteisössä Représentante autorisé dans la Communauté européenne Rappresentante autorizzato nell'Unione europea 유럽 공동체의 공식 대리인 In de Europese Gemeenschap gevestigde gemachtigde Autoriseret representant i EU Representante autorizado na Comunidade Europeia Auktoriserad representant i EU Avrupa Topluluğunda yetkili temsilci		

V02

Revision Date: 2023-01-13

© 2023 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

www.globalpointofcare-abbott

Hangzhou, 310018, P.R.China

Hangzhou Economic & Technological Development Area,

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Abon Biopharm (Hangzhou) Co., Ltd.
MT Promed Consulting GmbH