Manuale d'uso

User's Manual

Gebrauchshandbuch




Manual de uso

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




Ⓓ


I

II

III

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Welcome to the Stet Cube system for stethoscope head disinfection, by means of UV-C LED.

Please read this user's manual carefully in order to use the device safely and achieve its best performance.

The product was developed by egoHEALTH S.r.l. (www.egohealth.it) in collaboration with Light Progress S.r.l. (www.lightprogress.it).

For more details, please visit the website: www.stetclean.com

1. INTENDED USE

Stet Cube is a Class I medical device specifically designed for the disinfection of stethoscopes.

2. RECEIPT

Handle with care. Fragile material! Verify that the packaging is complete and intact and be careful not to damage the contents when using sharp blades or cutters to open the box. Take out the device from the box and immediately check if it has not been damaged during transport.

Before using this device, carefully read the safety WARNINGS and all other instructions described below.

Once out of the box, we recommend connecting Stet Cube to the power supply to recharge the battery, using the cable included in the box. The first time you charge the battery, we recommend leaving Stet Cube plugged in for at least 12 hours. In case of anomalies in operation, press the Reset button (6) (see section 4.4).

3. OVERVIEW

The disinfection process based on UV rays with wavelength 260-290 nm (UV-C) is activated a few seconds after the closure of Stet Cube, only if the head of a stethoscope has been placed inside.

Stet Cube is disabled automatically at the end of the disinfection time set or during operation, if the Stet Cube lid is opened.

The UV-C LED turns on when the Stet Cube lid is closed, which activates a micro-switch positioned inside the base via mechanical pressure applied by a lever (8) on the lid hinge, and when the head of the stethoscope is detected via an infra-red proximity sensor (3) as illustrated in figure C. This double check ensures absolute safety and deep hygiene to operator and patient, as it contains UV-C in a closed environment and enables the device only in the presence of the head of the stethoscope. The

transparent housings in the base (2) (4), and the magnets on the inside of the lid (7) ensure the correct positioning and securing of the stethoscope within the device whether it is used in a horizontal or vertical position (wall mounted).

4. USE AND OPERATION

The device features mechanisms designed to ensure simple and immediate operation.

4.1. Insertion of the stethoscope

Disinfection is enabled by inserting the head of the stethoscope into Stet Cube and then closing the lid, which activates the micro-switch (8) and the infra-red proximity sensor (3) (Figs. C/D, phases I → II → III).

4.2. Membrane disinfection

Remove any dirt from the membrane before using Stet Cube. To make the action of UV-C more effective, it is recommended to dry the membrane of the stethoscope, if wet. When the stethoscope is correctly placed, the blue light (10) indicates that the disinfection process was activated correctly, and it will flash throughout the whole period T1 of standard disinfection (see point 5). The UV-C LED is activated a few seconds after the blue light starts flashing to indicate that the device was placed properly. For efficient disinfection, the stethoscope should be placed within the transparent housings (2) (4).

4.3 Disconnecting the stethoscope

At the end of phase T2 (deep disinfection), the blue light (10) turns off at the same time as the UV-C LED (1). The stethoscope will be deeply disinfected and can be removed from Stet Cube, ready for use (Figs. C/D, reverse phases III → II → I).

The disinfection process, both standard and deep, may be discontinued at any time by removing the stethoscope from Stet Cube.

To ensure effective disinfection, we recommend that you remove the stethoscope under the conditions described in point 4.2.

4.4. Battery charging and reset button

The battery can be recharged through a micro-USB port (5). A yellow light (9) flashes during the charging process, remains steadily lit when the battery is fully charged and turns off once you disconnect the cable.

It is advisable to recharge the battery at least once a week in order to not compromise charge capacity.

A reset button (6), which can be operated via a fine tip, restores the standard conditions of the device.

4.5. Wall mounting

Stet Cube can be attached to a vertical wall with the wall bracket provided (11). For wall mounting, it is first of all necessary to fix the bracket to the wall at the correct height, using the screws provided.

Stet Cube can then be fitted to the bracket via the housings underneath the base (12). The system can be removed from the bracket at any time by sliding Stet Cube upwards. If fitted near to an electricity outlet, Stet Cube can also be recharged directly from its position on the wall. "Vertical" disinfection operations are described in figure D (phases I → II → III). It is recommended to always place the head of the stethoscope within the transparent housings (2) (4) to ensure it is held firmly in place even with the lid open.

5. OPERATING RANGES

Stet Cube ensures effective disinfection throughout the shelf life of the UV-C LED (more than 50.000 uses).

During hygienization phases, the system defines automatically a T1 step (standard), lasting 3 minutes + T2 step, adding an extra 2 minutes to disinfect stethoscope more deeply.

6. LIGHTSIGNALS

The operating lights indicate:

-Blue light = UV-C LED operation

-YELLOW light = battery operation

6.1 Summary table of light signals

UV-C LED Control		BATTERY control	
T1	T2		
		 Charging	 Battery charged
		 (for 10 s)	Battery charge level
		 (for 10 s)	Battery empty
	 (upon plug-in)	Malfunction	
	 (upon plug-in)	UV-C LED out	

 Flashing - n flashes

 Fixed light

7. LIMITATIONS

Stet Cube can be used with almost all stethoscopes.

The incomplete closure of the lid may, however, prevent the device from working

8. IMPORTANT WARNINGS

We recommend you follow the instructions below.

⚠ Keep away from heat, sunlight, moisture and water.

⚠ It is recommended to dry the membrane of the stethoscope before inserting it into Stet Cube, to avoid the formation of condensate and ensure optimal disinfection.

⚠ Keep the Stet Cube away from sources of electromagnetic waves and ionising radiation (MRI, x-ray, microwave, etc.).

⚠ Do not tamper with or force the closure mechanism, to avoid activating the UV-C LED with the lid open and/or the stethoscope placed incorrectly.

⚠ Should the UV-C LED be activated when the lid is open, never direct the UV-C LED towards skin and eyes.

9. DEVICE CLEANING

Preferably use a soft dry cloth. If necessary, clean organic residues or stubborn dirt with denatured alcohol.

10. DISPOSAL

Dispose of the device properly, keeping in mind that it contains a lithium battery.

10.1. Main reference standards

- EC Directive 93/42/EEC as amended (2007/47/EC)
- Law Decree 46 of 1997
- CEI 14971:2012
- CEI 60601-1:2010
- CEI 60601-1-2:2007

11. TECHNICAL SPECIFICATIONS

11.1. General features

Power supply:

rechargeable lithium-polymer battery;
3,7V - 1000mAh voltage $5V \pm 5\%$,
non replaceable.

Consumption UV LED: 120 mA

Battery autonomy: 150 consecutive uses,
5 days in standby mode.

Input connector: micro-USB port (micro-B).

11.2. Environmental conditions

Operating temperature: from 0°C to $+40^{\circ}\text{C}$.


Relative operating humidity: $< 90\% \text{ RH}$.

Storage temperature: from -20°C to $+50^{\circ}\text{C}$.

Relative storage humidity: $< 90\% \text{ RH}$.

12. SYMBOLS

	CE Mark
	Manufacturer
	Device Disposal
	See Instruction Manual
	Serial Number
	Humidity Limitation
	UV rays source
	do not wet

This document is authentic in the Italian language.

Licensors :
egoHEALTH S.r.l.
Via Fiorentina, 1,
53100 Siena (Italia)
www.egohealth.it

Licensee:
Light Progress S.r.l.,
Loc. San Lorenzo, 40
52031 Anghiari (AR) Italia
(+39) 0575 749255
www.lightprogress.com
www.stetclean.com
P.IVA IT 01208950517

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We herewith declare that

- *The above-mentioned device is manufactured in compliance with the requirements laid down in Directive 93/42/EEC as amended in **Annex VII**;*
- *The above-mentioned device is to be considered as belonging to **class I** in accordance with rule 12 of Annex IX of the said Directive;*
- *The entire documentation concerning this device is enclosed in the Technical Dossier archived at the headquarters of Light Progress Srl and shall be stored for a period of at least 5 years from the last date of manufacture of the product;*
- *Light Progress S.r.l. has notified the Competent Body that the said device was put on the market, in order to ensure after-sales monitoring;*

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This statement is supported by the registration and subscription request submitted by Light Progress Srl to the Ministry of Health on July 17th, 2018 (No. 1725675) as: "manufacturer of class I medical devices"


Anghiari, lì 17/07/2018

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Ministry of Health registration Nr. :	N°1725675
CND:	Z12011399
GMDN:	62171
FDA:	
Establishment Registration :	3010276038
DUNS No.:	43-566-3059
Device Classification Name:	DISINFECTOR, MEDICAL DEVICES
Product Code:	MEC
Regulation Number:	880.6992
Device Listing Nr:	D295706
UDI:	(01)08052780880312 (21)(...)
Legal Entity GLN:	8052780880008
GS1 company prefix	805278088

www.stetclean.com

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FABBRICANTE / MANUFACTURER: : Light Progress Srl – Loc. San Lorenzo, 40 520310 Anghiari (AR)