EN

DE

FR

IT

ES

HU

bloom ferritin test

ferritin

instructions for use

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REF FER-REF-2-2.X; (5 tests)

Table of contents

1. General information	3
1.1. Intended use.....	3
1.2. Terms and definitions	3
2. Bloom System	3
2.1. Overview of the Bloom System	4
2.2. Bloom Ferritin Test	4
2.2.1. Principle of the method	4
2.2.2. Test performance	4
2.2.3. Cross reactivity and Interferents.....	4
2.2.4. Components.....	4
3. Safety Regulations.....	4
3.1. General notes on safety	4
3.2. Limitations of use and contraindications.....	5
4. Operating steps	5
4.1. Preparation.....	5
4.2. Blood sample collection	5
4.3. Process sampling and report collection.....	7
4.4. Device cleaning.....	8
5. Storage	8
6. Troubleshooting	8
7. Symbol information	8
8. Contact information.....	8

1. General information

1.1. Intended use

The Bloom Ferritin Test is a single-use self-testing kit for quantitative in-vitro diagnostic examination of ferritin in capillary blood, using lateral flow technology. It is intended for people aged above 18, interested in tracking and optimization of health and non-acute medical conditions. To be used in a home environment or store setting.

1.2. Terms and definitions

- **Warning:** Points out a situation where a potentially hazardous condition exists that could result in a critical or catastrophic situation such as death or severe long-term injury if not avoided.
- **Caution:** Points out a situation where a potentially hazardous

condition exists that could result in a serious situation such as short-term injury or impairment requiring medical intervention to correct if not avoided.

- **Information:** Points out a situation where a potentially hazardous condition exists that could result in a not serious situation, such as a slight inconvenience, if not avoided.

2. Bloom System


The Bloom Ferritin Test works only in conjunction with the Bloom App and Bloom Lab, collectively known as the "Bloom System". To get started, download and install the Bloom App on the smartphone and ensure you have access to a Bloom Lab.

Warning

Do not use the Bloom System as a diagnostic device. Using the information provided by the Bloom System for diagnosis can lead to mistreatment and severe consequences. Results should be reviewed by a medical practitioner.

2.1. Overview of the Bloom System

The Bloom System is a medical testing system. It provides screening and monitoring of target markers and medical guidance. It combines lateral flow assay-based bodily fluid testing (Bloom Test strips), an established diagnostics technology, with an electronic reader for test strip analysis (Bloom Lab) and user-friendly software for displaying and interpreting results (the Bloom App).

2.2. Bloom Ferritin Test

2.2.1. Principle of the method

The Bloom Ferritin Test strip is a rapid immunochromatographic assay based on gold nanoparticles. When the blood is mixed with the provided dilution buffer and applied onto the test strip, the diluted blood flows along the test strip and binds to anti-Ferritin antibodies, which detect the ferritin in the blood and will create a color line detected by the Bloom Lab. The Bloom Ferritin Test strip also includes a control line, which will be detected by the Bloom Lab to verify whether the test has run correctly. Users will receive their results in the form of a report via the Bloom App.

2.2.2. Test performance

Reference Test			Sensitivity	100.0 %
	Positive	Negative		
Bloom	Positive	65	2	67
Ferritin	Negative	0	33	33
Test	Total	65	35	100
			Specificity	94.3%
			Accuracy	98.2%

The Bloom Ferritin Test has a quantifiable range between 10 - 500 µg/L.

2.2.3. Cross reactivity and Interferents

Unusually high titers of heterophile antibodies or rheumatoid factor (RF) in the blood may affect the test results.

2.2.4. Components

The Bloom Ferritin Test includes one set of instructions for use and enough material for either one or five tests (indicated on the outer package). Material for one test comprises of the following:

1		1 x Bloom Ferritin Test strip in aluminium pouch		
2		1 x Alcohol wipe for finger cleaning		Servoprax GmbH, Wesel, Germany
3		1 x Sterile safety lancet for finger pricking		Medipurpose Pte. Ltd., Singapore, Singapore STERILE R
4		1 x Plaster for incision patching		BSN Medical GmbH, Hamburg, Germany
5		1 x Capillary tube (lithium heparin) for blood collection		Sarstedt AG & Co., Numbrecht, Germany
6		1 x Buffer		
7		1x Disinfecting cleaning wipe for Bloom Lab		PDI Ltd., Flint, United Kingdom

3. Safety Regulations

3.1. General notes on safety

1. Do not use this Bloom Test if you are intoxicated or find yourself in an acute situation of emergency or distress. In case of doubt, always consult with your medical practitioner.
2. The Bloom Report displays the result of the Bloom Test and, together with the answers given in the questionnaire, offers additional information. It is not intended to diagnose or inform about the treatment of medical conditions. Never avoid, put off, change or start treatment based on the content of the Bloom Report. In order to receive a diagnosis and treatment, consult a medical practitioner. In case of health concerns, seek medical attention immediately.
3. With the exception of the measured test value, all information presented in the Bloom Report is the result of Bloom Diagnostics' research conducted with literature available at the time of development. While the information in the Bloom Report is provided with great care, all content provided is subject to change depending on scientific discoveries and advances, and might not apply globally.

4. The validity of the Bloom Report relies on the correctness of your answers. Once submitted, the answers cannot be changed.
5. Do not perform multiple Bloom Tests at the same time. Do not use the Bloom Test in an environment with high humidity, such as the bathroom.
6. Do not use the Bloom Test, if its packaging or any of its components are damaged or show signs of tampering, malfunctioning, being broken, defective, or contaminated.
7. Do not attempt to modify, deconstruct, damage, or tamper with any elements of the Bloom System as this can result in malfunction and will void the support service and product's warranty.
8. The Bloom System should not be used with any hardware or software that is not manufactured or developed by Bloom Diagnostics for the explicit purpose of being used with the Bloom System. This includes using generic power cables and uncertified or generic test strips.
9. The Bloom Lab and entire Bloom System should only be purchased from authorized retailers and never second-hand.
10. This device may be used by people with reduced mental capabilities provided that they are supervised and are fully aware of the subsequent risks of use.

3.2. Limitations of use and contraindications

1. Do not take this test if you are under 18 years of age, if you are, or may be, pregnant.
2. The test result can be influenced by physiological or pathological factors, which might compromise the interpretation of the test result.
3. Ferritin levels can be affected by inflammation, such as a present infection or an acute event, such as a recent heart attack or surgery. Liver disease and chronic conditions, such as autoimmune disease, can influence the result and might compromise its interpretation.
4. Other factors such as age, menopause status, diet and medication/supplementation intake can influence your test result.

4. Operating steps

Always follow the Bloom System's operating steps in a sequential order.

Caution: When using a Bloom Test containing five tests:

- Keep the test packaging until you have performed all five tests, in order to initiate a test.
- Ensure the test kit is only handled by one individual.

4.1. Preparation


Before you start, ensure that you have:

- checked the use-by-date has not been exceeded (located on the package in the format YYYY-MM-DD), Do not use an expired Bloom Test;
- a smartphone with a working internet connection and the Bloom App installed;

- a working Bloom Lab, set up according to its instructions for use;
- checked the Bloom Ferritin Test package components for completeness and damage;
- a clean, stable workspace and are seated while carrying out the test.

Step 1: Connect Bloom App and Bloom Lab


1. Log in to the Bloom App and tap on the "+" button in the menu bar to initiate a new test.
2. Swipe the phone across the light ring of the Bloom Lab to connect.

Information: Successful pairing is indicated by a white heartbeat pulse on the Bloom Lab.

Step 2: Register the Bloom Test

1. Open the Bloom Test packaging.
2. Swipe the phone across the lid of the Bloom Test packaging to register. Or register manually, by entering its serial number, which is indicated on the bottom of the packaging with "SN".

4.2. Blood sample collection


Caution: All blood samples should be treated as potentially infectious. Prevent the contamination of surfaces by performing the blood sample collection over a removable surface, such as a paper sheet or the opened Bloom Test packaging.

Warning
The accuracy of the test result heavily relies on the correct collection of the blood sample. After registering the Bloom Ferritin Test, follow all instructions with utmost care.

Step 1: Select and disinfect the finger

1. Warm up hands by washing them with warm water or rubbing them together. This will improve the blood flow.
2. Select and disinfect the middle or ring finger using the provided alcohol wipe.

Disinfect the finger.


Caution: Do not perform the blood sample collection in one of the following situations:

- On an arm where an infusion has been performed within the last 24 hours.
- On a finger of the side of the body where a mastectomy has been performed.
- On a swollen, infected, wounded or previously punctured finger.
- If you have a history of fainting at the sight of blood.
- If you have severely impaired dexterity, or lack of fine motoric skills.
- Do not use any object other than the sterile lancets provided in the Bloom Ferritin Test to prick the skin.
- The lancet should only be used on clean and dry middle or ring fingers. Skin should be healthy, without any dermatological issues or wounds. Never prick a finger twice.

Step 2: Use the sterile lancet on the disinfected fingertip


1. Remove the protective cap of the sterile lancet.


2. Press the red end of the lancet on the side of the cleaned fingertip to trigger it.

Step 3: Collect blood sample immediately after the puncture is created


1. *Massage hand from palm to finger to create a large blood drop.*
2. *Let the capillary tube soak up the blood drop.*


3. Fill the capillary tube up to the marked line.

Step 4: Apply the plaster to the puncture on the fingertip using light pressure to prevent ongoing bleeding

Step 5: Add the dilution buffer to the capillary tube immediately


1. Open the dilution buffer dispenser by twisting the lid.

2. Fill the capillary tube completely with the dilution buffer up to the white stopper.


Step 6: Mix the sample

Gently flip the capillary tube two times to mix your blood sample with the dilution buffer.

Step 7: Transfer the collected blood onto the Bloom Test strip immediately


1. Open the Bloom Test strip pouch.
2. Place the Bloom Test strip on a flat surface.

3. Empty the filled capillary tube onto the sample window of the Bloom Test strip by pushing down the colored piston. Ensure that the capillary tube is completely emptied onto the test strip

Step 8: Let sample sink in


1. Tap the "Confirm" button in the Bloom App.
2. Wait for absorption (see Bloom App timer).

Caution: Examination results can be affected if the Bloom Test strip is moved before the liquid sample has been absorbed, dropped or not handled as intended.

Step 9: Insert Bloom Test strip into Bloom Lab

1. Insert the Bloom Test strip into the Bloom Lab. You should feel the Bloom Test strip lock into place. The Bloom Lab will show a white heartbeat pulse on the LED ring when the test strip is inserted correctly.

Insert Bloom Test strip into the Bloom Lab.

2. Processing begins. The progress of the measurement is communicated by the LED ring of the Bloom Lab with a spinning white loader. The Bloom App shows the estimated time remaining.
3. Once processing has been completed, the Bloom App indicates the successful completion of the testing.

Caution:

- Place the Bloom Lab on a stable, flat surface and do not move the device during measurement.
- Do not remove the Bloom Test strip from the Bloom Lab during the testing procedure.
- Ensure that the Bloom Lab does not come into contact with blood.
- Make sure that the Bloom Test strip is clean and free from dust before inserting it into the Bloom Lab.

Possible errors:

- In case of a recoverable error, the LED ring of the Bloom Lab will show a red pattern, indicating that something went wrong.
- In case of a fatal error, such as a general error requiring service of the Bloom Lab or power adapter, the LED ring of the Bloom Lab will show a solid red ring, indicating that a general hardware failure error occurred.

Step 10: Dispose components

The used Bloom Ferritin Test components are potentially biohazardous. Please dispose of them carefully according to applicable laws and regulations.

4.3. Process sampling and report collection

The final step of the testing process is to complete the questions needed to compile the personalized Bloom Ferritin Report using the Bloom App. While the Bloom Lab is processing the test result, the Bloom App will prompt you to do the following:

Check the Bloom App

1. Answer the questions. When you are presented with a multiple-choice question, choose all answers that apply and then confirm to proceed.

2. After around 10 minutes, the results from the Bloom Lab will be transmitted to the Bloom App. The app will alert you that the Bloom Ferritin Report is ready.
3. View the personalized Bloom Ferritin Report.

Information: If the Bloom App disconnects from the Bloom Lab during the measurement process, the Bloom Lab completes the measurement process and saves the results. If the same Bloom App reconnects with the Bloom Lab before another measurement is initiated, the saved measurement value will be sent to the Bloom App and the Bloom Lab will subsequently delete the measurement value and any unique identifier of the Bloom App. After a 10-minute period other Bloom Apps will be able to connect to the Bloom Lab and perform new tests, which will automatically delete the previous results.

4.4. Device cleaning

Caution: Please use protective equipment, such as gloves, when cleaning the surfaces of the Bloom Lab to avoid contamination.

After performing the Bloom Ferritin Test, please make sure to leave a clean workspace. Thoroughly wipe all surfaces of the Bloom Lab after every use with the provided disinfection cleaning wipe. Please check gaps for stains and clean thoroughly when needed.

5. Storage

In order to avoid hazardous situations, the Bloom Ferritin Test should be stored with care at all times, this includes:

- The recommended storage temperature between +8°C and +30°C (46°F and 86°F)
- A dry indoors environment without direct or excessive heat and humidity below 80% (bathrooms not suitable)
- Out of reach for children to avoid injury through misuse of components like lancets.

6. Troubleshooting

For further information about the Bloom Lab, Bloom App, frequently asked questions and troubleshooting, please visit www.bloomdiagnostics.com/faq or report the error using your Bloom App (Settings - Contact us).

7. Symbol information

CE Mark

Manufacturer

Use-by date

Consult instructions for use

Fragile, handle with care

Serial number

Keep dry

Keep away from sunlight

Sterilized using irradiation

Do not use if package is damaged

Single use

Batch number

Contains sufficient material for one test

Contains sufficient material for five tests

Catalog number

Recycling

Humidity exposure limit
80%
0%

Temperature exposure limits between +8°C and +30°C (46°F and 86°F)

In vitro diagnostic medical device

Authorized Representative

Importer

8. Contact information

For additional information regarding technical support and local distributors, contact:

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Bloom Diagnostics AG
Bloom Diagnostics GmbH