
Lûx

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
User Manual

My intended use

I, Lûx, am only intended for use by healthcare professionals for medical diagnoses. I am a diagnostic light for non-contact examination of pupillary responses and for illuminating bodily cavities.

My built-in diagnostic mode makes it easier for you to evaluate pupillary reflexes thanks to a control sound after 0.5 seconds.

Please observe a minimum distance of 100 mm when examining pupillary reflexes and note that this is a temporary examination of a maximum of 58 s.

Inappropriate use

If you use me in a manner other than described in this user manual, this does not correspond to my purpose. My manufacturer is not liable for damage resulting from improper use. In such cases, you, the user, bears the risk alone.

Please note that I am not waterproof, so do not test my swimming ability.

Operation

To turn the light on, hold down the button **1** at my rear end. To turn it off, simply release the button.

To activate the diagnostic mode, slide the switch **2** down (towards the button).

Should my battery be low, use the supplied USB-C cable **3** and connect me to a power source. A small blue light indicates the charging process. When the light turns off, I'm fully charged and ready to go again.

Safety notice

I am not meant for invasive use.

Because of swallowable small parts (LEDs, USB-C cover), you should not put me in the hands of children under 3 years old!

Check that I am functioning properly before each use.

If you find any damage on me, you should take me out of use immediately.

If I am not working, have the repair carried out by qualified specialist personnel only.

Likewise, I should not be used in areas where there is a risk of fire or explosion, or in the vicinity of strong magnetic fields.

Please note that my light may be harmful. The risk of eye damage increases as the duration of irradiation increases.

ISO 15004-2: Group 2 (LED)

Caution – The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater is the risk of ocular damage. An exposure time with this instrument of longer than 2 mins 50 secs per eye per day at 100mm distance will lead to the guideline hazard value being exceeded.

Please be aware that using the instrument at a distance closer to the eye than 100mm will result in a reduction of the maximum daily exposure time.

ANSI Z80.36-2021: Group 2 (LED)

Caution – The light emitted from this instrument is potentially hazardous. The longer the duration of exposure, the greater is the risk of ocular damage. Exposure to light from this instrument will exceed the recommended maximum exposure (RME) of 2.2 J/cm², unless additional action is taken by the user to minimize exposure, after 40 seconds.

The risk of retinal injury at an exposure of 2.2 J/cm² is not high, but because some patients may be more susceptible than others, caution is

advised if this radiant exposure value is exceeded. However, because of a significant risk of injury at exposures exceeding 10 J/cm², the user should avoid exposures longer than 3 min and 5 Seconds.

Since prolonged exposure to intense light can cause eye damage, the use of the device for eye examinations should not be prolonged unnecessarily.

Infants and people with aphakia or diseased eyes are at greater risk of eye damage. The risk may also be increased if the person examined has already had contact with the same instrument or another ophthalmic instrument with a visible light source within the last 24 hours.

I have been tested according to the standard IEC 62471:2006 and am photobiologically compatible!

Materials

My elegant case is made of aluminium and held together with stainless-steel screws. You can attach me to your coat pocket using my stainless-steel clip. Inside is my circuit board as well as an ICR 10440 lithium battery. My button and the toggle switch are made of ABS.

Depending on how brightly you want me to shine, you can use me with three LEDs of different brightnesses. The luminosity ranges from 5,000–8,000, 8,000–10,000, and up to 10,000–14,000 mcd. The cover for the USB-C port is made of silicone.

Changing components

Replacing the LED

Am I shining too brightly or too weakly for you? Then simply change my LED!

You can also easily adjust my light radius. Two different

positions are provided for this in my housing where the LED will find a safe place.

Replacing the battery

If my battery stops working, you can easily change it yourself.

Procedure when changing components


Place me on a firm surface. Use the hex key provided to loosen the first screw. Now you can take off the clip and remove the second screw.

Take off the top half of my housing.

To get to the circuit board, you will need to remove the button and the indicator light.


Now remove the circuit board from the housing.


If you want to replace my LED, pull it out of the circuit board and insert the new one the same way. Please note that the longer contact is inserted into the opening on the right, as shown in the illustration.

The luminosity of the LEDs is indicated by a colour marker:


When inserting the circuit board, you can choose between two positions for the LED to vary my luminous radius: Pos. A = large cone of light, Pos. B = small cone of light.

To change the battery, take the plug out of the connector and plug in a new battery using the same principle.

The connection only fits in one way, so you can't go wrong when plugging in the new battery.

However, make sure that the battery cable is placed around the screw opening, otherwise you won't be able to close my housing.

After replacing the desired parts, reinsert the circuit board, replace the indicator light and button, and close the housing.

Finished! I'm ready to go.

Cleaning, disinfection and storage

Manual cleaning

Please clean me regularly. It's very simple: wipe me thoroughly with a damp cloth.

My manufacturer recommends cleaning wipes with 70% isopropyl alcohol or soapy water.

If you use other cleaning or disinfecting agents, please test their compatibility with the materials first. Otherwise, my pleasant outer appearance may get damaged. Frequent cleaning doesn't bother me. I like to be clean.

Sterilisation

You CANNOT sterilise me. If you do, I will break and probably come to a sad end.

Storage

I'm an undemanding guy. It is best to store me dry, in my packaging, and at room temperature. I have nothing against cool rooms either.

Maintenance and warranty

If my LEDs, screws or battery are defective, you can replace them.

The warranty for the entire product will be voided if replacement parts are used that were not procured from the DocCheck Shop.

If you use me as intended and follow this instruction manual, the warranty for manufacturing and material defects is two years starting from the date of purchase. Obvious abuse and accidental damage are excluded. Cosmetic damage, such as scratches on the surface, are also excluded from the warranty.

Check and functional test

Please carry out a visual inspection for damage and wear before each use.

Disposal

Please note that I have to be sent to a separate collection point for electrical and electronic equipment. If you have to dispose of the battery, this must be done at appropriate collection points. You should take a look at the specific disposal rules in your country for this.

Note to the user and/or patient:

All serious incidents related to me must be reported to DocCheck Shop GmbH and the competent authority of the Member State.

Manufacturer

DocCheck Shop GmbH
Vogelsanger Str. 66, 50823 Cologne, Germany
[doccheckshop.com](https://www.doccheckshop.com)

Explanation of symbols

Medical device

CE marking

Manufacturer data (address)

Production date

Item number

Batch

Contains no latex/natural rubber

Separate disposal of electrical and electronic appliances.


Separate disposal of used batteries

My hero
[doccheck.com](https://www.doccheck.com)

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