# RE-CROSS

Cross-linked Sodium Hyaluronate

by  
Reviscon

Ⓧ	Gebrauchsanweisung	2-5
Ⓧ	Instructions for use	6-9
Ⓧ	Instrucciones de uso	10-13
Ⓧ	Notice d'Utilisation	14-17
Ⓧ	Istruzioni per l'uso	18-21
Ⓧ	Instructies voor gebruik	22-25
Ⓧ	Instruções de Uso	26-29
Ⓧ	Инструкция по использованию	30-33
Ⓧ	Kullanma Talimatı	34-37

CE 1984

**Produkt:**

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**Product:**

Re-Cross by Reviscon is lightly cross-linked sodium hyaluronate viscoelastic solution for single injection providing longer durability intraarticular use.

Re-Cross by Reviscon consists of 4.0 mL of viscoelastic solution in a 5.0 mL single use syringe.

**Description:**

Re-Cross by Reviscon is non-animal structured intraarticular injection which is sterile, non-pyrogenic, clear, non-inflammatory lightly cross-linked sodium hyaluronate dissolved in buffered physiological saline solution. The sodium hyaluronate is obtained from by bacterial fermentation and it is cross-linked with a very low amount of crosslinker. Thanks to the unique purification in the production process, it provides highly purified product which allows safety to the patients.

Re-Cross by Reviscon contains lightly crosslinked sodium hyaluronate, sodium chloride, disodium hydrogen phosphate, sodium dihydrogen phosphate and water for injection.

	<b>Re-Cross by Reviscon</b>
1.0 ml contains	22 mg sodium hyaluronate
pH	6.8-7.6

**Properties and efficacy:**

All synovial joints especially the weight – bearing joints, contain viscoelastic sodium hyaluronate. This substance has lubrication and shock absorbing properties, allowing these joints to move normally and painlessly. In patients with degenerative joint disease (osteoarthritis), the viscoelasticity of the synovial fluid is significantly impaired, causing the mechanical stress on the joint and the breakdown of the articular cartilage to greatly increase resulting in limited and painful joint movement. Intraarticular administration of high purity sodium hyaluronate, which has very good viscoelastic properties, can improve the quality of the joint's lubrication. The lubrication and shock absorbing properties of this product reduce pain and improve joint mobility. This effect last at least for 6 months after the first injection.

**Indications:**

Re-Cross by Reviscon is indicated as a viscoelastic supplement or a replacement for synovial fluid in human joints. Re-Cross by Reviscon is well suited for treatment of the symptoms of human joint dysfunctions such as osteoarthritis;

- non-surgical treatment options for knee and hip osteoarthritis (OA)
- treatment of symptoms associated with Temporomandibular Joint Dysfunctions (TMJ).
- for the treatment of chronic shoulder pain associated with osteoarthritis or rotator cuff damage.
- intra- and postoperative intraarticular injection in patients undergoing wrist arthroscopy.

The actions of Re-Cross by Reviscon are lubrication and mechanical support.

**Contraindications:**

Do not administer to patients with known hypersensitivity to any of the materials contained in Re-Cross by Reviscon.

Intraarticular injections are contraindicated in cases of present infections or skin diseases in the area of the injection site to reduce the potential for developing septic arthritis.

**Side Effects:**

Following the use of Re-Cross by Reviscon, patients may experience local symptoms in the joint being treated (pain, sensation of heat, reddening and swelling). The following adverse events have been reported for similar products: mild or moderate arthralgia, in rare cases skin rash, aseptic joint effusions, pruritus and muscular cramps. Further adverse events that have been observed in very rare cases are: allergic reactions, anaphylactic shock, hemarthrosis, phlebitis, pseudosepsis, severe acute inflammatory reaction (SAIR), nasopharyngitis, joint stiffness, tendonitis, bursitis, fever and myalgia.

**Warnings and Precautions:**

•Re-Cross by Reviscon is intended for single use only. The reuse of the product creates a potential infection risks for patients or users.

•Sodium hyaluronate is manufactured by bacterial fermentation and rigorously purified. However the physician should consider the immunological and other potential risks that can be associated with the injection of any biological material.

•Do not reuse syringe or needle. Any repeat usage of the syringe or needle carry a risk of contamination and infection of the patient.

•Do not re-sterilize the pre-filled syringe. Performance will be impaired. The contents of the syringe are sterile but the exterior of the syringe is not sterile.

•Do not use if product inner packaging is opened or damaged.

•Do not use after the expiry date printed on the pack.

•Dispose of the syringe and needle in accordance with accepted medical practice and applicable national, local and institutional requirements.

•There is no evidence concerning the safety of Re-Cross by Reviscon in human pregnancy, lactation and the children under 18 years of age. Administration during pregnancy and lactation is at the discretion of the physician.

•Temporary pain or swelling may occur after intra-articular (IA) injection.

•As with any invasive joint procedure, patients are advised to avoid strenuous or prolonged (i.e. more than an hour) weight-bearing activities such as running or tennis within 48 hours following intra-articular injection.

**Interaction with other agents:**

Sodium hyaluronate is incompatible with quaternary ammonium compounds such as benzalkonium chloride solutions. Therefore, Re-Cross by Reviscon pre-filled syringes must never come into contact with surgical instruments rinsed with these solutions.

**Directions for use:**

**FOR INTRAARTICULAR INJECTION. FOR SINGLE USE ONLY.**

Only medical professionals trained in accepted injection techniques for delivering agents to joint spaces should inject sodium hyaluronate for this application. The required amount of product is injected into the joint space to be treated by using appropriate size and length of the needle. Common needle gauges for injections into the knee are 18-22 gauge. The physician decides the final needle selection for all indications. In order to avoid intraarticular infection strict aseptic injection technique has to be applied. It is recommended that an ice-pack be placed on the joint undergoing treatment for 5 -10 minutes in order to prevent pain and swelling. In the case of effusion accompanied by severe pain the fluid must be removed from the affected joint.

1. Twist and pull out

2. Insert and twist the cannula

**WRONG!**  
(Backstop opening at front)


**CORRECT!**  
(Backstop opening at back)

**Storage:**

Store Re-Cross by Reviscon between 2°C and 25°C. Protect from light and shocks. Do not freeze.

## Symbols:

	Manufacturer
	Do not reuse
	Do not resterilize
	Sterilized using steam (Product)
	Consult instructions for use
	Caution
	Keep away from sunlight
	Keep dry
	Do not use if package is damaged
	Use by date
	Latex free
	CE certified
	Storage temperature range
	Batch code

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VSY Biotechnology GmbH