GILCHRIST-BANDAGE

- (D) zur Immobilisierung der Schulter
- (GB) for immobilisation of shoulder
- (NL) voor het immobiliseren van de schouder
- (DK) til immobilisering af skulderen
- (F) pour l'immobilisation de l'épaule
- (I) per l'immobilizzazione della spalla
- (ES) para inmovilizar el hombro
- (P) para a imobilização o ombro
- (RU) Для иммобилизации плеча
- (SA) مثبت حركة الكتف

Gebrauchsinformation
Information for use
Gebruiksinformatie
Brugsinformation
Notice d'utilisation
Istruzioni per l'uso
Información de uso
Informações de uso
Инструкция применения
دواعى الاستعمال

D Anlegetechnik:

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GB Application:

Roll out the GILCHRIST-BANDAGE. Slip the elastic stockinet sleeve onto the injured arm allowing the attached shoulder padding to be positioned on the shoulder.

The shoulder padding should cover the shoulder joint. Ensure that the stockinet sleeve is placed as near as possible to the armpit.

Starting from the shoulder padding end, the bandage is then placed around the neck and onto the healthy shoulder.

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GB

One piece bandage ready for use for immobilisation of the upper arm and shoulder joint.

Range of application:

- Fractures in the upper arm region (e.g. subcapital fracture of the humerus)
- luxation of upper arm
- dislocation of the acromioclavicular joint
(acromioclavicular dislocation or shoulder joint dislocation): Tossi I or Tossi II
- postoperative care following osteosynthesis in the upper arm area
- postoperative care following surgical treatment of acromioclavicular dislocation: Tossi III

Contraindications:

With severe dislocation of the bone fragments which can only be reduced by operation.

Mode of operation:

In the event of a fracture, the bone fragments must be immobilised to ensure successful healing of the bone. Immobilisation also helps to reduce pain. Immobilising the upper arm in the case of torn ligaments in the area of the shoulder joint regenerates and scars the ligament tissue. Using GILCHRIST-BANDAGE for postoperative care stabilises the area which has been operated and reduces pain.

All serious incidents that occur in connection with the medical device must be reported to the manufacturer and to the Competent Authority of the Member State in which the user and/or patient is established.

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NOBAMED Paul Danz AG