The LumiraDX INR Quality Controls (hereafter referred to as "Quality Controls") are optional liquid quality controls to be used with the LumiraDX Instrument (hereafter referred to as "the Instrument") and the LumiraDX INR Test (hereafter referred to as "INR Test").

Read these instructions thoroughly before using the Quality Controls.

Inspect the Quality Controls packaging and contents for damage before use. Report any damage to LumiraDX Customer Services and do not use this kit if any damage is observed to the contents. Do not Perform Quality Controls testing if the Quality Controls material has not been fully reconstituted.

To ensure that you are using the Instrument, the INR Test and the Quality Control correctly, read the appropriate Platform User Manual, INR Test Product Insert and this entire pack insert. In addition, please watch the LumiraDX Platform training video available at [kcl.lumiradx.com](http://kcl.lumiradx.com). The Quality Controls are intended for professional use only.

The LumiraDX Logo, LumiraDX are trademarks of the LumiraDX Group.

#### Intended Use

The LumiraDX INR Quality Controls are intended for liquid quality control testing performed on the LumiraDX instrument when used with the LumiraDX INR Test Strip. The Quality Controls provide users with assurance that the device is performing within specification.

#### Summary and Explanation of the Test

The LumiraDX INR Quality Controls are an optional quality control for the instrument when used with the LumiraDX INR Test. The Quality Control material is in a device intended for medical purposes for use in a test system to estimate test precision and to detect systematic analytical deviations that may arise from reagent or analytical instrument variation and may be used for proficiency testing.

#### Reagents

Each Quality Control kit contains lyophilized human plasma with varied levels of coagulation factors, stabilizers, preservatives and a calcium chloride Reconstitution Solution. The Quality Control ranges are assigned by the LumiraDX Test Strip. The LumiraDX INR Test Strip performance is traceable to the WHO International Reference standard rITF/16.

#### Warnings and Precautions

For in vitro diagnostic use.

- The plasma control contains human source material that was tested and found nonreactive for HIV antibody, Hepatitis B Surface Antigen and Anti-HCV at the donor stage. This product, as with all human based specimens, should be treated as potentially infectious and handled with proper laboratory safety procedures to minimize the risk of transmission of infectious disease.
- All components of this kit can be discarded as Biohazard waste according to the local guidelines.
- Refer to the product safety data sheet for risk and safety phrases and disposal information.
- The product safety data sheet is available for users upon request.
- Requirements of the appropriate licensing or accrediting body should be incorporated into your quality control program.
- Exercise the normal precautions required for handling all laboratory reagents.

#### Storage and stability

• Store controls between 2°C and 8°C (36–46°F). DO NOT FREEZE.

• Allow the Quality Controls material to come up to room temperature (15–30°C; 59–86°F) for at least 10 minutes (but not above two (2) hours) before preparing its reconstitution.

• Unopened, lyophilized controls that are stored between 2°C and 8°C (36–46°F) can be used until the expiration date.

• Controls are stable for 30 minutes between 15°C and 30°C (59–86°F) after adding the Reconstitution Solution. Do not refrigerate after reconstitution.

#### Carton Contents

- 3 vials Level 1 Quality Control lyophilized plasma (normal INR levels) (Blue Cap)
- 3 vials Level 2 Quality Control lyophilized plasma (therapeutic INR levels) (Red Cap)
- 6 vials of Reconstitution Solution (White Cap)
- 6 large bulb pipettes
- 6 small bulb pipettes

#### INR Quality Controls Pack Insert

- Materials required but not provided with the Control Carton:**
- LumiraDX Instrument
  - LumiraDX INR Test Strips
  - LumiraDX Connect- if connectivity required (refer to LumiraDX Connect User Manual)
  - Quality Control Ranges Product Insert as included in the INR Test Strip Carton

#### 1. Getting ready to test

You will need the LumiraDX Instrument and the following supplies:

- LumiraDX INR Test Strip(s)
- LumiraDX INR Quality Controls Level 1 or Level 2
- 1 Vial of Reconstitution Solution
- Large bulb pipette
- Small bulb pipette

#### Performance characteristics

Quality Control precision was determined in measurement system analysis for the Quality Controls with the LumiraDX INR Test. The results were generated over multiple days by multiple operators and instruments. The results are shown in INR.

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**CE** CE mark applies to LumiraDX Instrument, Test Strips, Quality Controls and Connect Hub only.

#### 2. Preparing the Quality Controls

##### Quality Control Level 1

Mean	SD	%CV	N
1.25	0.081	6.5	235

##### Quality Control Level 2

Mean	SD	%CV	N
4.77	0.310	6.5	235

#### Limitations

This product is designed as a liquid control for quality control performed exclusively by the LumiraDX INR Platform. The quality control result is subjected to the limitations of the test Platform. Deviations may indicate potential problems with one or more components in the test Platform. The LumiraDX Instrument and INR Test Strips have on-board controls to detect errors and prevent false INR results when analysis is performed. Therefore deviations observed when testing with the LumiraDX INR Quality Controls may not invalidate previous results obtained from LumiraDX INR tests.

#### LumiraDX Customer Service

If you have questions call the LumiraDX Customer Services at

00800 58647239.

Any adverse results experienced with the use of this product and/or quality problems should also be reported to LumiraDX Customer Services at the above number or by email: [customerservices@lumiradx.com](mailto:customerservices@lumiradx.com).

#### Usage

The customer international free phone number will be free from most landlines in countries where the LumiraDX Platform is commercially available, but could be classified as an international call from a mobile and therefore chargeable. The charge will depend on the country you are calling from.

#### For Return Policy

If there is a problem with the LumiraDX INR Quality Controls you may be asked to return them. Before returning product please obtain a return authorization number from LumiraDX Customer Services. This return authorization number must be on the shipping carton for return. For ordinary returns following purchase, please contact Customer Services for terms and conditions.

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