## Technische Dokumentation - Tülip

### Akustische Performance von Tülip

#### Einführung

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## Acoustic Performance of Tülip

### Introduction

In addition to features such as durability and easy cleaning, the acoustic quality is particularly important in the development of a new stethoscope. In order to ensure the best quality in terms of auscultation for our stethoscope Tülip, we have collaborated with the Technical University of Cologne to develop a setup that allows for extensive acoustic testing.

### Requirements for a stethoscope in terms of auscultation:

- The instrument should convert a high proportion of mechanical energy from the body's interior into sound waves so that they can be perceived by the user
- The stethoscope must not attenuate transmitted energy, as the original energy of body sounds is already very low  
→ Further attenuation could mask important or characteristic sounds

### Body sounds and their frequency

*Table 1: Body sounds and their frequency range, own representation*

| Sounds   | Frequency range |
|--|-----------------|
| Gastrointestinal   | 97 - 1034 Hz    |
| Cardiac (murmurs, mitral stenosis)   | 22 - 281 Hz     |
| Cardiac (ventricular septal defect, atrial septal defect, aortic stenosis and pulmonic stenosis) | 205 - 775 Hz    |
| Respiratory  | 689 - 2584 Hz   |

## Materials & Methods

The approach of the test setup to demonstrate the quality of auscultation is based on the methods mentioned in Renoll, 2021 and Weiss, 2019. Stethoscopes are tested using an acoustic phantom that aims to replicate the acoustic properties of the human body. The figure shows a schematic overview of the test setup.


Figure 1: Schematic of the Test System, own representation

### Procedure:

- The measured frequency ranges from 20 Hz to 2 kHz (x-axis)  
→ Spectrum of frequencies of sounds in the human body that can be perceived with a stethoscope
- The y-axis represents the sound pressure level
- Each graph reflects the average of ten measurements
- Ballistic gelatin represents human skin
- 200g and 1kg weights represent the pressure exerted during auscultation
- A high sensitivity microphone is used to record the response of the diaphragm
- The recording takes place in an isolated chamber


Figure 2: Placing the chest piece on a block of gelatin, self-captured

Figure 3: Insertion of the earloops into the measuring device, self-captured

## Results


Figure 3: Graphic representation of the results from the acoustic testing of the tested stethoscopes, each with a 1kg weight, own representation

Table 3: Results of the tested stethoscopes depending on the frequency range, own representation

| Stethoscope  | Representation | <160Hz   | >160Hz  | >1kHz                                 |
|--|----------------|--|---|---------------------------------------|
| <b>Tülip</b>                                       | Green line     | Delta of 20 dB; max. SPL at 90 dB                | Curve remains relatively constant with occasional small peaks, SPL ranging between 70 and 85 dB | Curve drops sharply to 25 dB at 2 kHz |
| <b>Tülip with Cøpper</b>                           | Red line       | Delta of 18 dB; max. SPL at 83 dB                | Relatively constant, drops at 550 Hz  | Curve continues to decline            |
| <b>Littmann Classic III, Reference-stethoscope</b> | Grey line      | Delta of 35 dB; significant dips → Min. at 50 dB | Constant values with a peak at 90 dB; SPL ranging between 75 and 90 dB.                         | Curve flattens out significantly      |

## Analysis

The graph shows that Tülip has consistent values with minimal dips. Especially in the low-frequency range (<160Hz), the sound pressure level (SPL) of Tülip is at least 10 dB higher than that of the reference stethoscope. In the frequency range >160Hz, the reference stethoscope shows higher values compared to Tülip. The overall difference in the entire frequency range is on average only 2 dB. Higher frequencies are recorded when using the 1kg weight.

The graph of Tülip with Cøpper shows the lowest levels. In the low-frequency range, the values still remain above those of the Littmann Classic III. It is noticeable that the sound pressure level drops already at 500Hz.

## Interpretation

A high SPL serves as an indicator of the efficiency of converting mechanical energy into acoustic energy. A loss in magnitude in certain areas of the spectrum can result in critical sounds being missed during auscultation. Low levels can be caused by factors such as too stiff diaphragms, hose blockage, non-airtight stethoscope, or ill-fitting ear tips.

Tülip consistently delivers high SPL values with minimal dips, indicating that it can perceive all relevant sounds of the human body (as shown in Table 1). In particular, heart sounds can be auscultated better with Tülip compared to the reference stethoscope. Applying more pressure on the chest piece improves the quality of auscultation.

Auscultation is a highly subjective perception, where the experience of the user plays an important role. Therefore, the purely technical evaluation does not provide a complete picture of the quality for the user. However, the acoustic performance of the stethoscope is a significant starting point to ensure that the stethoscope supports the user's work to the best possible extent.

## References

1. **Rennoll V**, McLane I, Emmanouilidou D, West J, Elhilali M. Electronic Stethoscope Filtering Mimics the Perceived Sound Characteristics of Acoustic Stethoscope. *IEEE J Biomed Health Inform.* 2021 May;25(5):1542-1549. doi: 10.1109/JBHI.2020.3020494. Epub 2021 May 11. PMID: 32870803; PMCID:PMC7917155. <https://doi.org/10.1109/jbhi.2020.3020494>
2. **Weiss D**, Erie C, Butera J 3rd, Copt R, Yeaw G, Harpster M, Hughes J, Salem DN. An in vitro acoustic analysis and comparison of popular stethoscopes. *Med Devices (Auckl).* 2019 Jan 15;12:41-52. doi: 10.2147/MDER.S186076. PMID: 30697087; PMCID: PMC6339642. <https://doi.org/10.2147/mder.s186076>