

MICRO-TOUCH[®] Nitrile N30

Standard Thin

Lightweight nitrile, medical exam-grade gloves in a convenient, 30-count retail pack

- **Convenient pack size:** Perfectly sized for retail shelves or places that have limited storage with 30 gloves per box
- **Versatile protection:** Suitable for low-risk medical procedures, food handling or general household tasks
- **Easy donning:** An interior polymer coating makes donning easy
- **Heightened tactility:** Thin, lightweight design enhances tactility and dexterity
- **Reduced allergy risk:** Nitrile material prevents Type 1 (latex) allergic reactions

Recommended For

- Brief, low-risk examination procedures
- General patient care
- Protection from Type I latex allergy in HCW's or patients
- Protection from viruses
- Multipurpose use: Dish Washing, Cleaning, Pet Care & Grooming, Food Preparation, Gardening & many more!

FEATURED TECHNOLOGIES



MICRO-TOUCH[®] Everyday, Everywhere Protection

MICRO-TOUCH[®] premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH[®] Nitrile N30

Standard Thin

PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Blue	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.07 / 2.8
Finger Thickness (mm/mil)	0.11 / 4.3	Allergy Prevention	Latex (Type I)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	No
Sterile	No	Antistatic	Not Tested
Glove Length (mm/inches)	250 / 9.8	Product Segmentation	Standard Thin

ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	313043060	313043065	313043070	313043075	313043080

PACKAGING AND STORAGE

Packaging	30 gloves per dispenser; 50 dispensers per case; 1500 gloves per case
Shelf Life	5 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

Category III, CE 0493, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, Food Type: Rubber, ISO 9001, ISO 13485, ISO 10993-10:2002

Vulcanization Chemical Accelerators

Zinc Dibutyldithiocarbamate (ZDBC)

Only a very small number of users may be sensitive to this ingredient(s) and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email:
insidesalesus@ansell.com

Latin America & Caribbean
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil Tel: +55-11-3356-3100
Brazil Email:
luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email:
apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email:
protection@ap.ansell.com