

ENCORE[®] Latex Micro

S U R G I C A L G L O V E S

Latex, powder free micro surgical gloves offer exceptional tactile sensitivity and dexterity

- **Heightened dexterity:** Being 20% thinner than Ansell's standard powder-free latex gloves, ENCORE[®] Latex Micro offer exceptional level of tactile sensitivity and dexterity preferred for delicate procedures
- **Elevated ergonomic comfort:** Ansell's unique anatomic glove formers deliver surgical gloves designed for enhanced comfort and movement to help reduce hand fatigue
- **Double-gloving versatility:** This micro glove's brown color diminishes surgical glare while offering the wearer the option to use it as an inner glove when double gloving for easy breach detection
- **Enhanced donning aid:** A unique internal coating enables easy donning and doffing under dry or damp hand conditions, be it the first don or during intra-operative changes



CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

Recommended For

- All surgeries, especially delicate procedures in ophthalmology, micro, and cardiovascular
- Damp donning
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Prevention of eye fatigue

FEATURED TECHNOLOGIES



ENCORE[®] The Confidence to Perform

Dedicated to delivering optimized solutions engineered to provide a high degree of confidence to surgeons, O.R. nurses and surgical technicians seeking protection with a high level of dexterity, comfort and grip so they can perform at their best.

ENCORE[®] Latex Micro

SURGICAL GLOVES

PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Brown
Thickness	Micro
Cuff Style	Beaded with SUREFIT™ Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology
Allergy Prevention	None
Grip Level	Moderate
Double Gloving Recommendation	Outer glove or underglove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	• In process control before packing: 0.65 (Inspection level GI) • Manufacturing final release: 0.65 (Inspection level GI)
Sterile	Yes
Product Standards	ASTM D3577, ASTM D7160, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN 374-2 and -4, EN ISO 374-5:2016, ISO 10282
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Dibutyldithiocarbamate (ZDBC), Zinc Diethyldithiocarbamate (ZDEC)

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	330104055	330104060	330104065	330104070	330104075	330104080	330104085	330104090

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 844 494 7854
Email: insidesalesus@ansell.com
Sandel US Tel: 866 764 3327
Sandel email: sandel-cs@ansell.com

Central & South America
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil
Tel: +55-11-3356-3100
Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

India
Tel: +91 22 38124500
Email: info.india@ansell.com
Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com
China
Tel: +86(21) 3827 5005
Email: infochina@ansell.com