

1m4 Protect

Care product for stethoscope tubes

- ▶ Nourishes plastics and protects against material aging
- ▶ Restores flexibility to tubes and refreshes colors
- ▶ Water and dirt repellent
- ▶ Very economical
- ▶ Content: 50 ml
- ▶ For professional use only



Technical specifications

Dimension:	90 x Ø 40 mm
Weight:	75 g
HS-Code*:	34053000 DE I for CH: 100% VOC

Product data

Product	Retail packaging	Outer packaging	EAN	REF-No.
Protect	1 piece	25 x 1 piece	4745010564048	10002.00



Signal word: Danger

Hazard statements:

H225	Highly flammable liquid and vapour.
H304	May be fatal if swallowed and enters airways.
H315	Causes skin irritation.
H336	May cause drowsiness or dizziness.
H411	Toxic to aquatic life with long lasting effects.

Precautionary statements:

Prevention:

P210 Keep away from heat, hot surfaces, sparks, open flames, and other ignition sources. No smoking.
P273 Avoid release to the environment.

Reaktion:

P301 + P310 IF SWALLOWED: Immediately call a POISON CENTER/doctor.
P331 Do NOT induce vomiting.
P370 + P378 In case of fire: Use water spray, alcohol-resistant foam, dry chemical, or carbon dioxide to extinguish.
P391 Collect spillage.

Sales:

medichem Vertriebs GmbH • Sandhof 8-10 • 24768 Rendsburg • Germany
Tel.: +49 4331-72 18 5 • Fax: +49 4331-78 09 46 • E-Mail: info@medichem-online.de