



Disinfection on your desk

The first desk device stethoscope disinfection, specifically designed to ensure security and deep hygiene to patients thanks to the action of UV-C LED rays.




Always within reach

The stethoscope is the most used medical tool in the world. Despite this fact, only 10-20% of sanitary personnel regularly clean their stethoscope before using it on each patient. This bad practice exponentially increase the contamination of the tool after every use and the possibility of transmitting bacterias and associated illnesses.

STET CUBE, with its innovative and practical design conceived to be always ready to use on the studio desk, represents an effective and efficient solution to these problems.

Why an UV-C LED device



STETCUBE

TECHNOLOGICAL FEATURES

- innovative UV-C LED technology device
- automatic double level of treatment: 3 minutes for standard disinfection, 2 minutes more for a deeper disinfection
- compatible with every stethoscope type and dimension (pediatric, neonatal, cardiological...)
- microprocessor for irradiation and security controls
- special polycarbonate body
- status light indicator: ongoing treatment, treatment complete, low battery, malfunction
- rechargeable battery with standard micro USB cable
- high battery autonomy

