



SYREBOTM

Soft Rehabilitation Robotic Gloves

Gain Beautiful life with your hands!


Do you have the following problems?



Unable to take
objects alone



Unable to
eat/drink alone



Unable to dress
yourself alone



Unable to read
book/play phone alone

Syrebo™ Soft Rehabilitation Gloves is Your Best Companion!

Suitable for:

patients with hand dysfunction caused by hemiplegia /cerebral infarction/ stroke/ peripheral nerve injury of hand/ orthopedic surgery/ spinal cord injury/ hand burn/ children cerebral palsy, etc



● Passive Flexion & Extension

The rehabilitation glove will drive the affected hand to grasp the small ball/water cup and other objects.

Patients could intelligently adjust the time of flexion and extension.



● Innovative Mirror Training

With the help of syrebo™ rehabilitation glove, the healthy hand drives the affected hand to move synchronously, activates the mirror neurons, promotes the autonomic recovery of the brain and speeds up the rehabilitation process of hand function.





- **Active Game Training**

Game training could make rehabilitation more interesting, and directly activate the center nerve, accelerate the rehabilitation of stroke patients' brain motor intention control output.

It can be connected to App for operation (can be screened to TV)

- Intelligent Touch Button



Humanized interactive design

Elder could operate independently

Wireless portable, small size, big energy

Soft Pneumatic Exoskeleton Glove 🧤

Lightweight & Easy-wear 🧘

Adjustable Tightening Design 🔧

Safety & Comfort 🛡️



● Soft Hand Rehabilitation Robotic Gloves Series



C10



C11

	C10	C11
Passive Training	✓	✓
Mirror Training	✓	✓
Active Training		✓
Task-oriented Training	✓	✓

Model	Hand length (in.)	Hand width (in.)
#5 (child)	4.3~5.0	2.0~2.2
#7 (child)	4.7~5.3	2.2~2.4
#9 (child)	5.0~5.5	2.4~2.6
#12 (child)	5.4~6.0	2.6~2.8
S	5.7~6.5	2.8~3.1
M	6.3~7.3	3.1~3.7
L	7.1~7.8	3.6~4.1
XL	7.7~8.3	3.9~4.6

C11 could be connected to App for game training, currently only for Android.



• Patients' Feedback



Ying Bosong's son: 36 years old

Cause of illness: Right hand dysfunction due to cerebral infarction in March 2019.

Feeling of using: After more than 3 months of use, the hand has been significantly improved, and the grasping ability of the hand has also been significantly improved.



Wang Chunsen: 62 years old

Cause of illness: Cerebral hemorrhage in 2011, left hand hemiplegia, arm can not be raised, finger weakness.

Feeling of using: Through the introduction of a friend, I bought the rehabilitation robot gloves. After several months of rehabilitation training, my fingers, which have no strength for several years, do have some strength. when wear the gloves, I can hold some small objects. Now I have confidence in recovery.



Zhang Jian's father: 60 years old

Cause of illness: A brain hemorrhage occurred in April 2019.

Feeling of using: Now my father's hands are recovering very well, his fingers are not so stiff, and his normal straightening and clenching fists are also much more flexible. Compared with the patients who were sick at the same time, my father recovered better than them.




Han Chen: 38 years old

Cause of illness: In 2018, sudden hypertension caused cerebral hemorrhage.

Feeling of using: After using Syrebo™ Soft Rehabilitation Gloves for one and a half months, the immobile hands can now move, and simple grasps can be performed, the fingers can move, the hands can turn the wrists, and the effect is quite good. I am also full of confidence and expectation for this machine.

- Syrebo™ hand rehabilitation gloves series are applied in thousands of hospitals and families, and its effect is recognized by lots of institution



WHO WE ARE?

Shanghai SiYi Intelligent Technology Co., Ltd.

Shanghai SiYi Intelligent Technology Co.,Ltd, located in Shanghai International Medical Zone, is a high-tech enterprise specializing in the development of medical rehabilitation robots. SIYI INTELLIGENCE is the leader of software rehabilitation robot. We are committed to continuously creating a series of intelligent rehabilitation equipment from hospital to home, making unremitting efforts to comprehensively improve the quality of life of patients.



FDA certification



Clinical recognition
of hundreds of hospitals



More than 30 national
patent technologies



Factory specialized
production



Cutting-edge technology China's
first soft exoskeleton robot



Expert message:

"Hand dysfunction brings great difficulties to patients' life, but the rehabilitation process is often difficult. This is the status quo that can not be changed. To achieve the best rehabilitation effect, the most important thing is confidence, scientific rehabilitation training and perseverance. "



Hand Rehabilitation system

Shanghai SiYi Intelligent Technology Co.,Ltd

Office: 400-6886-275

Mobile&whatsapp: (+86) 17317983724

E-mail: wmsales@siyizn.com

Website: www.syrebo.com



WhatsApp



WeChat