

# MICRO-TOUCH<sup>®</sup> Nitrile Sterile

P e r f o r m a n c e   S e r i e s

## Nitrile, sterile, powder-free medical exam glove featuring extended cuff, packed as a pair


- **Barrier protection:** Packed as a pair of gloves per pack, these sterile gloves are made from a synthetic, nitrile formulation that provides high elasticity with superior barrier protection and puncture resistance
- **Sterile:** With a minimum 25kGy to provide sterility, they are an ideal solution for a range of sterile medical exam (non-surgical) procedures
- **Extended protection:** The 11.6" extended cuff provides additional protection
- **Secure grip:** Textured fingers provide extra grip and reliability
- **Reduced allergy risk:** Nitrile material prevents latex Type I allergic reactions
- **Specialized protection:** Tested and approved for use with chemotherapy drugs\*

\*see product packaging for details

### Recommended For

- Tested for use with chemotherapy drugs
- Sterile procedures outside of the O.R.
- Extra protection over the wrist and arm
- Sustained, high-risk examination procedures
- Protection from Type I latex allergy in HCW's or patients

### FEATURED TECHNOLOGIES



### MICRO-TOUCH<sup>®</sup> Everyday, Everywhere Protection

MICRO-TOUCH<sup>®</sup> premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH<sup>®</sup> Nitrile Sterile

Performance Series

## PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Blue	Cuff Length	Extended
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.11 / 4.3
Finger Thickness (mm/mil)	0.16 / 6.3	Allergy Prevention	Latex (Type I)
Available Sizes	S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Sterile	Yes	Antistatic	Not Tested
Glove Length (mm/inches)	290 / 11.6	Product Segmentation	Performance Series

## ORDERING INFORMATION

Size	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	6034151	6034152	6034153	6034154

## PACKAGING AND STORAGE

Packaging	1 pair per pouch 50 pairs per dispenser 4 dispensers per case 200 pairs per case
Shelf Life	5 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

## STANDARDS AND CERTIFICATIONS

ASTM D6978, ASTM D6319, Category III, CE 2797, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 13485, ISO 9001, ISO 10993-10:2002

## Vulcanization Chemical Accelerators

Zinc Dibutyldithiocarbamate (ZDBC)

Only a very small number of users may be sensitive to this ingredient(s) and hence may develop irritant and/or allergic contact reactions.

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email:  
insidesalesus@ansell.com

Latin America & Caribbean  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email:  
luvas.medicas@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email:  
apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email:  
protection@ap.ansell.com