



# Linktop Stethoscope

**Perform the professional heart and lung examination at the comfort of your home.**

## BENEFITS

### Examination quality control

The system helps carry out the examination correctly: a step by step tutorial will guide you. Excess noise will be monitored automatically, as it may impact the results.

### Easy contact with a doctor

After finishing the examination you can share the data with your doctor via e-mail or Social applications and consult them without leaving home.

### Wireless stethoscope

The stethoscope connects to a mobile device via Bluetooth without the need of device pairing.

### Noise control

Informs of excessive noise in the room.

### Automatic start/stop

Records sounds automatically once the stethoscope is put to the body correctly.

### Telemedicine

Enables remote consultation irrespective of where you are.

# SPECIFICATIONS

## Size

- **Product Size:**
  - 45 x 45 x 26 mm (H x W x D)
  - 1.77 x 1.77 x 1.02 inches (H x W x D)
- **Weight:** 78 g / 2.75 oz / 0.17 lb

## Requirements

- Android 6.0 or Higher
- iOS 11 or Higher
- The authorized application
- Bluetooth Smart
- Online connection

## The box contains:

- Linktop Stethoscope
- Magnetic USB charging cable
- Instruction manual
- The authorized 3.5mm Headphone (Optional)
- USB – mini jack recording cable (Optional)

The box content depends on the offer selected.

## Probe

- Chest piece Finish: Matte Finish
- Chest piece Technology: Single sided
- Chest piece Material: Kirsite
- Diaphragm Diameter: 43.8 mm / 1.7 inches
- Diaphragm Material
  - Polyurethane-Coated Silicone
  - Non-chill diaphragm cover
- Digital Electronic Filtering - Ambient Noise Reduction (ANR)

## User interface

- Product Color: White
- Suitable for: Adult, Infant, Pediatric
- Pickup frequency range: 20Hz-3KHz
- Eliminate 85% (on average) of ambient noise
- Amplify sounds up to 24x
- Frequency modes
  - Diaphragm
  - Bell
- Auscultation mode
  - Listen remotely via NexStetho App
  - Listen via 3.5mm headset
- Single key
  - Power On/Off
  - Switch Diaphragm/Bell mode
- Blue and Green LED indicators
- Bluetooth 5.0

## Power and Interface

- Power requirements: 3.3V~4.2V
- Battery capacity: 200mAh
- The power adapter: DC 5V/500mA
- Magnetic USB charging interface
- 1 x 3.5mm Headphone
- 1 x 3.5mm Recording

## Environmental requirements

- Working temperature: -10°C -40°C (14-113 °F)
- Relative humidity: 5-95%, non-condensing
- Storage temperature: -25°C~70°C

## Warranty

- 2-year limited hardware

## LEARN MORE

For more information on Linktop Stethoscope [www.linktop.com](http://www.linktop.com)

