



GIMA

FFP2 FILTERING MASK - black - GB,FR,IT,ES,PT

Code: 20735
Category: Face masks
Unit of sale: box of 20 pcs.
Minimum order: 1
Type: PPE
Class: III



EAN13: 8023279207354

Description: VALUE&SAFE FFP2 FILTERING MASKS WITH EAR LOOPS

- FFP2 without valve
- Languages: GB, FR, IT, ES, PT
- Colour: black
- Box of 20

Five-layer filtering FFP2 mask to protect yourself and others effectively.

Meeting FFP2 requirements of EN 149:2001+A1:2009 in accordance with the PPE regulation EU 2016/425.

- ergonomic design fit face closely
- non irritating, soft and skin-friendly inside non woven fabric
- adjustable nasal clip
- flexible and stretchable ear loop, painless after long-time wearing
- foldable in half for easy gathering
- each mask is individually packed in a transparent pouch

Technical Specifications:
Inhalation: 30 l/min% = 0.7
Inhalation: 95 l/min% = 2.4
Exhalation: 160 l/min% = 3.0
Penetration of filter material: = 6%
Size unfolded: 21 x 15.5 cm
folded: 10.5 x 15.5 cm