



Choose the right stethoscope to help you perform at your best.

3M™ Littmann® Stethoscopes feature adjustable headsets with soft-sealing eartips. Select Littmann stethoscopes are designed with tunable technology to give you the ability to hear different frequency sounds by simply adjusting the pressure applied to the chestpiece. Get the reliable acoustic performance, high quality and comfort you need with a Littmann stethoscope.



	Lightweight II S.E. Stethoscope	Classic II Pediatric and Infant Stethoscopes	Classic III™ Monitoring Stethoscope	Cardiology IV™ Diagnostic Stethoscope	Master Cardiology™ Stethoscope	CORE Digital Stethoscope
	Ideal for taking blood pressure readings and making limited patient assessments.	Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients, one long shift after another.		Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation. Useful in challenging and critical care environments.		High-performance auscultation. Hear more with amplification, noise cancellation and mobile connectivity.
Approximate Length	28"/71cm	28"/71cm	27"/69cm	22"/56cm, 27"/69cm	27"/68cm	27"/69cm
Approximate Weight	118g	105g, 95g	150g	167g, 177g	185g	232g
Warranty	2 years	3 years	5 years	7 years	7 years	2 years
Diaphragm	Tunable	Non-Tunable (Floating)	Tunable (2)	Tunable (2)	Tunable	Tunable (2)
Approximate Chestpiece Diameter	Double-sided 2.0 x 1.8"/5.2 x 4.5cm 1.4 x 1.3"/3.5 x 3.2cm	Double-sided Pediatric 1.4"/3.4cm 1.1"/2.7cm Infant 1.1"/2.8cm 0.7"/1.7cm	Double-sided Adult side 1.7"/4.3cm Pediatric side 1.4"/3.4cm	Double-sided Adult side 1.7"/4.3cm Pediatric side 1.4"/3.4cm	Single-sided 1.9"/4.8cm	Double-sided Adult side 1.7"/4.3cm Pediatric side 1.4"/3.4cm
Features	Teardrop-shaped chestpiece Fits under blood pressure cuff	Specially designed for the smallest patients	Convertible open bell	Convertible open bell Dual-lumen tubing	Dual-lumen tubing Best acoustic response in Littmann mechanical line	Up to 40x amplification* Switch from analog listening to digital Connects to Eko software on a smart device** to capture, store and share data



The power to learn, in the palm of your hand.

3M™ Littmann® Learning App – Download free version in the App Store® or Google Play™ Store

Need more help? Visit Littmann.com.au/FindMyStethoscope
Need more help? Visit Littmann.co.nz/FindMyStethoscope

1300 363 878 | www.littmann.com.au
0 800 80 8182 | www.littmann.co.nz

*At peak frequency, vs. analog mode
**Not included