

# GAMMEX® Latex Moisturizing

S U R G I C A L   G L O V E S

## Latex powder-free surgical gloves, featuring a built-in moisturizer to rehydrate the skin

- **Soft, supple skin:** GAMMEX® Latex Moisturizing surgical gloves feature a proprietary skin moisturizing technology, enabling a water-based coating with glycerol and dimethicone to retain moisture and rehydrate dry skin, improving the wearer's skin condition
- **Exceptional comfort:** These gloves are also designed using anatomic formers for added comfort and fit
- **Enhanced grip:** With a smooth finish and micro-textured external surface around the palm and fingers, these surgical gloves deliver excellent grip under dry or damp conditions
- **Easy donning:** Ansell's donning aid technology facilitates an easy first don, whether the wearer's hands are dry or damp
- **Optimal allergy prevention:** Additional cleaning processes ensure low residues of latex proteins and chemical accelerators, minimizing the risk of developing latex Type I and chemical Type IV allergies and sensitivities



**CAUTION:** Safe use of these gloves by latex-sensitized individuals has not been established - this product contains natural rubber latex, which may cause allergic reactions

### Recommended For

- Damp donning
- Extra comfort during long, sustained procedures
- Use as an underglove when double gloving

### FEATURED TECHNOLOGIES



### GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> Latex Moisturizing

SURGICAL GLOVES

## PRODUCT INFORMATION

<b>Material</b>	Natural Rubber Latex
<b>Powder Content</b>	Powder-Free
<b>Color</b>	Natural
<b>Thickness</b>	Standard
<b>Cuff Style</b>	Straight with SUREFIT™ Technology
<b>External Glove Surface</b>	Smooth with Micro-Textured Finish
<b>Internal Glove Surface</b>	Coated with HYDRASOFT™ Technology
<b>Allergy Prevention</b>	None
<b>Grip Level</b>	Moderate
<b>Double Gloving Recommendation</b>	Underglove
<b>Tested For Use With Chemotherapy Drugs</b>	No, not tested in accordance with ASTM D6978
<b>Freedom from Holes (Inspection level I)</b>	0.65 AQL
<b>Sterile</b>	Yes
<b>Product Standards</b>	AS/NZS 4179, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107
<b>Quality / Environmental Standards</b>	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
<b>Packaging</b>	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
<b>Shelf Life</b>	3 years
<b>Storage Instructions</b>	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
<b>Vulcanization Chemical Accelerators</b>	Xanthogen Polysulphides, Sodium di-n-butylthiocarbamate (SDBC), Zinc Diethyldithiocarbamate (ZDEC)

## ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
<b>Product Code</b>	330049055	330049060	330049065	330049070	330049075	330049080	330049085	330049090

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: insidesalesus@ansell.com  
 Sandel US Tel: 866 764 3327  
 Sandel email: sandel-cs@ansell.com

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: cslac@ansell.com  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: luvas.medicas@ansell.com

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: info@ansell.eu

**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: apac.medical@ansell.com

**India**  
 Tel: +91 22 38124500  
 Email: info.india@ansell.com  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: protection@ap.ansell.com  
**China**  
 Tel: +86(21) 3827 5005  
 Email: infochina@ansell.com